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Instruction for use English

Definition

edelweiss PEDIATRIC CROWNS are laser-sintered, pre-polymerized, highly filled nano-hybrid composite enamel shells. Their anatomical profiles are replicas of the natural deciduous tooth in the form of highly esthetic and improved physical properties that facilitates their direct use as individual crowns for anterior and posterior deciduous teeth.

edelweiss PEDIATRIC CROWNS are pre-fabricated, radio-opaque PEDI-ATRIC CROWNS that form the anatomical basis of the enamel layer of the tooth and together with edelweiss COMPOSITE can be used in cases of shallow or deep caries for simple or complex restorations of both anterior and the posterior deciduous teeth. The PEDIATRIC CROWNS provide complete protection of the deciduous tooth up until exfoliation while at the same time maintaining functional support through proper occlusal contacts with opposing teeth. Occlusion can be easily obtained by adjusting the PEDIATRIC CROWN and there is no need to adjust the opposing natural tooth.

Restoring the tooth is quick and easy and provides a bio-esthetic and bio-functional end result that mimics the natural dentition. Minimal preparation is required and is limited by the extent of the caries, hence there is no risk of iatrogenic pulpal exposure to the deciduous tooth.

PEDIATRIC CROWNS are cemented with the same material they are made of, which produces a stable intaglio of the entire restoration (between dentin/enamel + composite + PEDIATRIC CROWN). Edelweiss PE-DIATRIC CROWNS can also be cemented with conventional Glass Ionomer Cement or Resin Modified Glass Ionomer Cement (RMGIC).

Composition

Barium glass
Bis-GMA-based matrix made of dental resins
Pigments
Additives
Catalyst

Technical data

Filler content: 82 wt %
variation range of inorganic filler particles lies between 0.02-3 µm.

Features of edelweiss COMPOSITE

- Easy handling
- Excellent physical properties
- Easy and fast polishing
- Low shrinkage due to Nano-technology and high degree of fillers (82%)
- Excellent abrasion with natural behavior pattern
- Zinc and fluoride component provides antibacterial surface
- Natural fluorescence
- Radio-opaque

Indications

- Caries damage to tooth structure
- Replacing old defective restorations
- Restorations of all cavities from Class I to VI
- Mal-aligned teeth
- Tooth fractures
- Tooth discoloration
- Anatomical malformation
- Developmental defects
- Abrasion
- Erosion
- Cosmetic corrections
- Lifting of vertical dimension due to loss in the posterior region

Recommended application

- Sizing of PEDIATRIC CROWN The appropriate size of the PEDIATRIC CROWN is selected using the supplied sizing guide (edelweiss PEDIATRIC CROWN Sizing Guide). The sizing guide is placed over the deciduous teeth to be restored so that the mesio-distal margins becomes visible on the transparent sizing guide and the best matching (small, medium or large) edelweiss PEDI-ATRIC CROWN is then selected.
- Use of Rubber Dam is recommended whenever feasible.
- Tooth preparation

Tooth preparation is largely dictated by the extent of the caries and removed as indicated. edelweiss PEDIATRIC CROWNS are designed for a minimally-invasive preparation and it is recommended to preserve as much tooth structure as possible. The interproximal contacts should be reduced sufficiently to allow the selected crown to fit passively. This reduction should be performed using a tapered diamond or carbide bur on all planes of the tooth to follow the natural contours of the existing clinical crown. This would result in a preparation that is parallel to slightly tapered incisally/occlusally. edelweiss PEDIATRIC CROWN is tried over the prepped tooth and can be adjusted if necessary.

- Shade selection of dentin is done prior to tooth drying. Dentin shades (A0, A1, A2, A3, A3.5) selected using the appropriate VITA Shade Guide or edelweiss Shade Tabs. Prefabricated edelweiss PEDIATRIC CROWNS are translucent to mimic natural enamel.

5. Cementation of edelweiss PEDIATRIC CROWN

5.1. Cementation using edelweiss COMPOSITE

5.1.a. The inside of the crown is roughened very lightly with a diamond bur or sand-blasted with 25µm to 50µm aluminum oxide, rinsed off and air dried.

5.1.b. edelweiss VENEER Bond is applied to the inside of the crown using an applicator tip, lightly air dried to allow even spreading of the VENEER Bond and prevent pooling

5.1.c. Light-cure the edelweiss VENEER Bond for 10 seconds if using a 1200-2000 mW/cm² curing light, or 20 seconds if using a 800 mW/cm² curing light.

5.1.e. Etch the prepared tooth with 37% phosphoric acid etch gel for 15 seconds on enamel and 7 seconds on dentin. Rinse thoroughly for approximately 20 seconds and air blow gently with oil and water-free compressed air.

5.1.f. Dispense dentin adhesive onto an applicator tip and rub gently on the tooth surface for 20 seconds and light cure as per manufacturer’s instructions.

5.1.g. Cementing of edelweiss PEDIATRIC CROWN

The use of edelweiss Nano-hybrid composite is recommended to ensure esthetic and functional success. edelweiss COMPOSITE is applied to the edelweiss PEDIATRIC CROWN and positioned onto the prepared tooth, excess composite is removed and spot cured for 3 seconds. Final positioning checked and final cure achieved with 20 seconds around each wall and margins using 1200-2000 mW/cm² light curing unit. Polishing of margins achieved using green and white polishing cups.

5.2. Cementation using Resin Modified Glass Ionomer Cement edelweiss Following conditioning of the tooth surface, the inside of the crown is roughened using a diamond bur, rinsed off and air dried. When using RMGIC it is preferable to use edelweiss VENEER BOND on the inner surface of the crown. The cement is mixed as per manufacturer’s instructions, inserted into the crown and the crown positioned onto the prepared tooth. Excess cement is removed and light cured if necessary.

Precautionary Information

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known allergies to any of the ingredients. If irritation as a result of contact with oral soft tissue occurs, flush with large amounts of water and avoid usage. If allergic reaction occurs, seek medical attention as needed.

Contraindications

If known allergies to any of the components.

Side effects

Systemic side effects are not known. Individual cases of contact allergies have been reported for other products with a similar composition.

Interactions

Discoloration may occur during contact with cationic mouth washes as well as plaque disclosing agents or prolonged use of chlorhexidine.

Shelf-life and labeling

The expiry date and the batch description are given on the container.

Do not use after expiry date.

Important

To ensure optimal results, it is recommended to use in combination with all edelweiss components.

Note

May only be supplied to dentists and dental laboratories or on their behalf. Keep away from children! Only for use by dentists!

Date of issue of this package information leaflet 02-2020

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