**Eigenschaften des edelweiss Flowable****Nano-Hybrid Composite**

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**Instruction for use****English****Definition**

edelweiss POST & CORE system are laser-sintered and pre-polymerized highly filled nanohybrid composite post and core monoblocks for use when restoring endodontically treated teeth.

The post is a smooth surface tapered post that conforms to the internal post space design with a 1.4 mm and 1 mm diameter making it ideal for use in anterior and posterior teeth. The posts have a conical shape for perfect post space adaptation. The translucency of the post ensures uninterrupted light transmission for complete polymerization even in the apical region. The unique process of laser sintering allows for improved composite technology with superior mechanical properties.

The system consists of a prefabricated core attached to the post as a single unit. This unique and innovative monoblock design of the post and the core together in one unit provides complete stabilization within the root system. A core portion that corresponds to individual teeth in design is attached to the post. This is designed to simulate a crown preparation of the core material and thus saves the additional step of a core build-up. The attached core requires minimum preparation and can be adjusted for any type of restoration like veneer, inlay, onlay, crown or bridge. This makes core preparation easy and time saving. Reconstruction is quick and easy and provides an esthetic, functional and high quality outcome. The edelweiss POST & CORE ensures an highly esthetic result because of the inherent translucent post attached to an opaque core.

The POST & CORE is bonded to the root dentin with the same material it is made of, thus creating a single intaglio between root- dentine/cervical enamel + composite resin cement + POST & CORE. This allows for maximum retention of the post within the root canal.

**Composition**

Barium glass  
Bis-GMA-based matrix made of dental resins  
Pigments  
Additives  
Catalyst

**Technical data**

Filler content: 83 wt % Variation range of inorganic filler particles lies between 0.02-3 µm

**Features of edelweiss POST & CORE**

- Translucent tapered post allows for deep curing
- Opaque core allows for high esthetic restorations
- True monoblock effect enabling good retention
- Easy to use
- Cuts like dentin, can be easily adjusted in the mouth.
- Easy to customize
- Modulus of elasticity similar to natural tooth-preventing wedging effect
- Unmatched radiographic visibility: radio-opaque
- Superior adhesion to the root canal

**Features of edelweiss Flowable Nano-Hybrid Composite**

(recommended for POST&CORE cementation and core adjustment)

- Easy to use
- Excellent handling features
- Easy and quick application
- Superior physical features
- Antibacterial fillers
- Biocompatible

**Indications**

- Restoration of endodontically treated teeth
- Insufficient amount of remaining tooth structure to provide retention
- To provide support for tooth reconstruction
- Stabilization of restoration against pulling, bending and shearing forces
- Severe loss of clinical crown height as occurs in attrition or abrasion
- When the nonvital tooth is to be protected against fracture as the bridge abutment
- Dentin collar (ferrule) under 2-3mm
- When esthetics is especially important

**Recommended application**

The endodontically treated tooth should have good prognosis with no existing clinical symptoms.

Tooth preparation for endodontically treated teeth is generally done in three stages, namely:

Removal of the root canal filling material to the appropriate depth

Enlargement of the root canal post space

Preparation of the core for final restoration.

This can be achieved in a single step where the removal of gutta percha and the post space is created in a single step because of the unique qualities of the edelweiss POST DRILL. There is no need for a separate build-up of the core as the core is attached to the post.

1. Select the appropriate post and drill size by radiographic examination. Determine the required preparation depth.
2. Remove the coronal root filling material with Peeso Reamer drills or Gates Glidden drills (1000rpm-2000rpm) with a minimum of 5mm apical gutta percha seal intact or with the edelweiss POST DRILL.
3. Prepare the post space to the predetermined depth using the edelweiss POST DRILL (2000-5000rpm) with moderate pressure using water spray.
4. The corresponding size post is tried in the canal and should fit easily within the canal without binding.
5. The core portion of the POST & CORE can be adjusted and customized using the diamond drill and made to fit within the core space extra-coronally.
6. As the final step before cementing, clean or disinfect the root canal with a 2.5-5.25% sodium hypochlorite solution (NaOCl). Make sure no residues of temporary material or chemical disinfectant remains in the post space.
7. Rinse thoroughly with water and dry with paper tips.
8. Post preparation: Clean the surface of the POST & CORE with an alcohol swab and dry with oil-free air. Apply edelweiss VENEER Bond to the surface of the post and light cure for 20 seconds.
9. Post cementation:
- 9.1 Cementation using edelweiss Flowable composite: Etch the root canal walls with 35% phosphoric acid for 15 seconds, wash thoroughly with water and gently air-dry. Remove excess water from the post space using paper points. Apply primer-adhesive into the root canal with a thin micro brush/appliator tip. Remove excess primer adhesive with a paper point, gently air-dry and light-cure for 20 seconds or as per manufacturer's instructions. Insert the tip of the edelweiss Flowable syringe as deeply as possible into the root canal and starting apically, apply the cement, keeping the tip immersed in the cement; move the tip slowly upwards as the level of the paste rises. Keeping the opening of the tip immersed in the cement during the entire application process will minimize the entrapment of air bubbles/voids in the cement. Seat the post, twisting it slightly, and apply moderate pressure to hold it in position while removing excess cement. Light cure for 40 seconds.
- 9.2. Cementation with other types of Cements: edelweiss POST & CORE can be cemented with other types of luting cements like self-adhesive resin cements. Follow manufacturer's instructions for the respective cements.
10. Core adjustment: Finally, the core can be adjusted accordingly, to the final restoration.
11. Cleaning and Disinfection: Clean and autoclave drills before the first use and after use. Wash all blood or debris from drills under water. Clean drills using an ultrasonic cleaner in chlorhexidine gluconate solution. Remove the cleaning agent thoroughly using a brush under running water. Sterilize by autoclaving at 135°C for 10 minutes.

**Precautionary Information**

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known allergies to any of the ingredients. If irritation as a result of contact with oral soft tissue occurs, flush with large amounts of water and avoid usage. If allergic reaction occurs, seek medical attention as needed.

**Contraindications**

If known allergies to any of the components.

**Side effects**

Systemic side effects are not known. Individual cases of contact allergies have been reported for other products with a similar composition.

**Shelf-life and labeling**

The expiry date and the batch description are given on the container. Do not use after expiry date.

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## Important

To ensure optimal results, it is recommended to use in combination with all edelweiss components.

## Note

May only be supplied to dentists and dental laboratories or on their behalf. Keep away from children! Only for use by dentists!

Date of issue of this package information leaflet  
April, 2019

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**English**

**Definition**

edelweiss POST & CORE system are laser-sintered and pre-polymerized highly filled nanohybrid composite post and core monoblocs for use when restoring endodontically treated teeth. The post is a smooth surface tapered post that conforms to the internal post space design with a 1.4 mm and 1 mm diameter making it ideal for use in anterior and posterior teeth. The posts have a conical shape for perfect post space adaptation. The translucency of the post ensures uninterrupted light transmission for complete polymerization even in the apical region. The unique process of laser sintering allows for improved composite technology with superior mechanical properties.

The system consists of a prefabricated core attached to the post as a single unit. This unique and innovative monobloc design of the post and the core together in one unit provides complete stabilization within the root system. A core portion that corresponds to individual teeth in design is attached to the post. This is designed to simulate a crown preparation of the core material and thus saves the additional step of a core build-up. The attached core requires minimum preparation and can be adjusted for any type of restoration like veneer, inlay, onlay, crown or bridge. This makes core preparation easy and time saving. Reconstruction is quick and easy and provides an esthetic, functional and high quality outcome. The edelweiss POST & CORE ensures an highly esthetic result because of the inherent translucent post attached to an opaque core.

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**Features of edelweiss POST & CORE**

- Translucent tapered post allows for deep curing
- Opaque core allows for high esthetic restorations
- True monobloc effect enabling good retention
- Easy to use
- Cuts like dentin, can be easily adjusted in the mouth.
- Easy to customize
- Modulus of elasticity similar to natural tooth-preventing wedging effect
- Unmatched radiographic visibility: radio-opaque
- Superior adhesion to the root canal

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Tooth preparation for endodontically treated teeth is generally done in three stages, namely:  
Removal of the root canal filling material to the appropriate depth  
Enlargement of the root canal post space  
Preparation of the core for final restoration.

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1. Select the appropriate post and drill size by radiographic examination. Determine the required preparation depth.
2. Remove the coronal root filling material with Peeso Reamer drills or Gates Glidden drills (1000rpm-2000rpm) with a minimum of 5mm apical gutta percha seal intact or with the edelweiss POST DRILL.
3. Prepare the post space to the predetermined depth using the edelweiss POST DRILL (2000-5000rpm) with moderate pressure using water spray.
4. The corresponding size post is tried in the canal and should fit easily within the canal without binding.
5. The core portion of the POST & CORE can be adjusted and customized using the diamond drill and made to fit within the core space extra-coronally.
6. As the final step before cementing, clean or disinfect the root canal with a 2.5-5.25% sodium hypochlorite solution (NaOCl). Make sure no residues of temporary material or chemical disinfectant remains in the post space.
7. Rinse thoroughly with water and dry with paper tips.
8. Post preparation: Clean the surface of the POST & CORE with an alcohol swab and dry with oil-free air. Apply edelweiss VENEER Bond to the surface of the post and light cure for 20 seconds.
9. Post cementation:
- 9.1 Cementation using edelweiss Flowable composite: Etch the root canal walls with 35% phosphoric acid for 15 seconds, wash thoroughly with water and gently air-dry. Remove excess water from the post space using paper points. Apply primer-adhesive into the root canal with a thin micro brush/appligator tip. Remove excess primer adhesive with a paper point, gently air-dry and light-cure for 20 seconds or as per manufacturer's instructions. Insert the tip of the edelweiss Flowable syringe as deeply as possible into the root canal and starting apically, apply the cement, keeping the tip immersed in the cement; move the tip slowly upwards as the level of the paste rises. Keeping the opening of the tip immersed in the cement during the entire application process will minimize the entrapment of air bubbles/voids in the cement. Seat the post, twisting it slightly, and apply moderate pressure to hold it in position while removing excess cement. Light cure for 40 seconds.
- 9.2. Cementation with other types of Cements: edelweiss POST & CORE can be cemented with other types of luting cements like self-adhesive resin cements. Follow manufacturer's instructions for the respective cements.
10. Core adjustment: Finally, the core can be adjusted accordingly, to the final restoration.
11. Cleaning and Disinfection: Clean and autoclave drills before the first use and after use. Wash all blood or debris from drills under water. Clean drills using an ultrasonic cleaner in chlorhexidine gluconate solution. Remove the cleaning agent thoroughly using a brush under running water. Sterilize by autoclaving at 135°C for 10 minutes.

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