

AmpFi 2X qPCR MasterMix

DM-AmpqMM2500

Store at -20°C.

Description

AmpFi 2X qPCR MasterMix provides analysis of DNA samples for quantitative real-time in a convenient, reliable setup. This qPCR MasterMix contains Taq Polymerase and DNA Polymerase, providing rapid extension rates. This polymerase provides increased yields, and sensitivity, while optimizing reaction times by up to 70%, compared to wild-type Taq DNA polymerase.

The AmpFi 2X qPCR MasterMix possesses 5'-3' polymerase and 5'-3' exonuclease activities but does not have 3'-5' exonuclease activity. It generates amplicons with 3'-dA tails. qPCR products created with this MasterMix are compatible with TA cloning vectors. Additionally, the MasterMix includes a dye that is similar to both SYBR Green™ and EvaGreen™.

| Product Component | Quantity |
|-------------------------|------------------|
| AmpFi 2X qPCR MasterMix | 2500 rxn (25 ml) |
| ROX Reference Dye | 240 µl |

Protocol

The suggested quantity of ROX Reference Dye to be added to the MasterMix can differ based on the specific type of qPCR machine being used.

- No ROX equipment: Not needed.
- Low ROX equipment: 1 µl/1.25 ml or 22.5 µl/25 ml MasterMix.
- High ROX equipment: 11.5 µl/1.25 ml or 225 µl/25 ml MasterMix

1. Mix individual components before use and assemble reaction on ice.

| Component | Volume |
|--------------------------------|-------------------------------|
| AmpFi 2X qPCR MasterMix* | 10 µl |
| Forward Primer (10 µM) | 0.5 µl |
| Reverse Primer (10 µM) | 0.5 µl |
| Template DNA | Variable (100 ng genomic DNA) |
| Nuclease-free H ₂ O | Up to 20 µl |

*Reaction buffer contains 1.5 mM Mg²⁺

2. Gently mix the reaction components, and briefly centrifuge. Use thermocycling conditions:

| Step | Temperature | Duration (Standard) | Duration (Fast) | Cycles |
|---------------------|--|---------------------|-----------------|--------|
| Enzyme Activation | 95°C | 3 min | 3 min | 1 |
| Denaturation | 95°C | 15 sec | 1 sec | 40 |
| Annealing/Extension | 60°C | 1 min | 10 sec | |
| Melting Curve | Refer to specific guidelines for the instrument used | | | |

Notes

- The specialized buffer in this formulation enhances yields, sensitivity, and specificity beyond what is typically achieved with wild-type Taq polymerase.
- It's best to begin the qPCR immediately after preparing the reaction mixture. If that's not feasible, store the mixture on ice until you're ready to start the qPCR.
- For miRNA cDNA templates or other relevant applications, employ the standard thermocycling conditions.