

MACHEREY-NAGEL

User manual



- NucleoSpin® Inhibitor Removal








May 2021 / Rev. 02

DNA Clean up

Protocol at a glance (Rev. 02)

Method 5.1 for samples with moderate contaminations of humic substances or other types of inhibitors

NucleoSpin® Inhibitor Removal








1 Prepare sample		100 µL sample in a 1.5 mL or 2.0 mL tube
2 Adjust DNA binding conditions		300 µL IR1X Mix and incubate 1 min at RT 210 µL ethanol
3 Bind DNA		10,000 × g 30 sec
4 Wash silica membrane	 	500 µL IRW 10,000 × g 30 sec 300 µL IRW 10,000 × g 2 min
5 Prepare elution		100 µL BE RT, 1 min
6 Elute highly pure DNA		11,000 × g 1 min

DNA Clean up

Protocol at a glance (Rev. 02)

Method 5.2 for samples with considerable contaminations of humic substances
(intense brownish color)

NucleoSpin® Inhibitor Removal

1 Prepare sample		100 µL sample in a 1.5 mL tube
2 Adjust DNA binding conditions		500 µL IR1 mix
3 Bind DNA		10,000 × g 30 sec
4 Wash silica membrane	 	500 µL IRW 10,000 × g 30 sec 300 µL IRW 10,000 × g 2 min
5 Prepare elution		100 µL BE RT, 1 min
6 Elute highly pure DNA		11,000 × g 1 min

Contact MN

Germany and international

MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0
Toll-free: 0800 26 16 000 (Germany only)
E-mail: info@mn-net.com

Technical Support Bioanalysis

Tel.: +49 24 21 969-270
E-mail: tech-bio@mn-net.com

USA

MACHEREY-NAGEL Inc.
924 Marcon Blvd. · Suite 102 · Allentown PA, 18109 · USA
Toll-free: 888 321 6224 (MACH)
E-mail: sales-us@mn-net.com

France

MACHEREY-NAGEL SAS
1, rue Gutenberg – BP135 · 67720 Hoerdt Cedex · France
Tel.: +33 388 68 22 68
E-mail: sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifiée) au capital de 186600 €
Siret 379 859 531 00020 · RCS Strasbourg B379859531 · N° intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG
Hirsackerstr. 7 · 4702 Oensingen · Switzerland
Tel.: +41 62 388 55 00
E-mail: sales-ch@mn-net.com

Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Reagents, consumables, and equipment to be supplied by user	4
1.3	About this user manual	5
2	Product description	6
2.1	The basic principle	6
2.2	Kit specifications	6
2.3	Handling, preparation, and storage of starting materials	6
2.4	Elution procedures	7
3	Storage conditions and preparation of working solutions	8
4	Safety instructions	10
5	Protocols	11
5.1	Protocol for DNA clean up from samples with moderate contaminations of humic substances or other types of inhibitors	11
5.2	Protocol for DNA clean up from samples with considerable contamination of humic substances (intense brownish color)	13
6	Appendix	15
6.1	Troubleshooting	15
6.2	Ordering information	16
6.3	Product use restrictions / warranty	16

1 Components

1.1 Kit contents

NucleoSpin® Inhibitor Removal		
REF	740408.10	740408.50
NucleoSpin® Inhibitor Removal Column	10	50
Binding Buffer IR1*	10 mL	30 mL
Additive IRX*	1 mL	4 mL
Wash Buffer IRW* (Concentrate)	6 mL	12 mL
Elution Buffer BE**	13 mL	13 mL
Collection Tubes (2 mL)	10	50
Collection Tubes (1.5 mL)	10	50
User manual	1	1

1.2 Reagents, consumables, and equipment to be supplied by user

Reagents

- 96–100 % ethanol (for preparation of Wash Buffer IRW and adjustment of binding conditions for protocols 5.1)

Equipment

- Manual pipettors
- Centrifuge for microcentrifuge tubes
- Vortex mixer (e.g. Vortex-Genie 2 from Scientific Industries)
- Personal protection equipment (lab coat, gloves, goggles)

* For preparation of working solutions and storage conditions see section 3.

** Composition of Elution Buffer BE: 5 mM Tris/HCl, pH 8.5

1.3 About this user manual

It is strongly recommended reading the detailed protocol sections of this user manual if the NucleoSpin® Inhibitor Removal kit is used for the first time. Experienced users, however, may refer to the Protocol at a glance instead. The Protocol at a glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure. All technical literature is available on the internet at www.mn-net.com.

2 Product description

2.1 The basic principle

The NucleoSpin® Inhibitor Removal kit is designed for fast and efficient clean up of pre-purified DNA samples contaminated with PCR inhibitors. Samples might contain PCR inhibitors such as humic substances, heme, polyphenols, tannins, or dyes. Such inhibitors might originate from insufficient purification procedures or challenging sample materials such as soil, blood, fruit, brownish water, processed food (e.g. tea, coffee) or other.

Due to the diverse nature and amount of inhibitors that might be present in DNA preparations from problematic samples, the NucleoSpin® Inhibitor Removal kit provides two alternative procedures: method 5.1 is recommended for samples with slight contamination of humic acids or other PCR inhibitors. Method 5.2 is recommended for samples considerably contaminated with humic substances (strong brownish color).

2.2 Kit specifications

Kit specifications at a glance

Parameter	NucleoSpin® Inhibitor Removal
Technology	Silica membrane technology
Format	Mini spin column
Sample material	DNA solutions contaminated with PCR inhibitors
Sample amount	100 µL
DNA recovery	Typically >75 %
Elution volume	50–100 µL
Preparation time	15 min (6 preps)
Binding capacity	60 µg*

2.3 Handling, preparation, and storage of starting materials

DNA containing eluates should be kept on ice for short term storage and frozen at -20°C or below for long term storage.

*theoretical value

2.4 Elution procedures

In addition to the standard elution method, several modifications are possible to increase yield, concentration, and convenience.

- Convenient elution (standard elution): Elution can be performed by a single addition of 100 μ L Elution Buffer onto the column.
- High yield: Elution can be performed in two serial elutions of 100 μ L each, resulting in a total volume of 200 μ L.
- High concentration: Elution can be performed by application of 100 μ L Elution Buffer, which is then re-used in a second elution step, resulting in 100 μ L eluate with a high DNA concentration. Alternatively, the elution volume can be reduced down to 50 μ L. Please note that this typically will reduce the total amount of DNA recovered.

3 Storage conditions and preparation of working solutions

Attention: Buffer IR1 contains chaotropic salt. Wear gloves and goggles!

CAUTION: Buffer IR1 contains guanidine hydrochloride which can form highly reactive compounds when combined with bleach (sodium hypochlorite). Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

Before starting the NucleoSpin® Inhibitor Removal procedure, prepare the following:

Wash Buffer IRW: Add the indicated volume (see on the bottle or table below) of ethanol (96–100%) to IRW concentrate. Mark the label of the bottle to indicate that ethanol was added. Store Wash Buffer IRW at room temperature (18–25 °C) for up to one year.

NucleoSpin® Inhibitor Removal		
REF	10 preps 740408.10	50 preps 740408.50
Wash Buffer IRW	6 mL Add 24 mL ethanol	12 mL Add 48 mL ethanol

Preparation of buffer IR1X by supplementation of Buffer IR1 with additive IRX for use in procedure 5.1

Per preparation, combine 240 µL of buffer IR1 with 60 µL additive IRX. For convenience, a master mix sufficient for several preparations can be prepared.

Number of preps	Vol of IR1	Vol of IRX	Total Vol of IR1X
1	240 µL	60 µL	300 µL
2	480 µL	120 µL	600 µL
3	720 µL	180 µL	900 µL
4	960 µL	240 µL	1.2 mL
5	1.2 mL	300 µL	1.5 mL
10	2.4 mL	600 µL	3 mL
12	2.88 mL	720 µL	3.6 mL
20	4.8 mL	1.2 mL	6 mL

When preparing buffer IR1X for multiple preparations, it is recommended to prepare mix for one additional preparation to compensate for pipetting errors and attain sufficient mix for the planned number of preparations.

Attention: Additive IRX is viscous! Pipet slowly to avoid pipetting errors due to the viscosity of the additive.


Mix the solution by either pipetting up and down several times, incubation for several minutes on a rolling or inverting incubator or by moderate vortexing.

Note: The mixture contains detergent. Do not vortex strongly in order to avoid excessive foaming. Alternatively, you can prepare the master mix several hours in advance and wait until the foam dissolved. The mixture IR1X is stable for at least six months at 18–25 °C.

4 Safety instructions

The following components of the NucleoSpin® Inhibitor Removal kit contain hazardous contents. Wear gloves and goggles and follow the safety instructions given in this section.

Only harmful features do not need to be labeled with H and P phrases up to 125 mL or 125 g. *Mindergefährliche Eigenschaften müssen bis 125 mL oder 125 g nicht mit H- und P-Sätzen gekennzeichnet werden.*

Component	Hazard contents	GHS symbol	Hazard phrases	Precaution phrases
<i>Inhalt</i>	<i>Gefahrstoff</i>	<i>GHS-Symbol</i>	<i>H-Sätze</i>	<i>P-Sätze</i>
IR1	guanidine hydrochloride 50–66 % <i>Guanidinhydrochlorid 50–66 %</i> CAS 50-01-1	 WARNING ACHTUNG	302, 315, 319	264W, 280sh, 301+312, 330

Hazard phrases

- H 302 Harmful if swallowed.
Gesundheitsschädlich bei Verschlucken.
- H 315 Causes skin irritation.
Verursacht Hautreizungen.
- H 319 Causes serious eye irritation.
Verursacht schwere Augenreizung.

Precaution phrases

- P 264W Wash with water thoroughly after handling.
Nach Gebrauch mit Wasser gründlich waschen.
- P 280sh Wear protective gloves/eye protection.
Schutzhandschuhe/Augenschutz tragen.
- P 301+312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
BEI VERSCHLUCKEN: Bei Unwohlsein GIFTINFORMATIONSZENTRUM/Arzt anrufen.
- P 330 Rinse mouth.
Mund ausspülen.



The symbol shown on labels refers to further safety information in this section.
Das auf Etiketten dargestellte Symbol weist auf weitere Sicherheitsinformationen dieses Kapitels hin.

5 Protocols

5.1 Protocol for DNA clean up from samples with moderate contaminations of humic substances or other types of inhibitors

This protocol is especially recommended for samples containing inhibitors from e.g. blood (heme) or plant (polyphenols) or other and can also be used for samples with small amounts of inhibitors from water/soil (humic substances, slight brownish to yellow color).

Before starting the preparation:

- Check if buffer IR1X was prepared from IR1 and IRX according to section 3.
- Check if wash buffer IRW was prepared according to section 3.

Procedure

- 1 Starting material: Supply 100 μ L DNA solution (e.g. DNA solution containing moderate contamination of humic substances (slight brownish color)) or other types of inhibitors in a 1.5 mL or 2 mL tube (not provided).

Note: If your starting material consist of less than 100 μ L, fill it up to 100 μ L with Buffer BE.



**100 μ L sample
in a
1.5 mL
or
2.0 mL tube**

- 2 Add **300 μ L IR1X** and mix incubate for 1 min.

Note: Premix IR1 and IRX according to section 3 to obtain IR1X before starting the preparation.

Add **210 μ L ethanol** and mix.

Note: Do not premix IR1X and ethanol – a sequential addition of IR1X and ethanol is recommended.

Apply mixture onto **NucleoSpin® Inhibitor Removal Column** resting in a Collection Tube (2 mL, provided).



**300 μ L IR1X
Mix and incubate
1 min at RT
210 μ L ethanol**

- 3 Centrifuge for 30 s at 10,000 \times g.

Discard the flowthrough and reuse the collection tube.



**10,000 \times g
30 sec**

- 4** Add **500 µL Wash Buffer IRW** onto the column.
 Centrifuge for 30 s at 10,000 × g.
 Discard the flowthrough and reuse the collection tube.
 Add **300 µL Wash Buffer IRW** onto the column.
 Centrifuge for 2 min at 10,000 × g.
- Note: If flow through contaminates the column outlet upon removal of the assembly from the centrifuge or upon removal of the column from the collection tube, repeat the centrifugation step.*
- Discard the flow through with collection tube and put column into a fresh 1.5 mL collection tube (provided).
- 5** Add **100 µL Elution Buffer BE** onto the column.
Note: See section 2.4 for alternative elution procedures.
 Incubate for 1 min at room temperature.
- 6** Centrifuge for 1 min at 11,000 × g.
 Discard the column and use the purified eluate for further analysis.



500 µL IRW




**10,000 × g
30 sec**



**300 µL IRW
10,000 × g
2 min**



**100 µL BE
RT,
1 min**



**11,000 × g
1 min**





5.2 Protocol for DNA clean up from samples with considerable contamination of humic substances (intense brownish color)

This protocol is especially recommended for samples containing larger amounts of inhibitors originating from e.g. water / soil (humic substances) with moderate to strong brownish color.

Before starting the preparation:

- Check if wash buffer IRW was prepared according to section 3.

Procedure

- | | | |
|---|--|--|
| <p>1 Starting material: Supply 100 μL DNA solution (e.g. DNA solution containing humic substances (brownish color)) in a 1.5 mL or 2 mL tube (not provided).</p> <p><i>Note: If starting material is <100 μL, fill up to 100 μL with Buffer BE.</i></p> |  | <p>100 μL sample
in a
1.5 mL tube</p> |
| <p>2 Add 500 μL IR1 and mix.</p> <p>Apply mixture onto NucleoSpin® Inhibitor Removal Column resting in a Collection Tube (2 mL, provided).</p> |  | <p>500 μL IR1
mix</p> |
| <p>3 Centrifuge for 30 s at 10,000 \times g.</p> <p>Discard the flow through and reuse the collection tube.</p> |  | <p>10,000 \times g
30 sec</p> |
| <p>4 Add 500 μL Wash Buffer IRW onto the column.</p> <p>Centrifuge for 30 s at 10,000 \times g.</p> <p>Discard the flowthrough and reuse the collection tube.</p> <p>Add 300 μL Wash Buffer IRW onto the column.</p> <p>Centrifuge for 2 min at 10,000 \times g.</p> <p><i>Note: If flow through contaminates the column outlet upon removal of the assembly from the centrifuge or upon removal of the column from the collection tube, repeat the centrifugation step.</i></p> |  | <p>500 μL IRW</p> <p>10,000 \times g
30 sec</p> <p>300 μL IRW</p> <p>10,000 \times g
2 min</p> |
| <p>Discard the flow through with collection tube and put column into a fresh 1.5 mL collection tube (provided).</p> | | |

DNA Clean up

- 5** Add **100 μ L Elution Buffer BE** onto the column.

Note: See section 2.4 for alternative elution procedures.

Incubate for 1 min at room temperature.



100 μ L BE

**RT,
1 min**

- 6** Centrifuge for 1 min at 11,000 \times g.

Discard the column and use the purified eluate for further analysis.



**11,000 \times g
1 min**

6 Appendix

6.1 Troubleshooting

Problem	Possible cause and suggestions
No or low DNA recovery	<p data-bbox="468 308 829 330"><i>DNA washed off the silica membrane:</i></p> <ul data-bbox="468 347 967 395" style="list-style-type: none"> • make sure to prepare the Wash Buffer by adding the appropriate amount of ethanol. <p data-bbox="468 419 757 442"><i>Low initial DNA concentration:</i></p> <ul data-bbox="468 459 967 507" style="list-style-type: none"> • Provide higher initial DNA concentration or uses more sensitive analysis methods
	<p data-bbox="468 531 972 553"><i>Falsely quantification of DNA in the provided sample:</i></p> <ul data-bbox="468 571 978 802" style="list-style-type: none"> • Presence of PCR inhibitory substances (e.g. humic substances, polyphenols) in a sample will substantially influence DNA quantification by spectrophotometry and fluorescent methods causing a considerable over- or under-estimation of DNA, especially in colored samples or samples with unacceptable $A_{260/280}$ or $A_{260/230}$ ratios. Do not trust DNA quantification results of impure DNA solutions.
	<p data-bbox="468 826 792 849"><i>Insufficient PCR inhibitor removal:</i></p> <ul data-bbox="468 866 978 1042" style="list-style-type: none"> • The kit is designed to remove diverse PCR inhibitors like e.g. humic substances (brownish color) from DNA solutions. Because the chemical nature of PCR inhibitory substances is diverse, some inhibitors might not be effectively removed. If procedure 5.1 does not give satisfactory results, try procedure 5.2.
PCR inhibition	<p data-bbox="468 1066 631 1088"><i>Brownish eluate:</i></p> <ul data-bbox="468 1106 967 1177" style="list-style-type: none"> • If an insufficient decoloration of the sample is observed using procedure 5.1, try procedure 5.2 for DNA clean up.
No or insufficient decoloration of the sample	<p data-bbox="468 1201 575 1224"><i>Low purity:</i></p> <ul data-bbox="468 1241 978 1369" style="list-style-type: none"> • Quality ratio determination strongly depends on a sufficient amount of DNA measured. Make sure to use a sufficient amount of DNA that has been validated to enable a meaningful ratio determination with the photometric systems used.
Low $A_{260/280}$ or $A_{260/230}$ ratio	As above.
No improvement of $A_{260/280}$ or $A_{260/230}$ ratio	As above.

6.2 Ordering information

Product	REF	Pack of
NucleoSpin® Inhibitor Removal	740408.10/50	10/50
NucleoSpin® eDNA water	740402.10/50	10/50
Collection Tubes (2 mL)	740600	1000
Buffer BE (125 mL)	740306.100	1

6.3 Product use restrictions/warranty

NucleoSpin® Inhibitor Removal kit components are intended, developed, designed, and sold FOR RESEARCH PURPOSES ONLY, except, however, any other function of the product being expressly described in original MACHEREY-NAGEL product leaflets. MACHEREY-NAGEL products are intended for GENERAL LABORATORY USE ONLY! MACHEREY-NAGEL products are suited for QUALIFIED PERSONNEL ONLY! MACHEREY-NAGEL products shall in any event only be used wearing adequate PROTECTIVE CLOTHING. For detailed information please refer to the respective Material Safety Data Sheet of the product! MACHEREY-NAGEL products shall exclusively be used in an ADEQUATE TEST ENVIRONMENT. MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products in other fields of application. Application on the human body is STRICTLY FORBIDDEN. The respective user is liable for any and all damages resulting from such application. DNA/RNA/PROTEIN purification products of MACHEREY-NAGEL are suitable for IN VITRO-USES ONLY! ONLY MACHEREY-NAGEL products specially labeled as IVD are also suitable for IN VITRO-diagnostic use. Please pay attention to the package of the product. IN VITRO diagnostic products are expressly marked as IVD on the packaging. IF THERE IS NO IVD SIGN, THE PRODUCT SHALL NOT BE SUITABLE FOR IN VITRO-DIAGNOSTIC USE! ALL OTHER PRODUCTS NOT LABELED AS IVD ARE NOT SUITED FOR ANY CLINICAL USE (INCLUDING, BUT NOT LIMITED TO DIAGNOSTIC, THERAPEUTIC AND/OR PROGNOSTIC USE). No claim or representations is intended for its use to identify any specific organism or for clinical use (included, but not limited to diagnostic, prognostic, therapeutic, or blood banking). It is rather in the responsibility of the user or - in any case of resale of the products - in the responsibility of the reseller to inspect and assure the use of the DNA/RNA/protein purification products of MACHEREY-NAGEL for a well-defined and specific application. MACHEREY-NAGEL shall only be responsible for the product specifications and the performance range of MN products according to the specifications of in-house quality control, product documentation and marketing material. This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish to get an extra copy. There is no warranty for and MACHEREY-NAGEL is not liable for damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; defects in products or components not manufactured by MACHEREY-NAGEL, or damages resulting from such non-MACHEREY-NAGEL components or products. MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES


ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, REPRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS. In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications. This warranty is exclusive and MACHEREY-NAGEL makes no other warranty expressed or implied. The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty. Product claims are subject to change. Therefore, please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications. Last updated: 02/2020, Rev. 01

Please contact: MACHEREY-NAGEL GmbH & Co. KG
Tel.: +49 24 21 969-270
e-mail: tech-bio@mn-net.com

Trademarks:

NucleoSpin® is a registered trademark of MACHEREY-NAGEL GmbH & Co. KG

All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and cannot be seen as a kind of recommendation or assessment). Regarding these products or services we cannot grant any guarantees regarding selection, efficiency, or operation.



www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG
Valenciennener Str. 11
52355 Düren · Germany

DE	Tel.: +49 24 21 969-0	info@mn-net.com
CH	Tel.: +41 62 388 55 00	sales-ch@mn-net.com
FR	Tel.: +33 388 68 22 68	sales-fr@mn-net.com
US	Tel.: +1 888 321 62 24	sales-us@mn-net.com



A068515/0510.15