

NucleoType Plant PCR

Storage conditions and preparations of working solutions

The NucleoType Plant PCR kit should be stored upon arrival at +4 °C or -20 °C (recommended). The kit is stable for at least 12 months when stored at this temperature. The kit can be shipped at ambient temperature (18–25 °C) for up to 3 months. Short time exposure (up to 14 days) at temperatures up to 37 °C is tolerable.

Store all kit components at +4 °C or -20 °C (recommended) upon arrival and after first time usage. Store NucleoType HotStart PCR Master Mix in the dark, e.g. within the product box in a freezer (-20 °C; recommended) or fridge (+4 °C). Avoid prolonged exposure of the mix to light. Setting up PCR at average laboratory illumination is tolerable. Do not expose the mix to direct sunlight.

NucleoType HotStart PCR Master Mix (2x) is ready to use.

Prepare a primer mix according to the recommended concentration per primer as described.

Kit contents

NucleoType Plant PCR			
REF	25 preps 743202.25	100 preps 743202.100	500 preps 743202.500
Plant Transfer Tool (PTT)	25 pieces	100 pieces	500 pieces
NucleoType HotStart PCR Master Mix (2x)	125 µL	500 µL	2 x 1250 µL
User manual	1	1	1

Reagents, consumables, and equipment to be supplied by user

Reagents

- Primer for plant specific target of interest
- Water (PCR grade; for primer dilution and reaction fill-up)

Consumables

- Disposable pipette tips
- PCR tubes

Equipment

- Manual pipettes
- Personal protection equipment (lab coat, gloves, goggles)
- PCR cyclers
- Gel electrophoretic equipment or Bioanalyzer® for analysis of generated amplicons

About this user manual

It is strongly recommended for first time users to read the detailed protocol sections of the NucleoType Plant PCR kit before using these products.

All technical literature is available online at www.mn-net.com.

Please contact technical service regarding information about any changes to the current user manual compared with previous revisions.

Product description

The basic principle

Many plant genotyping methods are based on DNA purification from plant material, followed by PCR amplification of genes of interest. However, DNA purification from plant material is a time consuming and elaborate process.

The NucleoType Plant PCR Kit is designed for rapid plant typing experiments using plant leaf material as sample, without the need to purify DNA.

For plant leaf material, the NucleoType Plant PCR kit provides the patent pending Plant Transfer Tool (PTT) for sample take up, inhibitor inactivation, and transfer of the sample into the PCR. The sample harvest (pierce PTT into leaf) and transfer of plant material into the prepared PCR mix (short contact of PTT with PCR mix) takes less than five seconds. Due to the patent pending coating of the transfer tool, PCR inhibitors present in many plant materials are inactivated to a considerably extend during sample transfer.

Kit specifications

Kit specifications at a glance

Parameter	NucleoType Plant PCR
Technology	Direct PCR: Transfer of plant leaf aliquot with Plant Transfer Tool (PTT) directly into PCR Mix
Format	10 µL HotStart PCR (optional up to 50 µL)
Sample type	Plant leaf material from e.g., corn, soybean, wheat, <i>arabidopsis</i> , tobacco, cotton, grape wine, cress, as well as from kiwi, banana, and avocado fruit flesh
Preparation time	Sample preparation: < 1 min; PCR cycling: 30–90 min (cycler and target size dependent)
Amplicon size	Up to 1000 bp
Analysis	Gel electrophoresis: Approx. 30 min (40 samples); Bioanalyzer®: Approx. 40 min (12 samples)

Handling, preparation, and storage of starting material

The kit is designed to perform genotyping from fresh plant material like leaves. Plant material stored at 0–8 °C for several days or frozen material may also be used. However, kit performance with plant material stored for a long time might differ to the performance of fresh material.

Lysis, disruption, and transfer of sample material

In order to obtain reliable plant typing data, it is important to obtain a sufficient amount of DNA in a form suitable to serve as template for subsequent PCR amplification.

The NucleoType Plant PCR kit provides the Plant Transfer Tool (PTT) enabling an easy and patent pending sample uptake from plant leaves and transfer into the PCR mix. A special sample lysis step is not required. Upon insertion of the PTT into the plant leaf material, plant cells are disrupted and sufficient DNA will adhere to the PTT. During transfer of the sample into the PCR, inhibitors are inactivated due to the special coating of the PTT considerably. After sample transfer, it is possible to store the PCR mx for up to 2 hours at 4 °C–37 °C, due to HotStart function of the NucleoType HotStart PCR Master Mix.

Safety instructions

Use the product according to the user manual.

The product does not contain components requiring GHS hazard or precaution phrases.

PCR cycling parameters

Cycling conditions are depending on primer, target length, and PCR cycler setup. For several primer pairs with T_m ranging from 40–75 °C the following PCR programs have been used successfully.

PCR program 1 (three step program for typical endpoint PCR cycler)

Initial denaturation	95 °C	2 min	1 cycle
Amplification:	95 °C	15 s	40 cycles
	40–75 °C***	20 s**	
	72 °C	60 s**	
Extension	72 °C	1 min	1 cycle
Cooling	4 °C		
Total time		Approx. 70–100 min (total run time is annealing temperature and machine dependent)	

PCR program 2 (two step program for typical end point PCR machines)

Initial denaturation	95 °C	2 min	1 cycle
Amplification	95 °C	15 s	40 cycles
	60–72 °C*	60 s**	
	72 °C	1 min	
Extension	72 °C	1 min	1 cycle
Cooling	4 °C		
Total time		Approx. 60–70 min (total run time is annealing / elongation temperature and machine dependent)	

Note: For amplification of fragments smaller 1000 bp and/or for amplifications with a PCR machine with slow ramp rates (2 °C/s), the annealing time and extension time may be reduced stepwise, e.g. down to 15 seconds for annealing and 15 seconds for extension.

PCR program 3 (e.g., LightCycler® 1.5 in glass capillary)

Initial denaturation	95 °C	2 min	1 cycle
Amplification	95 °C	15 s	40 cycles
	40–75 °C***	15 s	
	72 °C	30 s	
Extension	72 °C	1 min	1 cycle
Cooling	20 °C		
Total time		Approx. 30–60 min (total run time is annealing / elongation temperature and machine dependent)	

Note: The LightCycler® is used herein solely as a fast cycling instrument, but not for quantitative PCR!

Analysis of PCR products

The PCR products (amplicons) can directly be analyzed using the following methods.

There is no need to add loading dye for gel electrophoresis because the PCR mix already contains a dye and suitable density.

There is no need to perform a proteinase digestion step prior to analysis of the amplicons.

- Gel electrophoresis: Apply the total PCR reaction onto an e.g., 1–2% agarose gel for analysis.
- Dye migration in
1% agarose gel: Approximately as 600 bp fragment
2% agarose gel: Approximately as 350 bp fragment
- Bioanalyzer® (Agilent): Use 1 µL with e.g. the Agilent DNA 1000 Kit.

* The optimal annealing/extension time is primer dependent. Only primers with melting temperature above 60 °C are recommended for this program.

** For initial testing an annealing/extension time of 20 and/or 60 seconds seconds is recommended!

*** Optimal annealing temperature is primer dependent and has to be determined empirically. A good starting point for testing is 50 °C. Optimally, a good annealing temperature for primer of your choice is determined with a temperature gradient cycler.

Protocols

Plant typing with plant leaf material using the Plant Transfer Tools (PTT)

A: *In situ* sampling (on-site sampling, greenhouse sampling)

Sampling may be executed *in situ*, i.e., the plant to be typed or parts thereof are not transferred to the lab.

1 Prepare Plant Transfer Tool

At the plant growing site, withdraw a Plant Transfer Tool (PTT) from the bag, grabbing it firmly on the brushed end.

Note: Do not touch the tip of the PTT with your finger or anything else except the sample in order to prevent any contamination and to prevent the active ingredient containing coating of the PTT tip to be wiped off.

2 Take up plant sample

Pierce the PTT into, or slightly through the plant leaf. Optimally, pierce the PTT into the leaf while the leaf is backed by a solid to semi solid support. Do not pierce the PTT more than approximately 1–3 mm into/through the leaf. See “Troubleshooting” section for additional recommendations.

3 Transfer sample into the PCR




Remove the PTT from the plant and dip the tip of the PTT briefly into the PCR tube (one touch, approximately 1 second, no stirring) containing the PCR mix. Alternatively, it is possible to transport the PTT with adhered plant material (which is often almost invisible) to the lab, followed by plant sample transfer into the PCR tube (one touch, approximately 1 second, no stirring) containing the PCR mix.

Note: If transporting the PTT with adhered plant material, it is important that the PTT tip carrying plant material does not touch anything in order to avoid removal of surface coating. The PTT can be placed tip up in, e.g., a tube rack.

4 Store sample or start PCR

Discard the PTT and close the lid of the PCR tube.

Note: After sample transfer, it is possible to store the PCR Mix for up to 2 hours in the dark at 4 °C–37 °C, due to HotStart function of the NucleoType HotStart PCR Master Mix. Avoid prolonged exposure of the PCR to light – do not expose to direct sunlight!



B: *Ex situ* sampling (lab sampling, sampling from cut off leaf material)

Plants, leaves, or parts of leaves can be transported to the lab for sampling.

1 Prepare Plant Transfer Tool

Withdraw a Plant Transfer Tool (PTT) from the bag grabbing it firmly on the brushed end.

Note: Do not touch the tip of the PTT with your finger or anything else except the sample in order to prevent any contamination and to prevent the active ingredient containing coating of the PTT tip to be wiped off.

2 Take up plant sample

Pierce the PTT into, or slightly through the plant leaf (see section A2).

3 Transfer sample into the PCR

Dip the tip of the PTT briefly into the PCR (one touch, approximately one second, no stirring), which was prepared in advance. Close the lid of the PCR tube and start PCR cycling.

Reaction setup for 10 µL PCR (single-plex or duplex)

A final PCR with 10 µL volume is the recommended standard reaction volume for the **NucleoType Plant PCR** kits. Due to the hot start technology of this product, the reaction setup can be performed at room temperature (18–25 °C).

Per reaction combine the following

- 5 µL NucleoType HotStart PCR Master Mix (2x).
 - 5 µL primer mix (each target primer with a concentration of 0.4 µM within the 5 µL primer mix, resulting in a final concentration of 0.2 µM in the PCR per primer).
- 10 µL final PCR volume, ready to receive the plant sample directly from the Plant Transfer Tool (PTT) or 1 µL from the seed lysate.

Note: The addition of silicone oil is not necessary and will impair removal of the liquid after the reaction. Therefore, the addition of silicone oil is not recommended.

Note: If desired, the final PCR volume can be scaled up by increasing all components proportionally.

Troubleshooting

Reduction of PCR volume

- *Reduction of initial 10 µL PCR volume during cycling*
Depending on the PCR tube size and quality, the initial 10 µL setup volume might shrink to approximately 8 µL during cycling. This is acceptable and does not impair typing performance. If volume reduction is more pronounced, use a smaller and/or tighter reaction tube.

To little or no amplicon detected

- *Unfavorable primer selection*
Make sure that the primers are selected well and are able to amplify the desired target from 1–10 ng of purified genomic DNA. Test different primer annealing temperatures.
- *Unfavorable storage conditions*
Avoid prolonged exposure of the mix to light – see section Storage Conditions.
- *Too much sample material in PCR*
Make sure to transfer plant leaf material with the Plant Transfer Tool (PTT) into the PCR. Typically, plant leaf material adhering to the PTT is faintly visible. Do not transfer bits and pieces of plant material.
- *No Plant Transfer Tool (PTT) used*
The substitution of the PTT by a common toothpick may cause reaction failure. The PTT is coated with an active component which inactivates several PCR inhibitors which are commonly present in plant leaves.
- *PTT not handled with care*
The PTT tip is coated with an active ingredient which might be wiped off in

case of too much contact of the tip with surfaces or fingers. Do not touch the PTT tip – grab it firmly on the brushed side!

- *PTT uses as a stirrer*
Do not “stir in” the sample into the PCR – a short, approximately 1 second dip in of the PTT tip into the PCR mix without stirring is sufficient for sample transfer!
- *Unfavorable PCR program*
Try to adjust annealing temperature and time as well as extension time. Note that PCR machines with rapid ramp rates require longer annealing and extension times than PCR machines with slow ramp rates, because there is less amplification time for the polymerase during ramping!
- *Insufficient uptake of sample material*
For very small leaves, its possible to circumvent insufficient sample uptake by using a support pad. The leaf should be backed up by a solid to semi solid support (e.g. rubber pad, piece of disposable plastic sheet) and the PTT is pierced into the leaf.
- For common leaves like tobacco or wine leaf, the PTT is inserted less than approximately 1 mm. Do not pierce the PTT in or through the plant leaf for more than 1–3 mm.

Amplicon does not have the correct size

- *Primer selection*
Make sure that the primers are selected well and are able to amplify the desired target from 1–10 ng or purified genomic DNA.

Amplicon number is not correct

- *Sensitivity of analysis method*
Make sure that the analysis method has enough resolving power to discriminate the two different sizes of DNA fragments.
- Use Bioanalyzer® instead of gel electrophoreses or increase electrophoresis time or gel concentration.
- Make sure that both primer pairs have a similar amplification efficiency. If this is not the case, titrate down the primer pair yielding an amplicon (use a smaller concentration for this primer pair).

Product use and restriction / warranty

NucleoType Plant PCR kit components are intended, developed, designed, and sold FOR RESEARCH PURPOSES ONLY. They are suitable for *in vitro* uses only. No claim or representation is intended for its use to identify any specific organism or for clinical use. MACHEREY-NAGEL does not warrant against damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; against defects in products or components not manufactured by MACHEREY-NAGEL, or against damages resulting from such non-MACHEREY-NAGEL components or products.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications.

For more detailed product use restriction/warranty please have a look at: www.mn-net.com

Please contact: MACHEREY-NAGEL Germany
Tel.: +49 (0) 24 21 969-270 • e-mail: TECH-BIO@mn-net.com

Trademarks / disclaimer:

Bioanalyzer® is a registered trademark of Agilent, Inc.
Eppendorf® is a registered trademark of Eppendorf AG.

LightCycle® is a trademark of ROCHE DIAGNOSTICS GMBH

All used names and denotations can be brands, trademarks or registered labels of their respective owner, even if they are not specially denoted. To mention products and brands is only a kind of information, i.e. it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment. Regarding these products or services we can not grant any guarantees regarding selection, efficiency or operation.