Precipitate small RNA

Add **500 µL** of **100 % isopropanol** (~0.8 vol) to the supernatant of the previous step. Incubate at **4 °C** for **30 min**.

Centrifuge for **15 min** at **12,000 x g**. Discard the supernatant.

Wash large RNA

Add **500 µL 75 % ethanol** to the pellet of the large RNA preparation. Centrifuge for **3 min** at **8,000 x g** and discard supernatant. Repeat washing. **Drying of the pellet is not necessary.**

Wash small RNA

Add **500 µL 70 % isopropanol** to the pellet of small RNA preparation. Centrifuge for **3 min** at **8,000 x g** and discard supernatant. Repeat washing. **Drying of the pellet is not necessary.**

Reconstitute RNA

Dissolve the RNA pellet in RNase-free water to obtain an RNA concentration of **1 µg/µL** for the large RNA fraction and about **0.1 µg/µL** for the small RNA fraction. Vortex the sample for **3 min** at **RT**.

6 Product use restrictions / warranty

NucleoZOL is developed, designed, and sold for research purposes only. We offer a one year warranty from date of delivery that our products will conform to applicable specifications set forth in the product specifications if not sold to persons set forth in § 13 of German Civil Code. In such a case, the provisions of the German Civil Code shall be valid. This warranty does not cover defects in or damage to the products which are due to improper installation or maintenance, misuse, neglect or any use other than ordinary commercial application. The seller shall not be liable for the goods being fit for a particular purpose unless otherwise agreed upon, to which the buyer intends to put them. This warranty is strictly exclusive. MACHEREY-NAGEL MAKES no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER; DIRECTLY OR INDIRECTLY; EXPRESS OR IMPLIED; INCLUDING; WITHOUT LIMITATION; AS TO THE SUITABILITY; REPRODUCIBILITY; DURABILITY; FITNESS FOR A PARTICULAR PURPOSE OR USE; MERCHANTABILITY; CONDITION; OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS OR THE SALE OF MACHEREY-NAGEL PRODUCTS. If not superseded by these terms the current terms and conditions for the sale of goods (export version) are additionally valid.

Patents:

US Patents 7,794,932 and 8,367,817 B2 and international patent

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany
Tel.: +49 24 21 969-0 · tech-bio@mn-net.com · www.mn-net.com

NucleoZOL

1 Contents

NucleoZOL		
REF	740404.6	740404.200
NucleoZOL reagent	6 mL	200 mL

2 Product description

NucleoZOL is a convenient reagent based on a combination of phenol with guanidinium thiocyanate which enables a quick isolation of total and small RNA from biological samples without liquid phase separation. NucleoZOL isolates pure and undegraded total RNA (small and large RNA) that is ready for downstream applications.

3 Storage conditions and preparation of working solutions

NucleoZOL can be stored at 15–25 °C and is stable for at least one year.

4 Safety instructions

NucleoZOL contains hazardous contents. Wear gloves and goggles!

GHS classification

Only harmful features must not be labeled with H and P phrases up to 125 mL or 125 g.

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5 Protocol

Please read the detailed protocol if using NucleoZOL for the first time. Download the detailed NucleoZOL user manual from the corresponding webpage www.mn-net.com/NucleoZOL. Experienced users may refer to this short instruction manual.

Before starting the preparation:

Reagents, consumables, and equipment to be supplied by the user:

- Reagents: RNase-free water, 75 % ethanol, 70 % isopropanol, 100 % isopropanol.
- Consumables: 1.5 mL, 2.0 mL, or 15 mL centrifuge tubes (depending on amount of sample to be processed per preparation), sterile RNase-free tips.
- Equipment: Manual pipettors, vortex mixer, centrifuge for microcentrifuge tubes, equipment for sample disruption and homogenization, personal protection equipment (e.g., lab coat, gloves, goggles), well ventilated working environment, RNase-free working environment.

5.1 Isolation of total RNA (small and large RNA in one fraction)

This protocol is used for the isolation of total RNA, containing small RNA (< 200 nt) and large RNA (> 200 nt) in one fraction. This protocol describes the isolation of RNA with 500 µL NucleoZOL. The procedure can be scaled up or down, dependent on the sample input.

Lyse sample

Tissue

Homogenize **50 mg tissue/500 µL NucleoZOL** with e.g., a rotor-stator-homogenizer.

Cells

Homogenize up to **5x10⁶ cells/500 µL NucleoZOL** by pipetting.

Liquids

Homogenize **200 µL liquid sample/500 µL NucleoZOL** by pipetting.

Precipitate contaminants

Add **200 µL water/500 µL NucleoZOL** to the lysate.

Mix vigorously for **15 s** and incubate **5 min** at **RT**.

Centrifuge for **15 min** at **12,000 x g**.

Precipitate total RNA

Transfer 500 µL of the supernatant to a new tube. Add **500 µL 100 % isopropanol**.

Incubate for **10 min** at **RT**. Centrifuge for **10 min** at **12,000 x g**.

Discard the supernatant.

Wash total RNA

Add **500 µL 75 % ethanol** to the pellet.

Centrifuge for **3 min** at **8,000 x g** and discard supernatant. Repeat washing. **Drying of the pellet is not necessary.**

Reconstitute RNA

Add an appropriate volume of water to reconstitute the total RNA to approx. 1 µg/µL. **Dissolve** RNA pellets by vortexing at room temperature for **3 min**.

5.2 Isolation of small and large RNA in two separate fractions

This protocol is used for the isolation of small RNA (RNA < 200 nt) and large RNA (RNA > 200 nt) in two separate fractions. This protocol describes the isolation of RNA with 500 µL NucleoZOL. The procedure can be scaled up or down, dependent on the sample input.

Lyse sample

Tissue

Homogenize **50 mg tissue/500 µL NucleoZOL** with e.g., a rotor-stator-homogenizer.

Cells

Homogenize up to **5x10⁶ cells/µL NucleoZOL** by pipetting.

Liquids

Homogenize **200 µL liquid sample/500 µL NucleoZOL** by pipetting.

Precipitate contaminants

Add **200 µL water/500 µL NucleoZOL** to the lysate.

Mix vigorously for **15 s** and incubate **5 min** at **RT**.

Centrifuge for **15 min** at **12,000 x g**.

Precipitate large RNA

Transfer 500 µL supernatant to a new tube.

Add **200 µL 75 % ethanol**. Incubate **10 min** at **RT**. Centrifuge for **8 min** at **12,000 x g**.

Continue with the pellet for large RNA isolation. Transfer supernatant to a new tube for miRNA isolation.