

# NucleoBond® RS purification workflow

## Protocol at a glance (Rev. 01)

	RS 10	RS 50	RS 100	RS 200
<b>1–3</b> <b>Cell cultivation and harvest</b>	4 °C, 15 min 4,500–6,000 × g			
	max. recom: ODV: 3,000 CDW: 8 g typ. Yield: ~10 mg	max. recom: ODV: 15,000 CDW: 40 g typ. Yield: ~50 mg	max. recom: ODV: 30,000 CDW: 80 g typ. Yield: ~100 mg	max. recom: ODV: 54,000 CDW: 145 g typ. Yield: ~190 mg
<b>4–5</b> <b>Cell lysis</b>  <i>(Important: Check Buffer LYS-EF for precipitated SDS)</i>	60 mL Buffer RES-EF	300 mL Buffer RES-EF	600 mL Buffer RES-EF	1080 mL Buffer RES-EF
	60 mL Buffer LYS-EF	300 mL Buffer LYS-EF	600 mL Buffer LYS-EF	1080 mL Buffer LYS-EF
RT, max 5 min				
<b>6 Neutralization</b>	60 mL Buffer NEU-EF	300 mL Buffer NEU-EF	600 mL Buffer NEU-EF	1080 mL Buffer NEU-EF
	Mix thoroughly until colorless Incubate 15 min on ice!			
<i>Maximum recommended pump speed</i>	5 mL/min	10 mL/min	10 mL/min	10 mL/min
<b>7 Equilibration of the column</b>	30 mL EQ-EF	120 mL EQ-EF	250 mL EQ-EF	400 mL EQ-EF
<b>8 Clarification of the lysate</b>	Clarify and filtrate the crude lysate to 0.45/0.22 µm			
<b>9 Binding</b>	Load cleared lysate on NucleoBond® RS column			
<b>10 1<sup>st</sup> Wash</b>	30 mL Buffer ENTO-EF	150 mL Buffer ENTO-EF	300 mL Buffer ENTO-EF	550 mL Buffer ENTO-EF
<b>11 2<sup>nd</sup> Wash</b>	30 mL Buffer WASH-EF	150 mL Buffer WASH-EF	300 mL Buffer WASH-EF	550 mL Buffer WASH-EF
<b>13 Elution</b>	30 mL Buffer ELU-EF	120 mL Buffer ELU-EF	250 mL Buffer ELU-EF	400 mL Buffer ELU-EF
<b>14 Precipitation</b>	21 mL Isopropanol	84 mL Isopropanol	175 mL Isopropanol	280 mL Isopropanol
	Mix thoroughly 4,500–15,000 × g 4 °C, 15 min			
<b>15 Washing and drying</b>	5 mL 70 % EtOH 4,500–15,000 × g RT, 5 min  Approx. RT, 10–15 min	15 mL 70 % EtOH 4,500–15,000 × g RT, 5 min  Approx. RT, 30–60 min	25 mL 70 % EtOH 4,500–15,000 × g RT, 5 min  Approx. RT, 45–90 min	50 mL 70 % EtOH 4,500–15,000 × g RT, 5 min  Approx. RT, 75–120 min
	Appropriate volume of TE-EF or H <sub>2</sub> O-EF.  It is recommended to choose the resuspension volume according to the requirements of the downstream application. (for var., @ 1 mg/mL)			
<b>16 Reconstitution</b>	8–10 mL	40–50 mL	80–100 mL	170–190 mL

# NucleoBond® RS purification workflow

## Protocol at a glance (Rev. 01)

