

## Exonuclease I, *E. coli*

Cat. Nos. X40501K, X40505K, and X40520K



Connect with Epicentre on our blog ([epicentral.blogspot.com](http://epicentral.blogspot.com)),  
Facebook ([facebook.com/EpicentreBio](https://www.facebook.com/EpicentreBio)), and Twitter ([@EpicentreBio](https://twitter.com/EpicentreBio)).

## 1. Introduction

Exonuclease I (Exo I) the product of the *sbcB* gene of *E. coli*, is an exodeoxyribonuclease that hydrolyzes single-stranded DNA (ssDNA) stepwise in a 3'→5' direction.<sup>1-3</sup> Hydrolysis generates deoxyribonucleoside 5' monophosphates and a terminal dinucleotide diphosphate.<sup>1</sup> The enzyme requires magnesium (optimal Mg<sup>2+</sup> concentration is 10 mM) and the presence of a free 3'-hydroxyl terminus.<sup>1</sup> Exo I is active under a wide variety of buffer conditions, allowing addition of the enzyme directly into most reaction mixes. Heat inactivation results from incubation at 80°C for 15 minutes.

Exo I is available in 1,000-, 5,000-, and 20,000-Unit sizes at a concentration of 20 Units/μl.

## 2. Applications

**Removal of residual ssDNA and oligonucleotides from reaction mixes.** Linear ssDNA and oligonucleotides can be selectively degraded from heterogeneous mixtures of nucleic acids in reaction mixes.

**Removal of ssDNA from nucleic acid mixtures.** Linear ssDNA can be selectively degraded from heterogeneous mixtures of nucleic acids with Exo I.

**Assay for regions of ssDNA.**<sup>4,5</sup> Use Exo I to assay for the presence of ssDNA containing a free 3'-hydroxyl end. This technique was used to characterize the endonuclease and helicase activities of purified *recBC* protein on circular fd phage DNA and duplex phage T7 DNA respectively.

## 3. Product Specifications

**Storage:** Store only at -20°C in a freezer without a defrost cycle.

**Storage Buffer:** Exo I is supplied in a 50% glycerol solution containing 50 mM Tris-HCl (pH 7.5), 100 mM NaCl, 1.0 mM dithiothreitol, 0.1 mM EDTA, and 0.1% Triton® X-100.

**Unit Definition:** One unit of Exo I results in the acid-solubilization of 10 nmol of nucleotides from calf thymus DNA in 30 minutes at 37°C.

**Quality Control:** Exo I is function-tested in a reaction containing 33 mM Tris-acetate (pH 7.5), 66 mM potassium acetate, 10 mM magnesium acetate, 0.5 mM dithiothreitol, 10 μg of denatured calf thymus DNA, and varying amounts of Exo I.

**Contaminating Activity Assays:** Exo I is free of detectable RNase, endonuclease, and double-stranded exonuclease activities.

## 4. Related Products

The following products are also available:

- Exonuclease III
- Lambda Exonuclease
- RNase-Free DNase I
- Mung Bean Nuclease
- OmniCleave™ Endonuclease
- Plasmid-Safe™ ATP-Dependent DNase

## 5. References

1. Lehman, I.R. and Nussbaum, A.L. (1964) *J. Biol. Chem.* **239**, 2628.
2. Kusher, S.R. *et al.*, (1971) *Proc. Natl. Acad. Sci. USA* **68**, 824.
3. Kusher, S.R. *et al.*, (1972) *Proc. Natl. Acad. Sci. USA* **69**, 1366.
4. Goldmark, P.J. and Linn, S. (1972) *J. Biol. Chem.* **247**, 1849.
5. Rosamond, J. *et al.*, (1979) *J. Biol. Chem.* **254**, 8646.

*Triton* is a registered trademark of Rohm & Haas, Philadelphia, Pennsylvania.

*OmniCleave* and *Plasmid-Safe* are trademarks of Epicentre, Madison, Wisconsin.

Visit our technical blog: [epicentral.blogspot.com](http://epicentral.blogspot.com)