6000TOCi & M800 Transmitter Quick Start Guide

English	3
Deutsch	17
Français	31
Italiano	45
Español	59
Česky	73
Dansk	87
Nederlands	101
Suomi	115
Magyar	129
Polski	143
Português	157
Svenska	171
Русский	185
中文	199
日本語	213
한국어	227
ไทย	241

6000TOCi & M800 Transmitter
30 473 835 Rev A

METTLER TOLEDO

Contents

1	Safety Instructions	4
2	Unpacking and Inspection	5
3	Dimensional Drawings	6
4	6000TOCi Sensor Installation and Start-up	6
5	Electrical Connection	8
6	Terminal Block (TB) definitions M800 Transmitter	8
7	6000TOCi Sensor	9
8	Menu Structure	9
9	Operation	10
10	General Setup	10
11	Disposal	12
12	Specifications	12
13	Service and Maintenance	13

Statement of intended use – The 6000TOCi is designed to measure Total Organic Carbon (TOC) concentration in high purity water applications. If this equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment against hazards may be impaired.

Warning!

Please refer to the 6000TOCi Operation Manual for safety instructions regarding installation and start-up. Follow all warnings, cautions, and instructions indicated on and supplied with this product.

1 Safety Instructions

This manual includes safety information with the following designations and formats.

Definition of equipment and documentation symbols and designations

WARNING: POTENTIAL RISK OF ELECTRIC SHOCK.

CAUTION: possible instrument damage or malfunction.

NOTE: Important operating information.

On the instrument indicates: There is alternating current present.

The following is a list of general safety instructions and warnings. Failure to adhere to these instructions can result in damage to the equipment and/or personal injury to the operator.

- The 6000TOCi Sensor should be installed and operated only by personnel familiar with the sensor and who are qualified for such work.
- The 6000TOCi Sensor must only be operated under the specified operating conditions.
- Repair or service of the 6000TOCi Sensor must be performed by authorized, trained personnel only.
- If this equipment is used in a manner not specified by the manufacturer, the protection provided by it against hazards may be impaired.
- Use only factory documented components for repair. METTLER TOLEDO accepts no responsibility for damage caused by unauthorized modifications to the sensor.
- Follow all warnings, cautions, and instructions indicated on and supplied with this product.
- Protective covers must be in place unless qualified personnel are performing maintenance.
- Protective covers must be in place at all times during normal operation.
- Install equipment as specified in this instruction manual. Follow appropriate local and national codes.
- Prior to shipping sensor back to the factory for repair or re-calibration, water **MUST** be drained from sensor to avoid damage due to freezing.

WARNINGS:

- Installation of cable connections and servicing of this product require access to shock hazard voltage levels. Switch or circuit breaker shall be in close proximity to the equipment and within easy reach of the OPERATOR; it shall be marked as the disconnecting device for the equipment.
- Main power must employ a switch or circuit breaker as the disconnecting device for the equipment. Electrical installation must be in accordance with the National Electrical Code and/or any other applicable national or local codes. Safety and performance require that this sensor be connected and properly grounded through a three-wire power source.

1 Safety Instructions continued

- Ozone gas (O_3) is generated inside the 6000TOCi Sensor enclosure during normal operation. The smell of ozone may be apparent when opening the front cover of the enclosure and caution should be taken when opening. Prolonged exposure to ozone gas is hazardous and may cause health problems.
- Electrical installation must be in accordance with the National Electric Code and/or any other applicable national or local codes.
- Safety and performance require that this instrument be connected and properly grounded through a three-wire power source.

NOTE! PROCESS UPSETS: Because process and safety conditions may depend on consistent operation of this sensor, provide appropriate means to maintain operation during sensor maintenance, replacement or sensor or instrument calibration.


2 Unpacking and Inspection

Inspect the shipping container. If it is damaged, contact the shipper immediately for instructions. Do not discard the container. If there is no apparent damage, carefully unpack the 6000TOCi Sensor. The box should contain the following items:

- 6000TOCi Sensor
- Installation Kit (see operation manual for list of contents)
- Operation Manual
- Quick Start Guide
- Product Insert Card
- Certificate of Test and Accuracy
- Declaration of Conformity

If items are missing or shipping container is damaged, notify METTLER TOLEDO immediately.

3 Dimensional Drawings

4 6000TOCi Sensor Installation and Start-up


Installation

Mount the 6000TOCi in an appropriate location. A minimum clearance dimension of 18" (45.7 cm) is required on the left side of the sensor for UV lamp change-out. 10" clearance above the sensor is required for CAL/SST. Two plumbing connections are provided on the right side of the enclosure labeled 'Sample Inlet' and 'Sample Outlet'.

4 6000TOCi Sensor Installation and Start-up continued

A third connection, a safety drain, is located on the bottom of the enclosure. Figure 1 shows a typical installed configuration. The installation kit provided includes installation tube assemblies, a sample conditioning coil and a high capacity inlet filter. Attach the high capacity inlet filter to the sample inlet connection as shown in Figure 2.

Figure 1


Connect the 0.125" (3 mm) PTFE tubing to a sample point isolation valve (not supplied). Warning: For inlet pressure above 85 psig (5.9 bar(g)), a pressure regulator (Thornton p/n 58 091 552) is required. Flush sample tubing before installing at the sensor. Close the isolation valve. Attach the other end of the PTFE tubing to the sample conditioning coil then attach the sample conditioning coil to the high capacity inlet filter. Connect the stainless steel drain tube (included) as shown in Figure 1 and direct outlet to an atmospheric drain.

Introduce Sample Flow

Slowly open the sample point isolation valve (not supplied), to start flow to the sensor. Allow approximately 3-5 minutes for the sensor to fill, and observe flow at the drain tube. Once there is flow, check that there are no leaks inside the sensor enclosure. Allow the sensor to initially flush with sample water (4-24 hours is recommended). Connect the M800 patch cable to the connector on the bottom left side of enclosure.

Figure 2


5 Electrical Connection

WARNING: MAIN POWER MUST EMPLOY A SWITCH OR CIRCUIT BREAKER AS THE DISCONNECTING DEVICE FOR THE EQUIPMENT. THE BREAKER SHOULD BE PROPERLY LABELED AND EASILY ACCESSIBLE FOR SYSTEM ISOLATION BY SERVICE PERSONNEL. THE CIRCUIT BREAKER SHALL HAVE A MAXIMUM OF 20 AMP RATING.

Two fittings are provided on the left side of the enclosure: the AC Power bulkhead fitting for power cable, and the M800 patch cable connection. Feed the power cable through the bulkhead fitting on the side-wall. Remove the two thumb screws protecting the TOC circuit board. Line and Neutral AC power connections are positioned on the top-left side of the TOC enclosure. The AC Ground connection is provided on a grounding plate located to the left of the Line and Neutral terminals. See Figure 3 for terminal wiring. See product label for specified fuse size.

Figure 3

6 Terminal Block (TB) definitions M800 Transmitter

6000TOCi Sensor Cable Connection

TOC Sensor uses 58 080 27X cables.


TB2 and TB4 – Terminal Assignment for 6000TOCi			
	TB2 (ISM Ch1,2)	TB4 (ISM Ch3,4)	6000TOCi
Term.	Function	Function	Sensor wire color
1	DI2+	DI6+	–
2	DI2–	DI6–	–
3	1-Wire_Ch1	1-Wire_Ch3	–
4	GND5V_Ch1	GND5V_Ch3	–
5	RS485B_Ch1	RS485B_Ch3	black
6	RS485A_Ch1	RS485A_Ch3	red
7	GND5V_Ch1	GND5V_Ch3	white
8	5V_Ch1	5V_Ch3	blue
9	24V_Ch2	24V_Ch4	–
10	GND24V_Ch2	GND24V_Ch4	–
11	1-Wire_Ch2	1-Wire_Ch4	–
12	GND5V_Ch2	GND5V_Ch4	–
13	RS485B_Ch2	RS485B_Ch4	black
14	RS485A_Ch2	RS485A_Ch4	red
15	GND5V_Ch2	GND5V_Ch4	white
16	5V_Ch2	5V_Ch4	blue

*For M800DP Transmitters equipped with Profibus DP communication capabilities, please refer to the M800DP Transmitter Operation Manual.

7 6000TOCi Sensor

The 6000TOCi Sensor housing provides four LED status lights and a single push button for local UV lamp control. See Figure 4.

Figure 4

Fault LED: Flashing when a Fault condition exists.


Error LED: Illuminates continuously when an Error condition exists.

Sensor status LED: Illuminates whenever AC power is on and the M800 and 6000TOCi are connected.

UV Lamp ON LED: When illuminated continuously, UV Lamp is on. When flashing, indicates sensor is in rinse and auto-balance mode.

UV Lamp ON key: Allows the user to turn UV lamp on or off at the sensor.

8 Menu Structure

9 Operation

Navigating the Menu Tree

To browse through the menus, touch the open arrows > and/or <.

To access a menu press the corresponding arrow ► in the same line.

Exit a Menu

Press  to exit the menu.

Press  to return to the Measurement Screen.

Return to the Measurement Screen

Press  to return to the Measurement Screen.

Press  to logout from the transmitter.

Entry of Data

The M800 displays a keypad for modifying values.

Press the  button and the transmitter will store the new value.

Press the "ESC" button to exit the keypad without changing data.

Selection Menus

Some menus require a selection of a parameter/data. In this case the transmitter displays a pop up window. Press the according field to select the value. The pop window will be closed and the selection will be stored

"Save changes" Dialog

If the M800 brings up the "SAVE CHANGES" dialog there are the following options.

"NO" will discard the entered values and exit the menu, "YES" will save changes made and "CANCEL" will return to the menu to continue configuring.


10 General Setup (PATH: /Guided Setup)

From the main screen, press the Configure icon () and select Guided Setup. From this screen, you will be able to select and configure the first measurement from the 6000TOCi that will be displayed by the M800, and configure and activate an analog output and/or set point to be associated with this measurement.

Analog outputs

By selecting Yes the linear 4-20 mA analog output Aout1 will be set up for measurement 1 when  is pressed. Selecting No will deactivate the analog output. Aout1 min, Aout1 max are the minimum and maximum measurement values for the 4 and 20 mA values respectively. Additional options may be configured from the analog output configuration menu.

10 General Setup continued

Set Points

On page 2 of the guided setup menu, a set point can be configured for that measurement. Once the Set Point option is set to Yes, then select the set point type you wish to activate. If Off is selected, the set point will be disabled. Once a set point type is selected, the necessary parameters for that set point can also be configured.


The following set point types can be selected:

- High (High value has to be set)
- Low (Low value has to be set)
- Between (High and Low value has to be set)
- Outside (High and Low value has to be set)

For Conductivity, the following additional options may be selected:

- USP (% safety margin below U.S. Pharmacopoeia limits)
- EP PW (% safety margin below European Pharmacopoeia limits for Purified Water)
- EPWFI (% safety margin below European Pharmacopoeia limits for Water for Injection)
- JP Cond (% safety margin below Japanese Pharmacopoeia limits)

After setting the set point value(s) a Relay can be selected from the relays available on the transmitter, and configured for that set point. The Relay delay is set to 10 seconds and the Hysteresis is set to 5%. These parameters can be configured from the Setpoints menu.

UV Lamp

Once flow has been established, the 6000TOCi is ready to begin taking TOC measurements. If sensor activation is desired at this time, set the UV Lamp to On. Allow approximately one minute for readings to begin. Additional time (approximately 4-24 hours) may be necessary for the readings to stabilize as the sample line rinses clean, and the TOC sensor reaches thermal stability.

Auto Start

The 6000TOCi sensor can be configured such that the UV lamp will automatically turn on after recovering from a power loss or fault condition. If automatic recovery is desired, set the Auto Start option to 'ON'. Auto Start can be configured in Configure (*), measurement parameter setting

It is recommended that Auto Start be set to 'ON'. The default setting is 'OFF'.

11 Disposal

When the sensor is finally removed from service, observe all local environmental regulations for proper disposal. Please note the UV lamp supplied with the 6000TOCI contains mercury.

12 Specifications

TOC Sensor

Measurement Range	0.05 - 2000 ppbC ($\mu\text{gC/L}$)
Accuracy	± 0.1 ppbC for TOC < 2.0 ppbC (for water quality > 15 M Ω -cm [0.067 $\mu\text{S/cm}$]) ± 0.2 ppbC for TOC > 2.0 ppbC and < 10.0 ppbC (for water quality > 15 M Ω -cm [0.067 $\mu\text{S/cm}$]) $\pm 5\%$ of measurement for TOC > 10.0 ppbC (for water quality 0.5 to 18.2 M Ω -cm [2.0 to 0.055 $\mu\text{S/cm}$])
Repeatability	± 0.05 ppbC < 5 ppbC, $\pm 1.0\%$ > 5 ppbC
Resolution	0.001 ppbC ($\mu\text{gC/L}$)
Analysis Time	Continuous
Initial Response Time	< 60 seconds
Update Rate	1 second
Limit of Detection	0.025 ppbC

Conductivity Sensor

Conductivity Accuracy	$\pm 2\%$, 0.02 - 20 $\mu\text{S/cm}$ $\pm 3\%$, 20 - 100 $\mu\text{S/cm}$; 0.1 Constant Sensor*
Cell Constant Accuracy	$\pm 2\%$
Temperature Sensor	Pt1000 RTD, Class A
Temperature Accuracy	$\pm 0.25^\circ\text{C}$

Sample Water Requirements

Temperature	0 to 100 $^\circ\text{C}$ **
Particle Size	< 100 micron
Minimum Water Quality	> 0.5 M Ω -cm (< 2 $\mu\text{S/cm}$), pH < 7.5 ***
Flow Rate	> 8.5 mL/min
Pressure	4 to 200 psig (0.3 bar(g) to 13.6 bar(g)) at sample inlet connection ****

General Specifications

Case Dimensions	11.9" [302.75 mm] W x 9" [229.8 mm] H x 5.7" [144.7 mm] D
Weight	11.0 lb. (5 kg)
Enclosure Material	Ignition Resistant Polystyrene Resin meeting UL 94V-0, Painted Aluminum
Ambient Temperature/Humidity Rating	5 to 50 $^\circ\text{C}$ / 5 to 80% Humidity, non-condensing
Power Requirements	100 - 240VAC, 50 - 60 Hz, 25W
Local Indicators	Four LED lights for Fault, Error, Sensor Status and UV Lamp ON
Ratings/Approvals	CE Compliant, UL and cUL (CSA Standards) listed, Conductivity and temperature sensors traceable to NIST, ASTM D1125 and D5391. Meets ASTM D5173 Standard Test Method for On-Line Monitoring of Carbon Compounds in Water by UV Light Oxidation

12 Specifications continued

Installation/Power/Enclosure

Inlet Connection	0.125" [3 mm] O.D. (6' [2 m] FDA compliant PTFE tubing supplied)
Outlet Connection	0.125" [3 mm] O.D. (6.5" [165 mm] fixed 316SS tube provided)
Inlet Filter	316SS, inline 60 micron
Wetted Parts	316SS/Quartz/PEEK/Titanium/PTFE/EPDM
Wall Mount	Standard, mounting bracket provided
Maximum Sensor Distance	300ft [91 m]

* Readout in equivalent S/m ranges selectable at M800

** Temperature above 70°C requires Sample Conditioning Coil (included)

*** For power plant cycle chemistry samples, pH may be adjusted by measurement after cation exchange.

**** Process pressure above 85 psig (5.9 bar(g)) requires optional High Pressure Regulator p/n 58 091 552.

Specifications subject to change without notice.

13 Service and Maintenance

The 6000TOCi Sensor is designed to minimize service and maintenance. There are virtually no moving mechanical components, therefore normal operating wear and tear is minimal. This reduces the amount of consumable components as well as the amount of time needed to maintain the sensor. Listed below are instructions on how to perform simple periodic maintenance, which includes UV Lamp change (every 4500 hours of operation), filter replacement (typically every 6 months), and general cleaning.

UV Lamp Replacement


WARNING: UV RADIATION HAZARD

Apply power to UV lamp only when installed in housing in accordance with operation manual. DO NOT remove UV lamp from housing unless power is off.

Always protect eyes and skin from direct exposure to UV light.

METTLER TOLEDO Thornton recommends replacement of the UV lamp inside the 6000TOCi Sensor after 4500 hours of operation, or 6 months of continuous use, not to exceed one year. This is a simple procedure that requires only a few minutes to complete. The following steps explain the proper procedure for the UV lamp replacement. Refer to the illustration below.

13 Service and Maintenance continued


CAUTION: Use of a UV lamp other than those provided by METTLER TOLEDO Thornton specifically for use with the 6000TOCi Sensor will affect performance and void the warranty of this product.

1. At the sensor, turn off the UV lamp by pressing the lamp control button (UV lamp ON LED will turn off). If the LED does not turn off, check that the Sensor Key Lock is in the off position in the M800. Refer to **Parameter Settings: Sensor Key Lock** function in the operation manual.
2. Once power to the UV lamp is off, remove the UV lamp replacement side access panel on the left side of the sensor using a Phillips head screwdriver.
3. Disconnect the power cable to the UV lamp.
4. Loosen the black UV lamp holding screw.
5. Slide the cable of the UV lamp through the side opening of the enclosure and gently slide the UV lamp out of the oxidation chamber assembly (stainless steel cylinder). Be careful not to let the UV lamp hit the quartz glass coil inside the chamber.
6. Use the gloves supplied with each replacement bulb. Hold the new lamp from the ends of the lamp. Do not touch the clear glass portion of the bulb. Slide the new UV lamp into the side opening of the enclosure and into the oxidation chamber opening until it stops. Do not use excessive force to insert the UV lamp as this may cause damage to the lamp or the internal components of the oxidation chamber.
7. Tighten the UV lamp holding screw until snug. Do not overtighten.

CAUTION: Do not overtighten holding screw.

8. Feed the power cable through the side opening of the enclosure. Re-connect the cable to the power connector.
9. At the M800, enter the TOC Maintenance menu at Path: Calibrate (L₁)/TOC Maintenance ► Replace Lamp button to register the lamp replacement.

10. After the lamp replacement has been completed, a TOC calibration must be performed. In addition, a System Suitability Test is recommended for applications where compliance with Pharmacopeia regulations is necessary.


Front Panel Cleaning

Clean the front panel with a damp soft cloth (water only, no solvents). Gently wipe the surface and dry with a soft cloth.

13 Service and Maintenance continued

Technical Support


For technical support and product information for the 6000TOCi, contact your local METTLER TOLEDO Sales Office or representative.

High Capacity Inlet Filter Replacement

The 6000TOCi Sensor includes a high capacity filter (shown below) containing a filter element that should be replaced (P/N 58 091 551, package of 2) after approximately 6 months of operation, depending on water quality conditions. Detailed instructions to replace this filter are included in the replacement package.

Once the filter replacement has been completed, the replacement must be registered with the 6000TOCi to reset the filter operating information. At the M800, enter the TOC Maintenance menu at Path: Calibrate (L*)/TOC Maintenance ► Replace Filter button to register the filter replacement.

REPLACING THE FILTER ELEMENT FOR THE HIGH CAPACITY INLET FILTER

Draining and Shipping Instructions

The 6000TOCi contains a glass coil which facilitates oxidation of the water sample. Water freezing within the coil during transportation will cause this coil to break. Therefore, it is important that all water be removed from the TOC sensor prior to shipment.

To drain the water, connect the TOC sensor to AC power. To put the flow control valve into the 'Drain' mode, press and hold the UV Lamp button for 5 seconds until the Sensor Status LED begins flashing. This mode will fully open the automatic flow control valve to allow draining the sensor more easily. Once the valve has been placed into the 'Drain' mode, use the supplied plastic syringe to push the water out of the sensor, blowing from the inlet, and pushing water out through the sensor outlet.

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© Mettler-Toledo Thornton 06/2018. Printed in United States.
- Switzerland** **Mettler-Toledo (Schweiz) GmbH**
METTLER TOLEDO, GmbH
METTLER TOLEDO Thornton, Inc.