

Flexible single-use TFF systems for full-scale production

A unique open system design provides complete flexibility to meet your exact process requirements, while every part of the TFF system features best-inclass components to optimize performance.



Engineered to optimize your process

The Agilitech Advantage:

- Low shear design maximizes yield
- Optimized TMP control streamlines your process
- Easily run processes with more viscous solutions
- Reduced cost of consumables

A magnetically levitated pump combined with a cubical mixing bag with magnetic mixer limits shear on the product throughout the filtration process, maximizing yield at every TFF step. The pump maintains optimal transmembrane pressure to optimize

Technical Specifications

Attribute	T750	T1000
General Skid Specific	cations	
Pressure rating	60 psi/4 bar	60 psi/4 bar
Temperature Range	2-50°C	2-50°C
Humidity	10-90%	10 -90%
(non-condensing)		
Utility Requirements*		
Skid Power	220-230 VAC	220-230 VAC
Mixer Power	100-120 VAC	100-120 VAC
Wetted	USP <88> Class VI	USP <88> Class VI
Materials	and animal free	and animal free
Membrane	$0.5 m^2 to 5.0 m^2$	2.5 m ² to 10.0 m ²
Holder Capacity		
Feed Flow Rate	5-50 LPM	10-100 LPM
Feed and Filtrate	Clamp on ultrasonic	Clamp on ultrasonic
Flow Measurement	flow sensor (±3%)	flow sensor (±3%)
Temp Control	2-50°C	2-50°C
Inlets	1 product	To recirc bag: 3 inlets,
	2 buffers	1 product, 2 buffers
		To feed pump: 3 inlets,
		1 recirculating bag,
		2 additional feed ports
Tubing Internal		
Diameter		
Feed and Retentate:	0.75 in	1 in
Permeate:	0.5 in	0.75 in
Inlets:	0.5 in	0.75 in

Specifications represent standard Agilitech models and may differ when customized for specific customer applications.

filter performance throughout the TFF process. Combined with an air modulated diaphragm pressure control valve, retentate pressure reaches set point almost instantly to achieve the desired TMP. These design features work to eliminate pump pulsations and pinching control valves that lead to more erratic pressure control.

The system maxmizes product recovery in the recirculation bag. Users simply disconnect the tubing from the TFF skid and move the recirculation bag to the next step for further processing. This results in fewer consumables and reduced cost, while eliminating project delays due to supply chain issues with consumables.

Attribute	T750	T1000
Outlet Ports	1 Filtrate	Permeate: 2 outlets; waste or recycle Retentate; 3 outlets waste, recycle or recirculating bag
Tank Size	Up to 650 L	Up to 1000 L
Minimum working volume based on: Agitator Top Agitator Bottom Bag Sensors/Sampling	200 L Bag 37 L 14 L 25 L	650 L Bag 82 L 31 L 55 L
Hold-Up Volume (excludes filtration device)	1.5 L	2.5 L
Dims (W x D x H)** Skid, mm Bag Holder, mm Sensor Performance	1219 x 787 x 1810 1219 x 825 x 1595	2032 × 1270 × 2108 1371 × 1066 × 1595
Pressure	0-6 psi (± 2%) 6-30 psi (± 3%) 30-60 psi (± 5%)	0-6 psi (± 2%) 6-30 psi (± 3%) 30-60 psi (± 5%)
Filtrate Conductivity	0-250 mS/cm (±0.25 mS/cm)	0-250 mS/cm (± 0.25 mS/cm)
Vessel Weight	0-650 kg (± 0.3%)	0-1000 kg (± 0.3%)
Flow	0 - 50 LPM (± 3%)	0-100 LPM (± 3%)

*Compressed air (oil free) @ 80 psi/5.5 bar

**Highest Operator interaction: 5 ft 2 in (T750); 5 ft 8 in (T1000)

For more information or product support: +1.661.381.7950

info@agilitech.bio www.agilitech.bio