pulse oximetry





SMARTsat technology – Now you have a choice! High performance SpO₂ technology for OEMs

The SMARTsat technology platform was developed over the last few years, in close cooperation with wellestablished research institutions and university hospitals.

The very latest and innovative signal processing technologies and algorithms enable precise measurements, even under very difficult physiological conditions. SMARTsat technology has extremely low power consumption requirements, and is a real alternative for OEM customers whom require highend pulse oximetry solutions for their applications.

In accordance with ISO Standards, the approved pulse oximetry sensors are calibrated and evaluated against dyshaemoglobin-free reference measurements, which were determined from CO oximeter data and do not contain saturation components of the haemoglobin fractions SaCO and SaMet.

SMARTsat® OEM I & II SpO₂ accuracy validation results

	A_{RMS} in the SaO ₂ ranges ²⁾						
Sensor Type	60 - 100%	70 - 100%	90 - 100%	80 - 90%	70 - 80%	60 - 70%	
SF7500VM Finger Clip Sensor	1.6	1.5	1.4	1.3	1.7	2.3	
SC7500VM ¹⁾ Silicone Soft Sensor	1.6	1.6	1.1	1.7	1.8	1.9	
W7500VM Wrap Sensor	1.6	1.4	1.0	1.2	1.8	2.8	
EP7500VM Ear Clip Sensor	2.2	2.1	1.0	1.7	2.8	3.5	
10-AP-VM Disposable Sensor - Adult	2.1	2.1	2.2	1.9	2.0	2.2	

SMARTsat® OEM III SpO, accuracy validation results

A_{RMS} in the SaO ₂ ranges ²⁾						
Sensor Type	60 - 100%	70 - 100%	90 - 100%	80 - 90%	70 - 80%	60 - 70%
SF7500VM Finger Clip Sensor	1.6	1.5	1.4	1.4	1.7	2.0
SC7500VM ¹⁾ Silicone Soft Sensor	1.7	1.6	1.5	1.8	1.7	2.4
W7500VM Wrap Sensor	1.7	1.6	1.3	1.5	1.9	2.4
EP7500VM Ear Clip Sensor	2.4	2.3	1.4	1.6	3.2	3.3
10-AP-VM Disposable Sensor - Adult	2.5	2.4	2.1	2.3	2.6	3.4

1) Accuracy statement applies to SC7500VM, SCM7500VM and SCP7500VM

2) As inherent to their functional principle, pulse oximetry measurements underlie statistical spread; therefore only two thirds of the measurement data are within the specific range of +/- ARMS



Typical applications

SMARTsat is designed for integration into:

- High-End SpO₂ monitoring devices
- Patient monitors
- Transport monitors
- Sleep screening devices
- Defibrillators
- Home care monitors
 - ... and more



Robust, flexible and hygenic

Viamed provides a range of high quality and comfortable sensors that incorporate the very latest technology



Soft Silicone Sensors

Due to their robust design, the Soft Silicone Sensors are ideal for use in the tough environment of rescue services, emergency services and hospitals.

The Silicone Sensors from the 4000 series are manufacturered from premium materials and their design allows for effective high-level disinfection. This reduces the risk of nosocomial infection associated with surface-borne pathogenic microorganisms.

Available in Adult (SC7500VM), Small Adult (Medium) (SCM7500VM) and Paediatric (SCP7500VM) sizes.



Finger Clip sensors are the ideal solution for ambulatory use, or for long-term monitoring. Advanced manufacturing technologies, materials and design elements are found in these sensors. (SF7500VM)

Ear Sensor

High-quality ear sensor, which is comfortable and easy to use. The ear sensor has an additional sensor holder that prevents sensor movement and prevents it from being dislodged. (EP7500VM)



Silicone Wrap Sensor

Silicone Wrap sensors are for universal applications, for use on infants to adults. Wrap tapes are provided with mounting loops, enabling a simple fixation of the skin friendly wrap sensor. (W7500VM)

Disposable Sensors

Disposable sensors with Plaster Tape – the soft material is highly flexible and adaptable. Available in adult (10-AP-VM), paediatric (10-PP-VM), infant (IP-10-VM) and neonatal (10-NP-VM) sizes.



Sensor Wraps Pack of 12



Multi-site Y Senso Single Patient Use



Hydrogel Pads For use with Y sensor



Extension Cables Two standard lengths: 1.2m (XT6500) and 2.4m (XT6501)

Specification			OEM II⁴	OEM III ⁴	EB ⁴		
Measurement Range	SpO ₂	0 - 100 %	0 - 100 %	0 - 100 %	0 - 100 %		
	Pulse Rate	20 - 300 bpm	20 - 300 bpm	20 - 300 bpm	20 - 300 bpm		
	Perfusion	0.02 - 20 % (no motion)	0.02 - 20 % (no motion)	0.1 - 20 % (no motion)	0.1 - 20 % (no motion)		
Accuracy	SpO ₂	60 - 100 % +/- 2 A_{RMS} (no motion) ¹ 60 - 100 % +/- 3 A_{RMS} (motion condition) ² 60 - 100 % +/- 2 A_{RMS} (low perfusion) ³					
	Pulse Rate	20 - 300 BPM +/- 2 A _{RMS} (no motion)					
Power requirements							
Input voltage (ripple < 1	0mV)	4.5 – 5.5 VDC	3.2 – 3.6 VDC	3.1 – 3.6 VDC	4.5 – 5.25 VDC		
Average power consum	ption	110 mW @ 5V⁵	70 mW @ 3.3V⁵	24 mW @ 3.3V⁵	106 mW @ 5V ^{5, 6} 47 mW @ 5V ^{5, 7}		
Maximum power consu	mption	131 mW @ 5V⁵	83 mW @ 3.3V⁵	29 mW @ 3.3V⁵	125 mW @ 5V ^{5, 6} 55 mW @ 5V ^{5, 7}		
Environmental condit	ions						
Operation	-:	25 to +60 °C					
Storage		-40 to +70 °C					
Relative humidity	1 1	15 % - 95 % (operation, non-condensing) 10 % - 95 % (storage, non-condensing)					
Serial communication	1						
Baud rate	9	9600 baud up to 230400 baud					
Sampling rate adjusta	ible 7	5 Hz, 150 Hz		75 Hz, 150 Hz, 300 Hz			
Applied Standards							
SMARTsat module OEM I, II, III IEC 60601-1:2005 (3rd Ed); IEC 60601-1 (2nd Ed) IEC 60601-1-2:2014 (4th Ed); IEC 60601-1-2 (3rd Ed) EN ISO 80601-2-61:2011; EN ISO 9919:2009 ISO 14971:2007; IEC 60601-1-6:2010; IEC 60601-1-11:2010; IEC 62471:2006							
SMARTsat sensors	ľ	ISO 10993-1:2009; ISO 10993-5:2009; ISO 10993-10:2010					
Dimensions (L x W x H)		OEM I 63.5 x 50.8 x 6.3 mm	OEM II 37.0 x 14.0 x 5.0 mm	OEM III 31.0 x 14.0 x 5.0 mm	EB 61.0 x 26.0 x 16.0 mm		
Part Numbers							
OEM Module		0012700	0012710	0012720	0012730		
SMARTsat development OEM module, Development pla sensor flex cable, USB cable & SoftFlap Finger Cilp SpO2 sen	nt kit atform (includes: PC software), sor SF7500VM	0012760	0012761	0012762			
SMARTsat sensor plug			0012745				
SMARTsat sensor plug - PCB mounting version		0012746					
SMARTsat sensor basic flex cable		0012749					
SMARTsat® sensor low	v noise flex ca	ble	0012750				

¹ A_{rms}: ± 1 Arms represents approx 68% of measurements

² Tested with all motion patterns Index II

³ Tested with Prosim 8 simulator

⁴ The baud rate has significant influence on the total power consumption of the SMARTsat® OEM I, OEM II and OEM III modules.



⁶EB with OEM II,

7 EB with OEM III

Reduced data transfer rates and/or higher sampling rates (150 Hz, 300 Hz) increases power consumption.



Viamed Limited · 15 Station Road · Cross Hills Keighley · West Yorkshire · BD20 7DT · United Kingdom Tel: +44 (0)1535 634 542 Fax: +44 (0)1535 635 582 Email: info@viamed.co.uk Website: www.viamed.co.uk

