

Multiple drug configurations to meet your needs.

Alere iScreen® Drugs of Abuse Test Card

Urine Drug Screening Device



For professional in vitro diagnostic use only


Knowing now matters.™

Alere iScreen® Drugs of Abuse Test Card Urine Drug Screen

Features and Benefits

- 1-12 configurations available for a fully customized urine testing program
- Simple procedure: collect, dip, and read results at 5 minutes
- Administrator remains in control of test device
- Try the Multipanel Drugs of Abuse Test Card for expanded specimen validity testing: CRE, GLUT, NIT, OX, pH, S.G.*


Test in 3 Easy Steps Single and a multidrug test card procedures vary slightly. See package insert for details.



1
Collect specimen and remove plastic cap at bottom.



2
Dip revealed pads into sample for 10 seconds.



3
Replace cap on device, lay flat, and read results in 5 minutes.



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings. For healthcare professionals including professionals at point of care sites. Immunoassay for in vitro diagnostic use only. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Drug screen options and cutoff levels

Test	Calibrator	Cutoff Level
Amphetamine (AMP 1000)	d-Amphetamine	1,000 ng/mL
Amphetamine (AMP 300)	d-Amphetamine	300 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cocaine (COC 300)	Benzoyllecgonine	300 ng/mL
Cocaine (COC 150)	Benzoyllecgonine	150 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9-COOH	50 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (mAMP 1000)	d-Methamphetamine	1,000 ng/mL
Methamphetamine (mAMP 500)	d-Methamphetamine	500 ng/mL
Methylenedioxymethamphetamine (MDMA - Ecstasy)	d,l-Methylenedioxymethamphetamine	500 ng/mL
Opiate (MOP 300)	Morphine	300 ng/mL
Opiate (OPI 2000)	Morphine	2,000 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000 ng/mL

*Specimen Validity Tests: Creatinine (CRE), Glutaraldehyde (GLUT), Nitrite (NIT), Oxidants/PCC (OX), pH (pH), Specific Gravity (S.G.)



Products 800.340.4029 | aleretoxicology.com