

**Warning**

Do not use any other sterilization method

APPROVED DISINFECTANTS

LIQUIDS

ORION Laboratories

70% Isopropyl Alcohol*

Do not use 70% Isopropyl Alcohol on STARBoard with an IBN product code lower than IBN001536*Jaychem Industries**2% Chlorhexidine / 70% Alcohol
(Pink)**Betadine**

Povidone Iodine (7.5%)

C-Tec

Nature Clean (Diluted 1:10 with water)

BODE Chemie GmbH

Kohrsolin FF

STERIS

Coverage Spray TB*

Do not use STERIS Coverage TB on STARBoard with an IBN product code lower than IBN001536*CaviCide**

by Metrex Research

DISSOLVABLE

Du Pont

Rely on Virkon Tablets

SURFACE WIPES

PDI Sani-Cloth Bleach**PDI Sani-Cloth AF3****Reynard Health Supplies**

Surface Disinfectant Wipes

Diversey

Oxivir Tb Wipes

CaviWipes

by Metrex Research

Clinell Wipes**Discoloration**

Some disinfectants may cause slight discoloration to soft blue components. This will not affect the strength and will remain fit for purpose.

**Warning**

Do not use any other sterilization method

APPROVED CLEANING INGREDIENTS

Please refer to your disinfectant labelling to ensure all ingredients are listed below before usage.

Acid Based:

- ≤10% Malic Acid CAS 6915-15-7
- ≤6% Sulfamic acid CAS 5329-14-6

Alcohol Based:

- ≤5% Benzyl Alcohol CAS 100-51-6
- ≤5% 2-Butoxyethanol CAS 111-76-2
- ≤10% Butyldiglycol CAS 112-34-5
- ≤70% Isopropyl Alcohol CAS 67-63-0
- ≤10% Tridecanol CAS 69011-36-5

Ammonium and Chloride Based:

- ≤5% Benzalkonium Chloride CAS 68424-85-1
- ≤0.28% Benzethonium Chloride CAS 121-54-0
- ≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
- ≤2% Chlorhexidine CAS 55-56-1
- ≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
- ≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6

Sulphate Based:

- ≤0.1% PHMB CAS 27083-27-8
- ≤55% Potassium Peroxymonosulfate CAS 70693-62-8
- ≤3% Potassium Persulfate CAS 7727-21-1
- ≤12% Sodium Alkylbenzene Sulfonate CAS 68411-30-3
- ≤5% Sodium Toluene Sulphonate CAS 12068-03-0

Others

- ≤10% Glutaral CAS 111-30-8
- ≤0.15% Hydrogen Peroxide CAS 7722-84-1
- ≤0.63% Sodium Hypochlorite CAS 7681-52-9
- ≤7.2% Povidone Iodine CAS 25655-41-8

Discoloration

Some disinfectants may cause slight discoloration to soft blue components. This will not affect the strength and will remain fit for purpose.