


RECOVERY PATCH

COMMON AREAS OF APPLICATION*



*Not clinically proven for all injuries



www.kttape.com

RECOVERY PATCH


KT Recovery Patches are specially designed to reduce swelling and inflammation and its accompanying pain.*

APPLY PATCH

With body part in stretched position, apply patch(es) over swollen or inflamed* area. Do not stretch patch during application.

AFTER

Affected area shows reduction in swelling or inflammation*



For more information, go to kttape.com/recoverypatch

BEFORE YOU START

Here are some important tips to get you the best possible results.



APPLY AT ONSET
Apply within 48 hours of initial onset of swelling or inflammation.



CLEAN SKIN
Clean dirt and oils from skin and allow to dry.



TRIM HAIR
Trim excessive hair for better adhesion.



SKIN ON STRETCH
Position the area of application so that the skin is on stretch.



DON'T STRETCH
Do not apply any stretch to this patch. Stretching patch may result in skin irritation or blisters.




DON'T RIP OFF!
Gently remove patch from skin; use baby oil if necessary.

HOW TO APPLY RECOVERY PATCHES

APPLY AS NEEDED - CAN BE WORN UP TO 3 DAYS


1 Pull tab to remove the inner backing paper.



2 With the body part on stretch, place the patch on the skin with the "KT" towards the heart.



3 Pull the backing paper frame away slowly and smooth the edges. Rub patch to ensure full adhesion.



ADDITIONAL APPLICATIONS

For greater coverage, use these other applications

END-TO-END



SIDE-TO-SIDE



CRISS-CROSS



VIDEO INSTRUCTIONS



www.kttape.com/instructions

CUSTOMER SUPPORT

KT Recovery Patches are not a replacement for professional medical care. Consumers are encouraged to consult a medical professional before using KT Recovery Patches. Discontinue use and seek medical consultation if you develop worsening skin conditions, or experience any new pain, tingling, or numbness. Do not use on abdomen if pregnant. Do not place KT Recovery Patches directly over cuts or open wounds. Individuals with known skin sensitivity due to active infection, radiation or other burns, dermatologic conditions, or a compromised lymphatic system, congestive heart failure, renal insufficiency, cancer, or deep vein thrombosis should consult their medical doctor before using this product. Warranties and remedies limited to product replacement cost.



Obelis s.a. - Bd Général Wahis 53- 1030 Brussels, Belgium
Tel - +(32) 2 732-59-54 Fax - +(32) 2 732-60-03 Email - mail@obelis.net

KT Health, LLC - 584 East 1000 South, Suite 4 American Fork, UT 84003 U.S.A. • www.kttape.com
Email • support@kttape.com • ©2023 KT Health, LLC. All rights reserved.

