

......

HOW TO REMOVE TAPE (DO NOT RIP OFF)





Begin by peeling tape up at a corner.

Hold skin down where tape is being pulled off with one hand while gently pulling tape off with the other hand. Remove slowly.

PREPARE Trust KT to get you ready for your next goal. PERFORM Confidence so you can be at your best. RECOVER Amplify your recovery and be ready for what's next.

TRUSTED KT PRODUCTS

Prepare, Perform, Recover, Repeat.

	PREPARE	PERFORM	RECOVER
Kinesiology Tape	~	~	~
Creams/Gels	~	\checkmark	\checkmark
Blister Relief	~	~	~
Chafing Relief	~	~	
Ice/Heat		~	~

VIEW ALL PRODUCTS @ KTTAPE.COM

WARNING: For external use only. Do not use this product if you have poor circulation, fragile skin, or are diabetic. Do not use if you are allergic to any component of this product. Do not apply to any infected areas. Product should not be used to treat already formed bilsters. Only apply each strip one time. Discontinue use if skin becomes irritated. Do not use product if paper backing has been damaged or removed. This product is not made with natural rubber latex. One or more of the components of this product is a dyed material. There is a possibility that some color migration may occur during use. READ ALL cautions AND INSTRUCTIONS PRIOR TO USE. NOTE: Store away from sunlight in a dry, cool environment between 20°C to 25°C (68°F to 77°F). Questions? Email us at support@kttape.com



NEW USER GUIDE



KT TAPE.COM