


CGM PATCH

KTTM TAPE

PRO EXTREME[®] WATER RESISTANT ADHESIVE



1 BEFORE YOU APPLY

Follow all directions and use the following tips for best results and enhanced water-resistance.



APPLY BEFORE ACTIVITIES
Apply patch overnight before exercising or exposing the patch to water.



CLEAN SKIN
Shower/bathe to remove dirt and oils from the skin. Clean the area with alcohol wipes and allow to dry.



TRIM HAIR
Remove excessive hair for better adhesion.





AVOID TOUCHING ADHESIVE
Avoid touching the adhesive side of the patch to prolong adhesion.



CAUTION: If you have skin sensitivities, cancer, or are pregnant, consult your doctor before use. Improper application or applying on areas of sensitive skin may result in skin irritation.

2 HOW TO APPLY



3 AFTER YOU APPLY

Follow these care tips to ensure best results.



DO NOT DRY WITH HEAT


Do not use a hair dryer to dry the patch when it becomes wet. Instead, gently pat with a towel.



TRIM LIFTING EDGES

Trim with scissors any edges that begin to peel back.

4 HOW TO REMOVE



FIRST

Begin by peeling the patch up at the edge.

SECOND

For easier removal, try applying baby oil to the patch.

THIRD

Gently pull tape back - not up. Remove slowly.



SUPPORT@KTTAPE.COM
TEXT US (385) 376-1811

LEARN MORE @ [KTTAPE.COM/CGMPATCH](https://www.kttape.com/cgmpatch)

