



A proudly Canadian company

# KN95 Protective Mask

Bacterial Filtration Efficiency  $\geq 95\%$

*Redacted Version*



## White and black colour options available



## Medical Standard GB19083-2010 & Commercial Standard GB2626-2006 Available

Table 2: KN95 specifications

	Medical KN95 (GB 19083-2010)	Commercial KN95 (GB2626-2006)
Labelling Information on respirator	<ul style="list-style-type: none"> <li>• Number and year of publication of standard</li> <li>• Type and grade of filter elements (i.e. KN95)</li> </ul>	
Nasal splints	Must be adjustable	N/A
Mask Harness	<ul style="list-style-type: none"> <li>• No less than 10 Newton (1kg) per strap</li> <li>• Earloops Are Allowed</li> </ul>	
Fit	The mask shall be designed as to provide a good fit and overall fit factor of the mask shall <b>not be less than 100</b>	The disposable facepiece structure shall ensure the tight fitting with face, and be free from deformation during the service life.
Filter performance – (must be ≥ X% efficient)	≥ 95%	≥ 95%
Test agent	Sodium Chloride (NaCl) Particles	Sodium Chloride (NaCl) Particles

## Medical Packaging Option

### Full Case Count

- 200 units/ case
- Masks individually wrapped
- Black and white options available
- 85 cases/ pallet

### Full Case Dimensions


- 11.5 x 10.5 x 10 inches



## Retail Packaging Option

- Case of 20 in a box
- Includes 4 packs of 5 masks
- Black and white options available





Sponsor:

## Sodium Chloride (NaCl) Aerosol Test Final Report

---

**Test Article:** Respirator Face Mask  
**Study Number:**   
**Study Received Date:** 03 Feb 2016  
**Test Procedure(s):** Standard Test Protocol (STP) Number: STP0014 Rev 07


**Summary:** This procedure was performed to evaluate particulate filter penetration as specified in 42 CFR Part 84 and TEB-APR-STP-0059 for requirements on a N95 respirator. Respirators were conditioned then tested for particle penetration against a polydispersed, sodium chloride (NaCl) particulate aerosol. The challenge aerosol was dried, neutralized, and passed through the test article at a concentration not exceeding 200 mg/m<sup>3</sup>. The initial airflow resistance and particle penetration for each respirator was determined.

According to 42 CFR Part 84.64, pretesting must be performed by all applicants as part of the application process with NIOSH. Results shown below are part of that pretesting and must be submitted to and accepted by NIOSH for respirator approval.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

**Results:** The NIOSH N95 filter efficiency as stated in 42 CFR Part 84.181 is a minimum efficiency for each filter of ≥95% (≥5% penetration). The test articles submitted by the sponsor conform to the NIOSH N95 criteria for filter efficiency.

Test Article Number	Initial Airflow Resistance (mm H <sub>2</sub> O)	Particle Penetration (%)	Filtration Efficiency (%)
1	6.9	1.02	98.98
2	6.8	1.23	98.77
3	6.8	1.24	98.76
4	7.0	1.08	98.92
5	6.6	1.12	98.88
6	7.3	1.12	98.88
7	7.1	1.03	98.97
8	6.6	1.13	98.87
9	7.4	1.06	98.94
10	6.9	1.17	98.83
11	6.8	1.18	98.82
12	7.1	1.25	98.75
13	7.5	1.14	98.86


  
 Study Director

Brandon L. Williams  
Study Completion Date

P.O. Box 571830 | Murray, UT 84157-1830 U.S.A. • 6280 South Redwood Road | Salt Lake City, UT 84123-6600 U.S.A.  
 www.nelsonlabs.com • Telephone 801 290 7500 • Fax 801 290 7598 • sales@nelsonlabs.com

Page 1 of 2

These results relate only to the test article listed in this report. Reports may not be reproduced except in their entirety. Subject to NLS terms and conditions at www.nelsonlabs.com.



Study Number   
 Sodium Chloride (NaCl) Aerosol Test Final Report

Test Article Number	Initial Airflow Resistance (mm H <sub>2</sub> O)	Particle Penetration (%)	Filtration Efficiency (%)
14	7.5	1.44	98.56
15	7.4	1.10	98.90
16	7.2	1.37	98.63
17	7.2	1.23	98.77
18	6.9	1.36	98.64
19	7.7	1.26	98.74
20	7.3	1.25	98.75

**Acceptance Criteria:** The filter tester must pass the "Tester Set Up" procedure. The airflow resistance and particle penetration of the reference material must be within the limits set by the manufacturer.

**Filter Test Procedure:** Prior to testing, respirators were taken out of their packaging and placed in an environment of 85 ± 5% relative humidity (RH) and 38 ± 2.5°C for 25 ± 1 hrs.

The filter tester used in this procedure was a TSI® CERTITEST® Model 6130 Automated Filter Tester that is capable of efficiency measurements of up to 99.999%. It produces a particle size distribution with a count median diameter of 0.075 ± 0.020 µm and a geometric standard deviation not exceeding 1.86 µm. The mass median diameter is approximately 0.26 µm, which is generally accepted as the most penetrating aerosol size. The reservoir was filled with a 2% NaCl solution and the instrument allowed a minimum warm-up time of 30 min. The main regulator pressure was set to 75 ± 5 pounds per square inch (psi). The filter holder regulator pressure was set to approximately 35 pounds psi. The NaCl aerosol generator pressure was set to approximately 30 psi and the make-up airflow rate was set to approximately 70 liters per minute (L/min).

The neutralized NaCl test aerosol was verified to be at 25 ± 5°C and 30 ± 10% RH by the acceptance of the manufacturer's reference material. The NaCl concentration of the test aerosol was determined in mg/m<sup>3</sup> by a gravimetric method prior to the load test assessment.

An entire respirator was mounted on a test fixture, placed into the test article holder, and the NaCl aerosol passed through the outside surface of the test article at a continuous airflow rate of 85 ± 4 L/min. In accordance with NIOSH policy, three respirators were challenged until 200 ± 5 mg of NaCl had contacted the filter. Based upon the load pattern of NIOSH Type 2, the initial penetration reading of the remaining 17 filters was recorded.

P.O. Box 571830 | Murray, UT 84157-1830 U.S.A. • 6280 South Redwood Road | Salt Lake City, UT 84123-6600 U.S.A.  
 www.nelsonlabs.com • Telephone 801 290 7500 • Fax 801 290 7598 • sales@nelsonlabs.com

Page 2 of 2

# KN95 Protective Mask



## Test Report

(2020) WSZ FHL

Product Name Protective Mask against haze

Applicant \_\_\_\_\_

Manufacturer \_\_\_\_\_

Test Type Entrusted inspection

Jiangsu Guojian Testing Technology Co.,Ltd  
检验专用章



## Test Report

[2020] WSZ FHL Page 1 of 2

Product name	Protective Mask against haze	Specification	Willow leaf L
Applicant/Add Tel	Brand		
Manufacturer/Add Tel			
Sample grade	KN95	Sample number	
Sample quantity	26 PCS	Receiving date of sample	2020.02.28
Test type	Entrusted inspection	Article number/Batch number/Style number	Batch number: 200203
Test date	2020/02/28-2020/03/02	Test site	Testing room
Sample state	Meet the requirements		
Test standard(s)	GB 2626-2006 (Respiratory protective equipment -Non-powered purifying particle respirator)		
Test items	Filter efficiency, Respiratory resistance		
Test conclusion	The sample upon examination, the projects meet the requirements of the Q/0285HNS 039 standard, the specific test results see page 2.		
Note	The Applicant requires the items inspected to be determined in accordance with Q/0285HNS 039 (protective mask). The sample information are provided by Applicant. This report is only responsible for samples received.		



Approver: 陈永 Reviewer: 万恒 Main surveyor: 杨莹



## Test Result

[2020] WSZ FHL Page 2 of 2

Number	Test item	Unit	Technical requirements	Test Result		Single Item decision	
				Test Result	Test Result		
1	Filter efficiency	—	≥95.0%	Untreated	98.9%	Qualified	
					98.8%		
					99.1%		
					99.0%		
					99.1%		
					99.2%		
				Pretreatment	98.7%		
					98.9%		
					99.2%		
					99.0%		
					98.7%		
					98.8%		
2	Respiratory resistance	Pa	≤175	Untreated	103.2	Qualified	
				Pretreatment	98.7		
				99.6			
	Expiratory resistance	Pa	≤145	Untreated	83.1		Qualified
				Pretreatment	80.2		
				81.2			
Note							



The End



PAGE 1 of 2

## CERTIFICATE OF FDA REGISTRATION

Certification No. [REDACTED]

Dear Official Correspondent:

This document provides notification of the registration number assigned to your establishment:

Establishment:

Address: [REDACTED]

Owner/Operator Number:

Registration Number: [REDACTED]

*Helen Chan*

General Manager  
CIC REGISTRATION LLC  
Email: [cc-086@icloud.com](mailto:cc-086@icloud.com)  
Web: <http://www.ctc-086.com/>



Validity: Nov.7, 2019 – Dec.31, 2020

Conclusion:

This certificate makes no other representations or warranties, nor does it make any representations and warranties to any person or entity other than the named certificate holder. CIC assumes no liability to any person or entity in connection with the foregoing. The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. CIC is not affiliated with the U.S. Food and Drug Administration.

FDA Official Website: <https://www.fda.gov/>



PAGE 2 of 2

## ANNEX TO CERTIFICATION

Listing Number	Premarket Submission Number	Product Code(s)	Device Name(s) / Proprietary Name	Activities
	Exempt	LRJ	Disinfectants, medical devices / Iodine swab	Manufacturer
	Exempt	KGX	Tape and bandage, adhesive / Medical Tape, Adhesive Bandage, Wound Plaster	Manufacturer
	Exempt	JES	FLOSS, DENTAL / Floss pick	Manufacturer
	Enforcement Discretion	OBR	Monofilament, over-the-counter / Face Mask	Manufacturer
	Exempt	FQM	BANDAGE, ELASTIC / Elastic Bandage, Gauze Bandage	Manufacturer
	Exempt	EMD	PACK, HOT OR COLD, DISPENSABLE / Warm Pad, Cold Pad, Cool Pad, Fever Relief Pad, Cool gel, Freeze gel	Manufacturer
	Enforcement Discretion	LRB	Pad, alcohol, device disinfectant / Alcohol Wipes, Alcohol prep pad	Manufacturer
	Exempt	KYF	APPLICATOR, ABSORBENT TYPED, NON-STERILE / Cotton Swab	Manufacturer
	Exempt	NRC	Dressing, wound, hydrophilic / Wound Dressing, I.V. Adhesive Dressing	Manufacturer
	Exempt	NAB	Gauze / sponge, nonresorbable for external use / Gauze Pad, Underpad, Cotton balls	Manufacturer
	Enforcement Discretion	OIO	First aid kit without drug / First aid kit	Manufacturer


THIS ANNEX IS ONLY VALID IF ATTACHED TO THE CERTIFICATION MENTIONED ABOVE.



## Certificate

The Certification Body of  
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization



has established and applies a quality management system for medical devices for the following scope:

**Manufacture and Distribution of Sterile Wound Plasters, Sterile Wound Dressings, Medical Adhesive Tapes**

Proof has been furnished that the requirements specified in

**EN ISO 13485:2016**

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2018-08-10

Certificate Registration No.: [REDACTED]

An audit was performed. Report No.: [REDACTED]

This Certificate is valid until: 2021-09-07

Certification Body



Date 2018-08-10



TÜV Rheinland LGA Products GmbH · Tillystraße 2 · 90431 Nürnberg  
Tel.: +49 221 806-1371 Fax: +49 221 806-3935 e-mail: cert-safety@de.tuv.com http://www.tuv.com/safety

## **Contact Information**

### **Phone Number**

**+1 (312) 450-5043**

### **Email**

**anatoliy@safedirectms.com**

### **Office Address**

**1100 - 33 Yonge Street**

**Toronto, ON**

**M5E 0A9**