Swift Basics AK-PK Slide Preparation Kit

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Welcome to the AK-PK Slide Preparation Kit Instruction! Whether you are a seasoned professional or a beginner hobbyist, Swift has the perfect microscope for you.

This AK-PK kit instruction manual will introduce you the components of slides and guide you how to prepare it. For more information about the AK-PK slides kit (product manuals, or customer service requests), please visit www.swiftoptical.com.

WARNING! This AK-PK set contains functional sharp needles, and sharp edges on scalpel and slicer. This kit is not to be used by children under 8 years old, and always should be used under the direct supervision of an adult. Keep small children away from experiments. Store out of reach of small children. Read the instructions before use and keep them for future reference. Inaccurate use of the kit can result in injury and harm to health. The preparation area should be kept clear of any obstructions, should be well lit and fresh, and near a clean water supply.

I. Product Parts

1	Blank Glass Slides	x12
2	Mini Slicer	x 1
3	Petri Dish	x 1
4	Hatchery (Micro Incubator)	x 1
5	Collecting Vials	x 4
6	Prepared Slide-Bee Wing	x 1
7	3X Magnifier	x 1
8	Blank Labels	x14
	Round Slide Covers	x12
	Cover Glasses	x12
9	Prepared Slide – Textile	x 1
10	Pair of Tweezers	x 1
11	Stirring Rod	x 1
12	Prepared Slide-Animal Fiber	x 1
13	Scalpel	x 1
14	Dissecting Needle	x 1
15	Pipette	x 1

II. Item Descriptions:

1. **Blank glass slides**– These are the clear slides that you will place prepared samples on for examination under your microscope.
2. **Mini Slicer**–It is used for scraping off bits of material for testing and to push down on soft samples to mash them flat.
3. **Petri dish** –This is a round, flat dish with a clear cover. It is used to grow and observe samples such as mold
4. **Hatchery (Micro Incubator)** –This is a place where the hatching of insects.
5. **Collecting vials** –These are little plastic bottles with tight-fitting lids. They are used to carry your samples from where you collected to where you have your microscope set up.
6. **Prepared Slide-Bee Wing** –This is a slide that have bee wing sample on then that have been prepared by professionals for you to examine.
7. **Magnifier** –This is useful for taking a close look at a sample before you examine it under the high-power magnification of your microscope. (Do not view sun with the magnifier because it will make your eyes injured)
8. **Blank labels**–These are little pieces of paper with sticky backs. You can write on them and stick them on your slides to record information such as when the sample was prepared
- Round Slide Covers**–These are little circles made of thin, clear plastic. They are used to cover very small samples on a slide when they are clean and dry they stick to the glass slide with a static electricity charge.
- Cover Glasses**–They are used to cover very small samples on a slide.
9. **Prepared Slide-Textile**–This is a slide that have textile samples on.
10. **Pair of Tweezers**–The tweezers are like little pinchers. They are used to pick up small samples and to handle samples that you don't want to touch with your hands.
11. **Stirring rod** – Use this rod to mix liquids until they are well blended. An example is when you mix salt in with water.
12. **Prepared Slide-Animal Fiber**–This is a slide that have animal fiber sample on.
13. **Scalpel** – A scalpel is a sharp blade that is used to cut very thin pieces of material so you can look at them with your microscope.

14. Dissecting Needle –A metal stick used to extract or separate sample material for viewing on a slide

15. Pipette–This is a plastic device that you can dip into a liquid to transfer a drop or two to a slide for examination.

III. Preparing a Slide:

1. Making a Wet Mount:

- Wipe the blank slide and cover glass clean.
- Thin the sample with a mini slicer. (Be extremely careful and operated under adult only).
- Then pick it up with tweezers, and put it on the center part of the glass slide.
- Next, add one drop of water on the sample with a dissecting needle. If the sample is clear, add one drop of methylene blue or eosin solutions. (Be extremely careful and operated under adult only)
- Gently put the cover glass on it, take care not to let any air bubbles in it.
- Remove any excess water or dyeing solutions with blotting paper.
- Then, it is ready for observation.

2. Making a Permanent Mount:

- Wipe the blank slide and cover glass clean.
- Continue as above but before covering the slide with the cover glass, add few drops of gum media (Canada balsam or transparent adhesive glue with a dissecting needle to the slide).
- Put down on the cover glass with tweezers or fix it in place, and leave it to dry for about a day.

NOTE: Remember to wash your hands immediately after doing the preparations and dispose the dyeing solutions down the drain not into a sink.

IV. Warranty Information and Customer Support

Our microscopes are manufactured to meet ISO 9001 standards. Swift warranties are as follows:

- Five (5) Year Warranty for Microscopes: Microscopes come with a five (5) year warranty against manufacturing defects. Does not cover normal wear, routine maintenance, add-on accessories, damage resulting from repair by unauthorized parties, accident, alteration, shipping, misuse or abuse is not covered.
- One (1) Year Warranty for Electrical and Video components. Does not cover light bulbs, batteries, fuses, or electrical cords.

All warranties start from the original date of purchase. Swift provides the repair or replacement of warranted parts for free, including labor, during the warranty period. Proof of original purchase is required. Buyers are responsible for shipping to and from our warehouse for warranty services. The warranty does not cover damages resulting from normal wear and tear, abuse, or unauthorized repairs.

Please visit www.swiftoptical.com for online download of instruction manuals and relevant software.

*For customers living outside the United States, Swift Optical Instruments, Inc. will provide standard warranty service. Both inbound and outbound shipping costs (including duties and taxes) is the responsibility of the consumer.

For more information or to submit a repair request, please contact our customer support department:

North America

Tel: (+1) 877-967-9438, option 1 (USA)

Email: customersupport@swiftoptical.com

Europe

Tel: (+44) 208-638-8819 (UK)

Email: customersupport@swiftoptical.com

*In order to serve Europe customers efficiently, we recommend to call the hotline during working hours 8 AM – 4 PM (London Time).

Disclaimer: We are constantly working to improve our instruments and to adapt them in response to customer feedback. These improvements occasionally involve small modifications to the mechanical structure and optical design of our microscopes. Therefore, some descriptions, illustrations, and specifications in this instruction manual may vary slightly from the microscope you receive.

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