Swift EP Series Eyepiece Camera Quick Start Guide

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Thank you for purchasing this EP series eyepiece camera. This camera is a delicate instrument and should be handled with care. When it not in use, please store it in a safe place with the dust cap attached. Do not let the camera get in contact with water or expose to extreme temperatures (below 0 or above 50 degrees Celsius).

Should you have any questions, please do not hesitate to contact our customer service and you will get full support.

I. Product Components

Standard Configuration

Eyepiece camera

Application CD

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Optional Parts

Following adapters are not included in the standard package. If you have any needs, please contact our customer service.

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II. Install Swift EasyView Application in Your Computer


1. Please run application CD and there are two Directories Windows and Mac.

For windows users: Please browse the CD and open the file Swift EasyView under directory Windows, choose Swift EasyView x64 (for 64-bit Windows system) or Swift EasyView x86 (for 32-bit Windows system) according to your computer model. Double click the Swift EasyView file and be ready to install it accordingly.

For MAC users: Please browse the CD and open the Directory Mac to find the Swift EasyView file. Double click the .dmg file to install the application and then you can find Swift EasyView in your Mac application list.

Note: If there's no CD-drive on your computer, please visit www.swiftoptical.com to download the software from support menu.

III. Use the Eyepiece Camera on the Microscope Eyepiece Tube


1. Install the Swift EasyView application from the CD included or www.swiftoptical.com

downloaded.

2. Remove the dust cap from the camera.
3. Choose the eyepiece tube adapter according to your microscope eyepiece tube size (If your microscope eyepiece tube is in 30mm or 30.5mm internal diameter, you need an additional $\Phi 23.2\text{mm}-\Phi 30\text{mm}$ or $\Phi 23.2\text{mm}-\Phi 30.5\text{mm}$ eyepiece tube adapter. They are not included in standard package. Please contact our customer support to purchase them)
4. Attach the eyepiece camera (or camera-adapter set) to your microscope's eyepiece tube.
5. Connect the one end of Mini USB cable into the Mini USB2.0 port on the top of camera and plug the other end into the computer USB2.0 port.
6. Open the application. If necessary, please refer to help manual under software help menu for the operation of the Swift EasyView.

IV. Use the Camera on a Trinocular Microscope with Eyepiece Adapter

1. Install the Swift EasyView application from the CD included or www.swiftoptical.com downloaded.
2. Remove the dust cap from the camera.

3. Choose the eyepiece tube adapter according to your microscope eyepiece tube size (If your microscope eyepiece tube is in 30mm or 30.5mm internal diameter, you need an additional $\Phi 23.2\text{mm}-\Phi 30\text{mm}$ or $\Phi 23.2\text{mm}-\Phi 30.5\text{mm}$ eyepiece tube adapter. They are not included in standard package. Please contact our customer support to purchase them
4. Attach the eyepiece camera (or camera-adapter) set to your microscope's trinocular eyepiece tube.
5. Connect the USB cable onto the camera and then plug it into the computer.
6. Open the application. If necessary, please refer to help manual under software help menu for the operation of the Swift EasyView.

V. Troubleshooting

1. Previewing fails, blank screen.
 - 1) Check whether there is source of light or not, and whether the exposure value of video attribute is normal.
 - 2) Check whether video device is connected, and video device driver is installed in Device Manager.
 - 3) Check the dialog box of "select video device" to make sure the correct target video device is selected.
 - 4) Close the software and restart it.
 - 5) Close the software, pull the cam out and plug it in again, then restart the software.
 - 6) Close the software, restart computer and the software.

Please operate in line with the steps above, till you can preview normally. But if it still cannot work after you have followed all of the six steps, please contact our customer service.

VI. Warranty information and customer support

Our microscopes are manufactured to meet ISO 9001 standards. Swift warranties are as

follows:

- Five (5) Year Warranty for Microscopes: Microscopes come with a five (5) year warranty against manufacturing defects. Does not cover normal wear, routine maintenance, add-on accessories, damage resulting from repair by unauthorized parties, accident, alteration, shipping, misuse or abuse is not covered.
- One (1) Year Warranty for Electrical and Video components. Does not cover light bulbs, batteries, fuses, or electrical cords.

All warranties start from the original date of purchase. Swift provides the repair or replacement of warranted parts for free, including labor, during the warranty period. Proof of original purchase is required. Buyers are responsible for shipping to and from our warehouse for warranty services. The warranty does not cover damages resulting from normal wear and tear, abuse, or unauthorized repairs.

Please visit www.swiftoptical.com for online download of instruction manuals and relevant software (Manuals are downloadable in the product details page under “Swift Basics” products and Software are downloadable under “Support” menu).

*For customers living outside the United States, Swift Optical Instruments, Inc. will provide standard warranty service. Both inbound and outbound shipping costs (including duties and taxes) is the responsibility of the consumer.

For more information or to submit a repair request, please contact our Customer Support department:

North America

Tel: (+1) 877-967-9438, option 1 (USA)

Email: customersupport@swiftoptical.com

Europe

Tel: (+44) 208-638-8819 (UK)

Email: customersupport@swiftoptical.com

*In order to serve Europe customers efficiently, we recommend to call the hotline during working hours 8 AM – 4 PM (London Time).

Disclaimer: We are constantly working to improve our instruments and to adapt them in response to customer feedback. These improvements occasionally involve small modifications to the mechanical structure and optical design of our microscopes. Therefore, some descriptions, illustrations, and specifications in this instruction manual may vary slightly from the microscope you receive.

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