# Thermo® Dewars

## Cryogenic Transfer Vessels

Operation Manuals and Parts List LT509X34 (7000803) Rev. 7

Visit us online to register your warranty  
[www.thermoscientific.com/labwarranty](http://www.thermoscientific.com/labwarranty)

**Thermo**  
SCIENTIFIC

Catalog numbers	Model Numbers	Volume
TY509X1	803	5 liter
TY509X2	804	10 liter
TY509X3	805	20 liter
TY509X4	806	30 liter

**MANUAL NUMBER LT509X34 (7000803)**

7	ECNM-000831	02/24/21	Removed Warranty caution from preface page ii, added Safety Precautions in Preface page iii, WEEE Compliance in page 7-1 and removed EUT statements in page 4-1. Updated Specifications for Environmental conditions in page 4-1.	TFI
6	41238	01/27/17	Combined languages into one file	slb
5	40578	2/18/16	Matched symbols and labels to unit	ccs
4	30559	11/13/13	Added environmental conditions to pg 4-1	ccs
3	28501	4/13/12	Removed CE & declaration of conformity reference on pg 4-1	ccs
2	27797	9/21/11	Added new manufacturer icon and euro rep.	ccs
1	27457/SI-10618	7/21/11	Added CE reference number - pg 4-1	ccs
0	--	5/13/10	Transfer to Marietta (was LT509x34 10/22/09)	ccs
<b>REV</b>	<b>ECR/ECN</b>	<b>DATE</b>	<b>DESCRIPTION</b>	<b>By</b>

**Important** Read this instruction manual. Failure to read, understand and follow instructions in this manual may result in damage to the unit, injury to operating personnel, and poor equipment performance. ▲

**Caution** All internal adjustments and maintenance must be performed by qualified service personnel. ▲

**Warning** Nitrogen gas can cause suffocation without warning. Store and use liquid nitrogen only in a well-ventilated place. As the liquid evaporates, the resulting gas displaces the normal air in the area. (The cloudy vapor that appears when liquid nitrogen is exposed to the air is condensed moisture, not the gas itself. The issuing gas is invisible.) In closed areas, excessive amounts of nitrogen gas reduces the concentration of oxygen and can result in asphyxiation. Because nitrogen gas is colorless, odorless, and tasteless, it cannot be detected by the human senses. Breathing an atmosphere that contains less than 19.5% oxygen can cause dizziness and quickly result in unconsciousness and death. Therefore, the use of oxygen monitoring equipment is strongly recommended. ▲

**Warning** Contact of liquid nitrogen with the skin or eyes may cause serious (freezing) injury. ▲

**Warning** Do not roll. Always keep liquid nitrogen vessel in an upright position. Exercise caution to prevent spilling and splashing. ▲

Material in this manual is for information purposes only. The contents and the product it describes are subject to change without notice. Thermo Fisher Scientific makes no representations or warranties with respect to this manual. In no event shall Thermo be held liable for any damages, direct or incidental, arising out of or related to the use of this manual.

**Intended Use:**

The Thermo Fisher Scientific cryogenic transfer vessels are intended to store biological samples at cryogenic temperatures. These products are comprised of a vessel and a lid, with or without a level monitor.

**Unintended Use:**

These products are not intended for use within the patient environment.

When this manual is translated into other languages, the US English version is binding.

© 2021 Thermo Fisher Scientific. All rights reserved. Printed in U.S.A.

Important operating and/or maintenance instructions. Read the accompanying text carefully.

Potential electrical hazards. Only qualified persons should perform procedures associated with this symbol.

Extreme temperature hazards. Only qualified persons should perform procedures associated with this symbol.

Potential biological hazards. Proper protective equipment and procedures must be used.

This symbol when used alone indicates important operating instructions which reduce the risk of injury or poor performance of the unit.

**CAUTION:** This symbol, in the context of a CAUTION, indicates a potentially hazardous situation which if not avoided could result in minor to moderate injury or damage to the equipment. ▲

**CAUTION:** This indicates a situation which may result in property damage. ▲

Before installing, using or maintaining this product, please be sure to read the manual and product warning labels carefully. Failure to follow these instructions may cause the product to malfunction, which could result in injury or damage.

**WARNING:** Unauthorized repair of your freezer will invalidate your warranty. Contact Technical Service at 1-800-438-4851 for additional information. ▲

## Do You Need Information or Assistance on Thermo Scientific Products?

If you do, please contact us 8:00 a.m. to 6:00 p.m. (Eastern Time) at:

1-740-373-4763

1-800-438-4851

1-877-213-8051

<http://www.thermofisher.com>

[service.led.marietta@thermofisher.com](mailto:service.led.marietta@thermofisher.com)

[www.unitylabservices.com](http://www.unitylabservices.com)

Direct

Toll Free, U.S. and Canada

FAX

Internet Worldwide Web Home Page

Tech Support Email Address

Certified Service Web Page

Our **Sales Support** staff can provide information on pricing and give you quotations. We can take your order and provide delivery information on major equipment items or make arrangements to have your local sales representative contact you. Our products are listed on the Internet and we can be contacted through our Internet home page.

Our **Service Support** staff can supply technical information about proper setup, operation or troubleshooting of your equipment. We can fill your needs for spare or replacement parts or provide you with on-site service. We can also provide you with a quotation on our Extended Warranty for your Thermo Scientific products.

Whatever Thermo Scientific products you need or use, we will be happy to discuss your applications. If you are experiencing technical problems, working together, we will help you locate the problem and, chances are, correct it yourself...over the telephone without a service call.

When more extensive service is necessary, we will assist you with direct factory trained technicians or a qualified service organization for on-the-spot repair. If your service need is covered by the warranty, we will arrange for the unit to be repaired at our expense and to your satisfaction.

Regardless of your needs, our professional telephone technicians are available to assist you Monday through Friday from 8:00 a.m. to 6:00 p.m. Eastern Time. Please contact us by telephone or fax. If you wish to write, our mailing address is:

Thermo Fisher Scientific (Asheville) LLC  
401 Millcreek Road, Box 649  
Marietta, OH 45750

International customers, please contact your local Thermo Scientific distributor.

# Table of Contents

<b>Section 1</b>	<b>Safety Information</b> .....	<b>1-1</b>
<b>Section 2</b>	<b>General Description</b> .....	<b>2-1</b>
<b>Section 3</b>	<b>Operation</b> .....	<b>3-1</b>
	Filling Instructions .....	3-1
	Withdrawal Device .....	3-2
	Installation .....	3-2
	Operation .....	3-4
	Removal of Liquid Discharge Device .....	3-4
	Liquid Withdrawal .....	3-5
	Measuring Liquid Nitrogen Quantity .....	3-5
	Filling Instructions (Thermo 8127) .....	3-5
	Low Pressure LN2 Supply Filling Method (8127) .....	3-6
	Funnel Filling Method .....	3-6
	Liquid Withdrawal for Model 8127 .....	3-7
<b>Section 4</b>	<b>Specifications</b> .....	<b>4-1</b>
<b>Section 5</b>	<b>Valve Maintenance and Repair</b> .....	<b>5-1</b>
	Ordering Procedures .....	5-1
<b>Section 6</b>	<b>Warranty Information</b> .....	<b>6-1</b>
<b>Section 7</b>	<b>WEEE Compliance</b> .....	<b>7-1</b>

# Section 1 Safety Information

Your Thermo® Dewar has been designed with function, reliability, and safety in mind. It is the user's responsibility to install it in conformance with local electrical codes. For safe operation, pay attention to the Warnings, Cautions and Notes throughout the manual.

## Warnings

Liquid nitrogen is extremely cold; it boils at  $-196^{\circ}\text{C}$ . To avoid injury caused by frostbite:

1. Use extreme care whenever handling liquid nitrogen, liquid nitrogen storage or transfer vessels, or any objects which have come in contact with liquid nitrogen.
2. Leave no areas of skin exposed.
3. Use in a well ventilated area.
4. Wear adequate PPE where necessary: face shield, cryogenic gloves and apron.
5. Do not wear pants with cuffs.
6. Never overfill liquid nitrogen vessels.
7. Always keep liquid nitrogen vessel in an upright position.
8. Do not tightly seal liquid nitrogen containers or prevent nitrogen gas from escaping.
9. Use extreme care to prevent spilling and splashing of liquid nitrogen during transfer.
10. Immediately remove any clothing or safety attire on which liquid nitrogen has been spilled.
11. Read all filling instructions carefully.
12. Release pressure and remove liquid nitrogen before working on a vessel.
13. Get immediate medical attention for any frostbite injuries caused by liquid nitrogen.

14. Use this product only in the way described in the product literature and in this manual. Before using it, verify that this product is suitable for its intended use. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

**Warning** The venting of nitrogen vapors will create a dilution of the air's oxygen concentration necessary to support life. Exposure to this diluted atmosphere can cause asphyxiation or even death. To avoid personal injury or death, DO NOT store or use liquid container in areas that have poor ventilation. Use in a well ventilated area. Place liquid containers outdoors or in a well ventilated area. Failure to comply with this warning may cause serious personal injury or even death. ▲

**Warning** Extreme temperature and pressure are associated with the Thermo 8127. To avoid serious personal injury, observe caution when removing parts or fittings. Wear eye protection and insulated, loose fittings gloves when removing parts or fittings. Failure to comply with this warning may result in serious personal injury. ▲

**Caution** DO NOT use regulators, valves, gauges, hoses, etc., that have been used in compressed air service. Failure to observe this warning could cause serious damage to your container and possible personal injury. ▲

## Section 2 General Description

Mechanical failure is the major cause of problems associated with liquid nitrogen storage vessels. The vessel necktube supports the full weight of the inner vessel and the liquid nitrogen it contains. A side or corner blow to the vessel causes the inner vessel to swing in a pendulum motion, damaging the neck. Any storage vessel which has been exposed to an accident, dropped, or lowered to hit on one corner, is more susceptible to failure. Model 8127 has a support (top and bottom), but mishandling may still cause problems. Carefully inspect your new Thermo vessel prior to use. Check for signs of damage which may have occurred during shipment. It is advisable to fill all new units with liquid nitrogen and watch the liquid nitrogen loss rate for a few days. If there are any problems, call Technical Services as soon as possible.

The on-line warranty registration address is listed on the front of this manual. The unit should be registered within 30 days to make warranty effective. This information must include the serial number located on the handle of the unit.

## Section 2

### General Description

## Section 3 Operation

**Warning** Liquid nitrogen is extremely cold; it boils at  $-196^{\circ}\text{C}$ . Extreme caution should be exercised when handling cryogenic liquids. Always wear proper safety attire when transferring liquid ni-trogen into or out of Thermo cylinders. Contact with a cryogenic liquid or cold gas may cause frostbite to unprotected areas of the body. Protect eyes and skin when transferring liquid. Stand clear of vent during filling. Loose fitting, insulated gloves and long sleeves are recommended for arm protection. Cuffless pants should be worn outside boots or overshoes to shed spilling liquid. Goggles or face shields should be worn if the possibility of splashing liquid exists. Failure to observe this warning may lead to severe burns or eye injuries. ▲

**Caution** Never overfill your Thermo dewar with liquid nitrogen. Filling the tank up to or above the bottom of the necktube may cause immediate or premature vacuum failure. ▲

### Filling Instructions

To avoid damage to your Thermo cryogenic vessel (Thermo<sup>®</sup> 5, 10, 20 and 30), which may result in premature vacuum loss, it is important that the following procedure be used during the addition of liquid nitrogen to a warm vessel and on subsequent additions.

1. Add only a small amount of liquid nitrogen (5-10 liters) to new or warm vessels.
2. Allow this small amount of liquid nitrogen to sit in the covered vessel for a minimum of 2 hours. This will limit stress caused by the sudden temperature change associated with adding liquid nitrogen to a warm vessel.
3. Fill your vessel to the desired level.
4. If you are filling your Thermo dewar from a pressurized source, make sure that the source tank is at a low pressure (22 PSI or below).
5. If the transfer hose is used for extracting liquid nitrogen from a pressurized liquid source, always use a phase separator at the end of the hose (see Parts List for transfer hoses and phase separators).
6. Remember to always wear proper safety attire over clothing, including: face shield, cryogenic gloves and apron.

## Withdrawal Device

The liquid discharge device provides liquid nitrogen for filling cold traps and smaller dewars without having to scoop out or pour liquid nitrogen from the vessel. The liquid discharge device utilizes the normal evaporation of the cryogenic liquid to pressurize a container by plugging the neck of the vessel and forming a hermetic seal. The pressurized vessel can then decant product through the liquid discharge device. Figure 1 shows a drawing of the liquid discharge device. The primary components of the liquid discharge device are listed below:

**Rubber Stopper:** Provides a seal against the walls of the necktube. Tightening the wingnut expands the rubber stopper.

**Pressure Gauge:** Indicates the pressure inside the container.

**Relief Valve:** Relieves pressure inside the container automatically when excessive pressure is reached. The relief valve is a safety device and cannot be adjusted.

**Vent Valve:** Provides means to vent product. The vent valve can be closed to pressurize the vessel or opened to depressurize the vessel.

**Discharge Valve:** Allows product to flow out of the vessel through the discharge device spout.

**Safety Cable:** Prevents discharge device from accidentally exiting completely from vessel.

## Installation

Perform the following steps to install the liquid discharge device on a vessel:

1. Attach the discharge device spout to the liquid discharge device (refer to Figure 1 for location). A transfer line can be used in place of the discharge device spout.
2. Verify that the rubber stopper is dry and free from grease and other contaminants. Check for nicks and gouges that may impair normal operation.
3. Verify that the cable is in good working condition.
4. Verify that the vent and discharge valves turn freely.
5. Loosen wingnut.
6. Lower discharge device into vessel. As the discharge device is lowered, the nitrogen in the vessel will boil until the dip tube has cooled down. The boiling action is normal. Continue lowering until the rubber stopper fits snugly.

### Installation (continued)

7. Tighten the wingnut until it contacts the washer. Tighten wingnut another 1½ to 2 turns (EXT RBR versions 4-5 turns).

**Caution** Tighten the wingnut finger-tight only, do not use any tools. Tightening the wingnut expands the rubber stopper against the necktube of the vessel. Excessive tightening may damage the necktube. ▲

8. Clip the safety cable to the vessel handle.
9. Close the vent valve and the discharge valve. Pressure will slowly begin to build. Several hours will be required to achieve sufficient pressure to operate the liquid discharge device. Pressure buildup can be accelerated by gently tipping the vessel at 45° angle.

**Caution** Excessive pressure can rupture the vessel. ▲


Figure 3-1. Components

## Operation

Perform the following steps to withdraw liquid from the vessel using the liquid discharge device:

1. Turn the discharge valve fully open and then turn back the valve a half turn. Do not leave the valve in the fully open position.
2. Close the discharge valve to end liquid withdrawal.

**Warning** Cryogenic liquids are extremely cold and can cause severe burns similar to frostbite. See CGA Bulletin No.P-12 Safe Handling of Cryogenic Liquids. (This bulletin may be ordered from Compressed Gas Ass'n., 4221 Walney Road, 5th floor, Chantilly, VA 20151.) ▲

**Warning** Nitrogen vapors in air may dilute the concentration of oxygen necessary to support or sustain life. Exposure to such oxygen deficient atmosphere can lead to unconsciousness and serious injury, including death. ▲

**Warning** Modifying safety valves on discharge devices causes potentially hazardous conditions. ▲

## Removal of Liquid Discharge Device

Perform the following steps to remove the liquid discharge device from the vessel:

1. Open the vent valve to remove all pressure.
2. Verify that the pressure gauge reads 0.

**Warning** Vessel must be completely vented before removing the discharge device. Any attempt to remove the discharge device before the vessel is completely vented may result in serious personal injury. ▲

3. Release wingnut.
4. Rock discharge device from side to side to free stopper.
5. Pull discharge device from vessel.
6. When stopper has passed rim of vessel, unclip safety cable.
7. Store liquid discharge device in a clean and dry area.
8. Cover vessel.

## Measuring Liquid Nitrogen Quantity

The procedure below pertains to Thermo 5, 10, 20, and 30 vessels.

1. Use a wooden dipstick. Never use a hollow tube or plastic dipstick to measure liquid nitrogen.
2. The liquid level will be indicated by a frostline which develops when the dipstick is removed.

## Liquid Withdrawal

1. Liquid withdrawal for Thermo 5, 10, 20 and 30 is always performed by pouring or utilizing a withdrawal device (for withdrawal devices, see Parts List). The withdrawal device pressurizes to approximately 4 PSI; the pressure forces liquid up the withdrawal tube and out the valve.
2. Always wear proper safety attire.

**Warning** Liquid nitrogen is extremely cold; it boils at  $-196^{\circ}\text{C}$ . Extreme caution should be exercised when handling cryogenic liquids. Always wear proper safety attire when transferring liquid nitrogen into or out of Thermo cylinders. Contact with a cryogenic liquid or cold gas may cause frostbite to unprotected areas of the body. Protect eyes and skin when transferring liquid. Stand clear of vent during filling. Loose fitting, insulated gloves and long sleeves are recommended for arm protection. Cuffless pants should be worn outside boots or overshoes to shed spilling liquid. Goggles or face shields should be worn if the possibility of splashing liquid exists. Failure to observe this warning may lead to severe burns or eye injuries. ▲

## Filling Instructions (Thermo 8127)

To avoid damage to your dewar, it is important that only the two following methods be used to fill your vessel. Failure to follow each step may result in damage which is not covered under the warranty. Read all instructions carefully. Should questions arise, contact Technical Services.

## **Funnel Filling Method**

The Thermo 8127 cylinders can be filled by removing the brass plug on the top center of the tank and inserting a funnel through the hole. Liquid nitrogen can then be poured directly into the cylinder through this hole. Be sure to stand clear of the vent during filling and always wear proper safety attire while handling liquid nitrogen.

1. Open vent valve completely, releasing any pressure built up inside the cylinder.
2. Remove brass plug located on the top center of the cylinder.
3. Insert funnel into hole.
4. Pour liquid nitrogen into cylinder until level gauge reads 7/8 full or until liquid nitrogen begins spitting from the vent valve.
5. Reinsert and tighten the brass plug.
6. Close vent valve completely.

## **Low Pressure LN2 Supply Filling Method (8127)**

The Thermo 8127 liquid nitrogen cylinders can be filled from a pressurized source of liquid nitrogen by attaching a transfer hose to the liquid withdrawal valve on the Thermo cylinders. The liquid nitrogen source pressure must not exceed 45 PSI. Read all instructions carefully before filling.

1. Attach transfer hose from liquid nitrogen source to the liquid withdrawal valve on the Thermo cylinder.
2. Open withdrawal valve on Thermo cylinder completely.
3. Open withdrawal valve on liquid nitrogen source. Liquid nitrogen source pressure must not exceed 45 PSI, since damage to Thermo cylinder gauges will occur. Optimum feed pressure is 35 PSI.
4. Open vent valve on Thermo cylinder until the pressure gauge on Thermo cylinder reads 22 PSI.
5. Continue to fill until Thermo cylinder weight is 180 pounds for the 8127. If the scale method is not possible, fill until liquid nitrogen begins spurting from vent valve.
6. Shut liquid nitrogen source valve completely.
7. Shut liquid nitrogen withdrawal valve on Thermo cylinder completely.
8. Carefully remove transfer hose from the Thermo cylinder. Some liquid nitrogen will remain in the hose under pressure after filling.
9. Close vent valve completely.

**Caution** Overfilling cylinders may result in damage to the level and pressure gauges. If overfilling should occur, remove excess liquid nitrogen by opening the liquid nitrogen withdrawal valve immediately. ▲

## Liquid Withdrawal for Model 8127

Thermo 8127 are to be used only for low pressure liquid withdrawal. The primary relief valve is factory set at 22 PSI. The secondary safety relief is set at approximately 35 PSI. Transferring liquid at higher pressures increases the flash-off rate of the liquid and adds to the risk of sparking.

To transfer liquid, attach the transfer hose or withdrawal spout to the liquid connections (see accessory sheet). Slowly open the liquid valve to flow the liquid. The liquid will vaporize at first until the transfer line or withdrawal valve cools down. If using a transfer hose to extract liquid from the Thermo 8127 into an open dewar, a phase separator is recommended on the end of the transfer line.

Transfer pressure should be kept to a minimum. The normal evaporation of the liquid will usually maintain enough pressure for transferring.

**Caution** If the liquid or vent valves are not closed, the container can become contaminated after being emptied. ▲

**Section 3**  
Operation

## Section 4 Specifications

Vessel Specs for Thermo® Dewars and Cylinders				
	TY509X1 (803)	TY509X2 (803)	TY509X3 (805)	TY509X4 (806)
Catalog Number (Model Number)	TY509X1 (803)	TY509X2 (803)	TY509X3 (805)	TY509X4 (806)
Diameter - inches (cm)	8.8 (22.4)	10.3 (26.2)	14.5 (36.8)	17.0 (43.2)
Height - inches (cm)	18.2 (46.2)	21.5 (54.6)	24.7 (62.7)	24.0 (61.2)
Neck Diameter - inches (cm)	2.2 (5.6)	2.2 (5.6)	2.0 (5.0)	2.5 (6.4)
Full Weight - pounds (kg)	17.0 (7.7)	31.0 (14.1)	55.0 (25)	82.0 (37.3)
Empty Weight - pounds (kg)	8.0 (3.6)	12.0 (5.5)	19.0 (8.6)	27.0 (11.6)
Static Evaporation Rate (liters/day)*	0.15	0.18	0.18	0.22
Static Holding Time (days)	33.3	55.6	111.0	145.0
Liquid Nitrogen Capacity (liters)	5.0	10.0	21.0	32.0

*\* Actual working time may vary due to current atmospheric conditions, container history, manufacturing tolerances, and any individual patterns of use.*

### Environmental Conditions

Operating temperature: -29°C to +60°C

10% to 95% relative humidity, non-condensing.

Altitude limit: 2,000 meters.

Storage: -29°C to +60°C; 10% to 95% relative humidity.

Intended for indoor use only.

**Section 4**  
Specifications

## Section 5 Valve Maintenance and Repair

When a defective valve is suspected, follow this procedure to repair it:

1. Void the tank of liquid product and release any pressure that is in the container.
2. If the vent valve must be repaired, allow it to warm-up before disassembling.
3. Remove the valve handle screw, washer, retainer cap and spring assembly.
4. Remove the valve handle and teflon thrust washer.
5. Unscrew the bonnet to remove the stem and stem seal.
6. Pull out the body insert and plug assembly.
7. Clean the seat.
8. Replace parts as needed and reassemble in reverse order.

### Ordering Procedures

Refer to the Specification Plate for the complete model number, serial number, and series number when requesting service, replacement parts or in any correspondence concerning this unit. All parts listed herein may be ordered from the Thermo dealer from whom you purchased this unit, or can be obtained promptly from the factory. When service or replacement parts are needed, check first with your dealer. If the dealer cannot process your request, then contact our Technical Services Department. Prior to returning any materials, contact Technical Services for a “Return Materials Authorization” number (RMA). Material returned without an RMA number will be refused.

## **Section 5**

### Value Maintenance and Repair

## THERMO FISHER SCIENTIFIC STANDARD PRODUCT WARRANTY (LN<sub>2</sub> Vacuum)

The Warranty Period starts two weeks from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period.

During the first year, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor included. *LN<sub>2</sub> Vacuum Integrity is covered for five years.* Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo postage paid and replacement parts are shipped FOB destination.

**THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY.**  
Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed instruction manuals carefully detail equipment installation, operation and preventive maintenance.

If equipment service is required, please call your Technical Services Department at 1-800-438-4851 (USA and Canada) or 1-740-373-4763. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application. Outside the USA, contact your local distributor for warranty information.

Rev. 0 5/10

## **THERMO FISHER SCIENTIFIC INTERNATIONAL DEALER WARRANTY (LN2 VACUUM)**

The Warranty Period starts two months from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period. Dealers who stock our equipment are allowed an additional six months for delivery and installation, provided the warranty card is completed and returned to the Technical Services Department.

During the first two years, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor excluded. *LN2 Vacuum Integrity is covered for five years.* Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters, reagents, tubing, and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo postage paid and replacement parts are shipped FOB destination.

**THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY.** Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed operation manuals carefully detail equipment installation, operation and preventive maintenance.

Contact your local distributor for warranty information. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application.

Rev. 0 7/10

## Section 7 WEEE Compliance

### WEEE Compliance

**WEEE Compliance.** This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with the following symbol. Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, and this product should be disposed of or recycled through them. Further information on our compliance with these Directives, the recyclers in your country, and information on Thermo Scientific products which may assist the detection of substances subject to the RoHS Directive are available at [www.thermofisher.com/](http://www.thermofisher.com/)

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Spanish

**thermoscientific.com**

© 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Thermo Fisher Scientific (Asheville) LLC**  
401 Millcreek Road  
Marietta, Ohio 45750  
United States

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