

BLADE CLEANING SYSTEM



MICRO-SCIENTIFIC BLADE CLEANING SYSTEM

An efficient and safe method for cleaning tissue and other debris from the surfaces of diamond and other delicate blades prior to sterilization. Blades can be cleaned without scratching, dulling, chipping or breaking.

- Neutral pH detergent
- Phosphate free
- Gentle scrubbing pads
- Corrosion inhibitors
- Biodegradable solution
- Easy as 1-2-3

Contents:

- Detergent [Blue]
- Rinse Solution [Clear]
- Scrubbing Pads [x3]

ORDER NO.: M60070

M60070-84637B-0219



An efficient and safe method for cleaning tissue and other debris from the surfaces of diamond and other delicate blades prior to sterilization. Blades can be cleaned without scratching, dulling, chipping or breaking.

- **Neutral pH detergent**
- **Corrosion inhibitors**
- **Gentle scrubbing pads**
- **Easy as 1-2-3**

BLADE CLEANING SYSTEM



Order No.	System Contents	Physical Properties	
M60070	Detergent, Rinse Solution Scrubbing Pads (x3)	Form:	Liquid, & Scrubbing Pads (x3)
		Color:	Clear (Rinse) & Blue (Detergent)
		Odor:	None
		pH Concentration:	8.5

Uses			Temperature
<input type="checkbox"/> Presoak	<input type="checkbox"/> Automatic Washer		Usage: Room temperature
<input checked="" type="checkbox"/> Manual	<input type="checkbox"/> Automatic Endoscope Reprocessor		Storage: Cool, dry area
<input type="checkbox"/> Ultrasonic	<input type="checkbox"/> Cart Washer		
<input type="checkbox"/> Evacuators			

Product Features	Consumer Benefits
Neutral pH:	Product safe and compatible with all metals.
Environmentally-Friendly:	Biodegradable, phosphate-free solution. Recyclable plastic tray.
Corrosion Inhibitors:	Reduces likelihood of rust, tarnish or corrosion.
Scrubbing Pads:	Gentle on delicate blades. Blades can be cleaned without scratching, dulling, chipping or breaking.
Easy As 1-2-3:	Clean, rinse, and rinse. A sparkling clean knife, ready for sterilization.

Procedures	Directions for Use
Directions:	Insert the blade directly into the saturated cleaning pad located in the large compartment. Gently move the blade back and forth in a cutting fashion to remove surface debris. Repeat until clean. Repeat the above procedure in the upper left-hand rinse compartment to remove the cleaning solution. Repeat step two in the upper right-hand rinse compartment as a final rinse. Follow blade manufacturer's rinsing and sterilization instructions. Dispose of blade cleaning system daily.