

Transparent Male Catheter Model


Instructions

Please read these instructions carefully before using.

How to Use

(See drawing on the right)

- ① Hang irrigator(water tank) to the exclusive stand that comes with the model, and fix.
 - ② Connect the longer tube coming out of the bladder to the connector at the top of the irrigator. (When doing this, set the stand some distance away from the model so that the tube does not dangle between the stand and bladder.)
 - ③ Connect the shorter tube coming out of the bladder to the connector below the clamp.
 - ④ Fill irrigator(water tank) with approximately 500cc of water. Next, open the clamp and water will enter the bladder. (It will take some time before the water becomes full.)
 - ⑤ In the beginning it may be difficult to insert, so use Nelaton catheter(size Eg.8 etc.) and spread exclusive lubricating oil thoroughly on the tip. Confirm that the catheter has been inserted inside the bladder.
 - ⑥ Practice catheterization using balloon-tipped catheter.
- ※After inserting the catheter, if it does not drain(urinate), add a small amount of water (approximately 100cc) to raise water pressure.




【Caution Regarding Use of Lubricant】

- When inserting catheter, always **spread lubricant thoroughly**.
- It is not possible to insert when the catheter angle differs. Do not insert forcibly.
- When there is water on the tip of the catheter, wipe it off before applying lubricant.
- Cannot be used for the human body.

【Caution When Using】

- **The external genitalia and bladder are fixed. Never attempt to detach them.**
- **Never attempt to detach the tube attached to the bladder.**
- The transparent portion of this product is susceptible to scratches.
If hit against something or strongly rubbed, it may cause damage or blurs, so please be careful.
- When it becomes dirty, wipe with soft cloth using kitchen detergent or wipe with cloth only.
※Please do not use alcohol, benzine, a thinner, or polishing powder.
- Any stain in the external genitalia should be wiped gently with a soft cloth using kitchen detergent, and then, wiped with a cloth moistened with water.
Further, before storing, slightly cover with powder using a brush or cloth etc.
- It will break if dropped, so please handle with care.

- ① Prepare a tray etc. (something to catch the drainage), close the clamp, and pull out the tube from the connector.
- ② As shown in the illustration, tilt the transparent model sideways until the shorter tube coming out of the bladder comes directly below.
- ③ Water remaining inside the bladder can be drained through the shorter tube coming out of the bladder.
- ④ Water remaining inside the irrigator will be drained by opening the clamp.



【How to Store】

- Pull out the tube connected to the connector at the top of the irrigator.
- Pull out the tube connected to the connector below the clamp.
- Water remaining inside the bladder can be drained through the shorter tube coming out of the bladder by tilting the body.
- Wipe off water remaining inside the irrigator using a slender stick wrapped with cloth.
- Please understand that a small amount of waterdrops will remain inside the tubes.
- Fold the irrigator, and put inside a plastic bag.
- Put exclusive stand and model back into the hard case.
- Do not store under direct sunlight or where it is hot and humid.