

THE SCIENCE BEHIND AGERELIEF

COMPRESSION SHOULDER BRACE



COMPRESSION THERAPY

Enhances blood circulation, funneling nutrient-rich blood to rotator cuff tendons, expediting the healing process.



STABILIZATION

Maintains the shoulder's natural alignment, fostering an ideal healing environment.



SUPPORT

Secures the shoulder, minimizing movement to prevent further injury and aid tendon recovery.

TIPS FOR BEST RESULTS



Consistently wear for a minimum of 14 days



Ensure a snug fit, not too tight, not too loose



Combine with a healthy lifestyle for optimal recovery

Tips for Use

- 01 Wear a thin layer of clothing under the brace for comfort.
- 02 Start with wearing the brace for short periods to allow your body to adjust.
- 03 Do not over-tighten; the brace should feel supportive, not restrictive.

Cleaning and Care Instructions

- 01 Wash with mild detergent
- 02 Air dry away from direct heat
- 03 Do not bleach or iron



HOW TO WEAR THE AGERELIEF

COMPRESSION SHOULDER BRACE

4 EASY STEPS TO WEAR YOUR BRACE



STEP 01

Slip your arm through the brace, ensuring the supportive pad sits comfortably on your shoulder.

Adjust the strap around your bicep to a snug fit that does not restrict circulation.



STEP 02

Pull the chest strap around your back and secure it in front for even compression.



STEP 03

Make final adjustments to all straps to ensure a comfortable fit with adequate support.



STEP 04