

Full Body Vibration Heated Massage Pad



How to use:

Spread the massage pad flat on the sofa or bed and other flat places, the DC head of the power adapter and the DC head of the massage pad are firmly docking, and the power plug of the power adapter is firmly docking into the electrical socket. At this time, the number on the plate of the hand control is to show "OFF", indicating that the massage pad has been powered on normally, and you can enter the next operation:

Neck Controller Instructions:

Switch

Press the switch button, one key to start up, the whole machine enters automatic mode and presses the second time to shut down the whole machine.

Direction Key

The default one-minute forward and reverse is turned on at startup, press the button the first time, open the corotation function, and press again to turn on the inversion function.

Force Adjustment Key

The machine has a total of 3 levels of force adjustment, the machine starts running in default medium, and presses the button for the first time, for low force gear operation. Press the key again for high-force gear operation. Press again to resume mid-range operation.



Full Body Vibration Heated Massage Pad

Heating key

Turn on the heating function by default, and press the button the first time to turn off the heating function. Press this button again to start heating again.

Energy Saving Key

Energy-saving design, press this button to enter energy-saving mode.

Smart Mode Key

In the power-on state, this button is pressed 3 times to experience 3 different intelligent massage modes, and 4 times is switched to the first intelligent massage mode, and so on.

Pad Controllers Instructions:

1. Switch and Time: Pressing this button will display "10", massage pad to start working, at this time in automatic mode, and default 10 minutes timing. Press again to display "20", and start 20 minutes of timing. Then press again to display "30", and start 30 minutes timing. Finally, press again to display "OFF" The massage pad will stop working.

2. Mode: Press this button to select 9 different massage modes of whole body vibration massage, "F1"- "F8"- "AU" a total of 9 modes to select.

3. Force Adjustment: Press this button to increase or decrease the vibration massage force, a total of 9 gears of force adjustment, namely "S1"-"S9".

4. Speed: Press this button to speed up or slow down the vibration massage frequency, a total of 9 gears of speed and slowness adjustable, namely "L1" - "L9".

5. Positioning: Press this button to control the massage area by part, and "P1" -"P4" will be displayed, with a total of 4 different areas adjustable.

6. Heating: Press this button to turn on the heating function of the massage pad, display "H1"-"H3", a total of 3 gears are adjustable when the button is pressed to display "HO" can turn off the back heating function of the massage pad.

Waist Pillow Gasbag Traction Instructions:

Press the gasbag ball to inflate the waist pillow.

Press the deflation button, massage the waist pillow Gasbag starts to deflate.

(Adjustable according to the user's own comfort need)