



RM533  
RM533B Bariatric



RM535  
RM535B Bariatric



RM531LG

## Day Chairs

# **Warranty**

---

All products carry a 12 month warranty from the date of purchase against manufacturing defects.

Misuse or Damage is not covered.

This warranty does not affect your statutory rights.

---

# Assembly / Adjustment



RM Seating products come ready assembled or may require minimal assembly.

## **Assembly of Chair**

**Step 1:** Slide legs on to each side of seat base.

**Step 2:** Adjust seat to comfortable depth positioning.

**Step 3:** Tighten screws to secure the frame.

**Step 4:** Slide backrest onto locating poles.



## **Adjustment of the legs** (if applicable)

**Step 1:** Using your thumb, push down the Push-Pin and slide the leg up or down until the Push-Pin clicks back into the location hole.

*(In some cases there maybe a push-pin on either side of the leg.)*

**Step 2:** Ensure the Pin/s is completely through the hole, and test to ensure the two pieces do not come apart.

**Step 3:** Ensure all legs are adjusted to the same height.

# Maximum User Weight

<b>Product Code</b>	<b>Description</b>	<b>Weight Limit</b>
RM533	High Back Day Chair	200kg
RM533B	High Back Bariatric Day Chair	250kg
RM535	Low Back Day Chair	200kg
RM535B	Low Back Bariatric Day Chair	250kg
RM531LG	Foot Stool	n/a

It is important to never exceed this limit and always ensure that the product weight capacity is sufficient before using the product. Failure to do so could result in injury or damage the product.

---

## Care & Maintenance

All of the product components are manufactured from high quality materials. However, it is very important that you only clean your Day Chair using a non-abrasive cleaner or mild detergent with a soft cloth. Abrasive cleaners and/or abrasive cleaning pads will seriously damage your product beyond repair.

All chairs and footstools should be only be used for their intended purpose.

All chairs and footstools should be inspected regularly to ensure moving components (push-pins, locking screws) remain secure. If in doubt, do not use and immediately contact your supplier for support.

Your safety is of the utmost importance.