



Features

- Light Weight
- Easy to fold with single push-button
- Interchangeable legs from tips to wheels
- Maximum User Weight Capacity of 130 Kg (290lbs)

Safety Instructions

- A risk assessment should always be carried out to ensure the product is suitable.
- Ensure all components are properly and securely installed.
- If normal function of the castors is not observed, return the Walker to the place of purchase for adjustment or replacement of parts.
- Ensure the user has sufficient strength to safely operate and control the Walker.
- Do not exceed the maximum recommended working load.
- Avoid storing the Walker in direct sunlight; handles or parts exposed to sunlight will heat up and may cause injury to the skin.
- When using outdoors, use only on designated footpaths or walkways; do not use on roads, streets or highways.
- Do not use the Walker on any incline greater than 9 degrees, or on a slope with water, ice or snow covering.
- Use caution when traversing kerbs or obstacles as this may cause the Walker to tip over and cause bodily harm or damage to the Walker

Assembly Instruction

1. Remove all contents from carton: Frame and Legs
2. Open each side of the frame until you hear a click
3. Attach the 4 Legs

Finding the Right Height

Elbow should be at 30 Degrees when standing and holding onto the walker handles.

