490 mm.

10 mm.

20 mm.

156 mm.

TRIM GUIDES

TRIM LINES

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
Brand: Neutónica | Model: BNST-1 | Made in: China


QUICK START GUIDE.

Read the instructions and/or indications carefully.

Once you've completed your daily routine and you're in bed ready to go to sleep:

- 01** Join both Neutonica Bands on the side of your preference. The rear band should have the widest curve upwards and the smallest curve downwards, as in the following image:
- 

- 02** Place the Neutonica Band around your head and fasten it on the right side.
- 


- 03** Using the nose and the center of the back band as reference, open and fasten each side as necessary so that both parts are centered and fitted securely, but without pulling too tightly.
*It must not completely cover your ears.
- 

- 04** Close your eyes and the cushion will be on your nose and closed eyes, but it shouldn't be tight or feel uncomfortable.

Now, with your right hand behind your ear, you will feel the **on button**.

Press and hold for two seconds, you will hear an on sound and you will be ready to relax.

Remember that the Neutonica Band will turn off automatically once it has finished carrying out its purpose or manually by holding the button for two seconds.

490 mm.

10 mm.


ABOUT ITS USE.

- The initial placement is different on each person. Once you find the ideal fit for you, it is not necessary to open the Neutonica Band. You can remove it and put it on from above your head without having to open it each time.
- Each time you use it you will hear a different track. The different durations are created to simulate the natural time of sleep latency. If the audio finishes and you haven't been able to fall asleep shortly after it's finished, we recommend you remove the Neutonica Band, do an activity out of bed for no longer than 15 to 20 minutes, and afterwards, retry another cycle with the Neutonica Band. With repeated use, your brain will get used to sleeping before the Neutonica's Band cycle ends within a time frame ranging from 2 weeks to 2 or 3 months, depending on each person and the cause of their insomnia.
- Once it is on, you can briefly press the button and you will hear the next audio track. If you notice that one of the tracks helps you sleep or relax better, you can use the same track choosing it each time you use the Neutonica Band.
- If you want to turn the Neutonica Band before it's finished, you must keep the button pressed until you hear the shut-down confirmation sound.
- Remember that another recommended use for the Neutonica Band is during the day between different activities, taking a relaxing moment with a complete cycle of the Neutonica Band. For example: before or after a work meeting, at the end of a work day, after suffering an intense emotional event, etc.


CONNECTION INSTRUCTIONS AND BATTERY CHARGE.

* Uses a factory rechargeable battery. Do not try to replace it.
For more information go to the website mentioned in this guide.

1.- Remove the device from the band.

2.- Connect the USB cord (included) to the device port and connect the other end to an electric outlet.

3.- When the light is red color, it means that the battery is charging, when the light is green color, it means that the battery is completely charged.

IMPORTANT NOTE.

We are pleased to be able to help you with any situation that has encouraged you to buy a Neutonica Band. If it is for lack of sleep, wishing to sleep better or simply as a gift, we also want to remind you some relevant details about its functioning.

The Newtonic Mask is not a medical device, therefore, it does not interfere with any specific medical conditions. It is completely safe in this sense. However, its function is not to create a sleep state spontaneously and definitely. Sleep is a very complex mechanism, whose functions are partially known by current science. Given that the Newtonica Band fosters the ideal conditions to achieve physiological sleep naturally within the average sleep latency time, please do not feel discouraged if you used it one day and did not obtain perfect sleep (especially if you have been suffering from insomnia for weeks, months, or years). The Newtonica Band is individually proven, however its ongoing and routine use always related to the sleeping habits we have widely spread on our social media, enhances the Newtonica Band's effect. If you desire an ongoing positive outcome in the long term, the Newtonic Mask can help you if you use it every day after a routine of going to bed and constantly for at least 21 days, observing the recommendations. **Take on the 21-day challenge and you will notice the difference.**

"This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation."

TO LEARN MORE ABOUT THE NEUTONICA BAND'S USE.
Visit our website: www.bandaneutonica.com/support

TRIM GUIDES

TRIM LINES