



Brighten your future !

# Best Choice, Best Works

## Contents

- Oil Spray for handpiece
- Scan Spray for CAD/CAM
- Occlusion Checker
- Die Spacer

- Wax Spray
- Crucible Coating Spray
- Stone Die Hardener & Thinner
- Acrylic Resin Liquid for Tray
- Temporary Tooth Liquid
- Stone Separator




» Oil spray for handpiece

## ● Denu Handpiece Lubricant

**Cleans well with no clotted residue**

- Cleaning with botanical oil
- Extension of lifespan for bearing products
- Effective lubrication



**Characteristics**

- The product provides excellent cleaning high & low speed handpieces.
- Its oil does not clog and is less sticky and does not leave any residue. Its alcoholic ingredient provides superior disinfection and high volatility making the oil smell disappear in 2-3 seconds after spraying.
- The product can be used for all kinds of handpieces including Kavo, W & H, NSK and Bien Air.
- Its high compression air in spray container enables highly accurate spraying, which creates minute and uniform coating film on the surface of bearings inside handpiece.

**Instructions**

1. Shake well before using it .
2. After decomposing air turbine handpiece from tubing coupler, tighten spray nozzle to handpiece.
3. Lubricate for 2 ~ 3 seconds till oil come out.
4. Even if you do not use, it is recommended to shertce it once a day for about one or tow secong to keep good



Contents : 520ml×1

## ● Denu Easy Scan

Spray for preparation of surfaces for friction testing of telescopic crowns and CAD/CAM scanning.



220ml(250g)

**Characteristics**

- Scanning the model, it helps accurate and clear results delivered
- Applying thinly and coating evenly, it provides excellent scan results
- Containing 2-3 times more volume(net 250g) than normal products
- Suitable for indirect application(model)

**Instructions**

- Shake it well for about 1 minute before use.
- Thin application over several times is more effective than one-time thick application.

**Cintents**

- Spray 250g(220ml)×1



## ● Denu Occlusion Checker

**Aerosol for articulation marking to check the occlusion of inside of prosthesis in dental laboratory**

- Spray with thin layer
- Easy to wash with water
- Paste to the metal checker

### Instructions

- Shake the product well about 10-20 times before using.
- Press spray button 1-3 times to avoid the checker left in the nozzle.
- Keep enough distance from the model when spraying.

### Storage

- Use outdoors or in a well-ventilated area.
- Keep children away from it.
- Do not expose to direct sunlight or heat.

### Contents

- Spray 150g x 1



150g

## ● Denu Die Spacer

**Spray for providing a space between crown and underlying tooth preparation**

### Characteristics

- One application makes a space of 3-5 $\mu$ m thick. Two applications makes a space of about 10 $\mu$ m thick.
- It dries fast within 10 seconds and its color is not smeared.
- It resists being peeled off when steamed.
- It is used to make inlays, onlays, crowns bridge prosthetic appliances.

### Instructions

- Shake enough before use.
- Adjust concentration as much as desired using thinner.
- Apply 1 or 2 times using a brush on the bottle cap.
- Apply 0.5mm short from margin.
- Dry for 30 to 60 seconds after application and apply a separating agent.

### Packing

- Bottle : 50ml(1.6 oz) / Color : Red



50ml

## ● Denu Wax Agent

**Smooth surface / No bubbles**

### Characteristics

- Provides a smooth surface between investment and wax in the ring furnace by spraying Denu Wax Agent on wax to create a coating film on it, which enables metal alloy or gold to be shaped excellently.
- Allow no bubbles between investment and wax in the ring furnace.

### Instructions

- Shake the spray well before use.
- Lightly spray 2 or 3 times on the wax.
- Shake the spray once a day for 1 or 2 seconds to keep good quality even when not in use.

### Contents

- Spray 420ml x 1



420ml

## ● Denu Crucible Coat Agent

**Crucible Coating**

- Protecting valuable gold from disappearing into the crucible during casting.

### Characteristics

- Applicable for all metal alloys
- No self-pollution against metal
- No need to use sand powder or asbestos
- Prevention of a molten metal from crucible deposition
- Prevents overheating of crucible by inducing even heat conduction inside molten metal.

### Contents

- Spray 420ml x 1



420ml

## ● Denu Stone Die Hardener & Thinner

Keep the liquid thin. Add thinner as much as needed.

### Denu Stone Die Hardener

- The product significantly protects the stone or plaster model, particularly its margins.
- Once coated, there are no scratches by nails or carving knives.

#### Instructions

- Denu Stone Die Hardener is applied to the surface on a stone or plaster model. After applied with 1 to 2 coats, the surface is set and hardened in 3 minutes.
- It protects margins and hardens, waterproofs and seals stone or plaster. But, it does not affect dimensions of die. It is easy to use with a brush on its bottle cap.

#### Contents

- Bottle : Hardener(50ml) ×1



### Denu Thinner For Stone Die Hardener

Denu Thinner is used to maintain the proper thin consistency over the life of resin.

#### Instructions

Add a few drops of the thinner.  
Add the thinner as needed to keep resin thin.

#### Contents

- Bottle : Thinner(50ml) ×1



## ● Denu Acrylic Resin Liquid

Liquid for fabrication of self-curing acrylic resin for tray

- Used to fabricate self-curing acrylic resin for denture repair, temporary tooth and rebasing resin.
- Produces less heat and smell and enables fast setting.
- Usable for fabrication of temporary tooth for crowns or bridges in dental laboratories.

#### Instructions

- Check the expiration date, impurities and chemical change before use.
- Mix powder and liquid at the ratio of 20g(powder) to 5ml(liquid). Mix them until the dough does not stick to the mixing container.
- Make tooth model with the dough and remove extra resin.
- Make a handle by using extra resin.
- Leave the tooth model for about 10 minutes for self-curing.
- Trim the resin by using burs or other instruments.

#### Storage

- Keep the product in a well-ventilated place after closing the container cap.
- Do not keep the product over 25°C(77°F).

#### Contents

- Can : 1,000ml×1



## ● Denu Stone Sep

Separating agent for dental stone or plaster

#### Characteristics

- Denu Stone Sep enables clear separation between stone or plaster models without generating impression during pin working.
- It seals stone or plaster pores completely and allows no build-up layer.
- The surface becomes hard, smooth and waterproof.
- Use the product as much as needed using the empty moist sprayer provided

#### Instructions

- Spray to moist stone. Spray the product on the surface of stone.
- Leave it for around 10 seconds after spraying.
- Wash redundant liquid with running water or blowing air.
- Surface will become sealed and separated.

#### Contents

- Set : Can 500ml×1



500ml (Can)



**denu**  
DENTAL FOR YOU

Good Choice, Best Works



#201, IL-SHIN B/D, 24 Tongil-ro, Jung-gu, Seoul, Korea  
Tel: +82-2-777-3510, E-mail: [hdikorea@hanmail.net](mailto:hdikorea@hanmail.net)  
Web site: [www.hdent.com](http://www.hdent.com)