

Indirect Pulp Capping Placement

Dentistry courtesy of Dr. Ross Nash



1
Indirect pulp capping. Infected dentin is removed and affected dentin remains allowing for incremental layers of TheraCal LC to be placed.



2
TheraCal LC directly applied to preparation in 1mm incremental layers and manipulated into a smooth surface.



3
Each increment is light-cured for 20 seconds.



4
Surrounding enamel is etched using BISCO's SELECT HV® ETCH.



5
BISCO's ALL-BOND SE® is applied and light-cured.



6
Final prep after core build-up using BISCO's CORE-FLO™ DC.

Ordering Information

4 SYRINGE PACKAGE

4 Syringes TheraCal LC (1g ea.), Accessories, Instructions

H-33014P

REFILLS

1 Syringe TheraCal LC (1g), Accessories, Instructions

H-3301P

50 Black Disposable Syringe Tips (22 Gauge)

X-80621N

MC-2334TC



The
**SECURITY
BLANKET**
for
**Pulpal
Protection**

TheraCal LC®




Seal and Protect with TheraCal LC
Pulp Capping Material and Liner



Rx Only

HIGH PERFORMANCE CHEMISTRY

Unique Hydrophilic Matrix



CALCIUM RELEASE

- Hydroxyapatite formation^{1, 2}

HEALING & INSULATION

- Liner insulates the pulp³, providing virtually no post-operative sensitivity
- Promotes healing¹

EASE OF USE

- Light curable, flowable-like handling

RADIOPAQUE

- Radiopaque properties allow for easy detection on x-rays

TheraCal LC Stimulates Apatite Formation

TheraCal LC is a light-cured, resin-modified calcium silicate filled liner designed for use in direct and indirect pulp capping, as a protective base/liner under composites, amalgams, cements, and other base materials. It can be used as a replacement for calcium hydroxide, glass ionomer, RMGI, IRM/ZOE and other restorative materials. TheraCal LC performs as an insulator/barrier and protectant of the dental pulpal complex.

The proprietary formulation of TheraCal LC consists of tricalcium silicate particles in a hydrophilic monomer that provides significant calcium release* making it a uniquely stable and durable material as a liner or base.

- Calcium release^{1*} stimulates hydroxyapatite and secondary dentin bridge formation^{2, 3}
- The alkaline pH promotes healing and apatite formation^{2, 4}
- Significant calcium release¹ leads to protective seal^{5, 7, 8}
- Forms a protective barrier that insulates the pulp^{5, 6}
- Moisture tolerant¹ and radiopaque - can be placed under restorative materials and cements

Direct Pulp Capping Placement

Dentistry courtesy of Dr. Mark Cannon




1 Hemostasis achieved prior to TheraCal LC direct pulp capping placement.



2 TheraCal LC applied directly to exposed pulp and light-cured in 1mm increments.



3 Etched, rinsed, and bonding applied.



4 Continued restoration of the tooth.



5 Final restoration.

PHYSICAL PROPERTIES

	Shear Bond Strength* (MPa)	Radiopacity (mm Al)	Calcium Release (24 h)
TheraCal LC	4.4(2.9)	2.6	213 (µg/cm ²)
Prisma® VLC Dycal®	0.9(0.9)	0.8	NA

* Visibly moist bonding.

Prisma and Dycal are registered trademarks of Dentsply Caulk.

¹ ADA definitions for direct and indirect pulp capping at <http://www.ada.org/en/publications/cdt/glossary-of-dental-clinical-and-administrative-ter>

² Apatite-forming Ability of TheraCal Pulp-Capping Material, M.G. GANDOLFI, F. SIBONI, P. TADDEI, E. MODENA, and C. PRATI J Dent Res 90 (Spec Iss A):abstract number 2520, 2011 (www.dentalresearch.org)

³ Selcuk SAVAS, Murat S. BOTSALI, Ebru KUCUKYILMAZ, Tugrul SARI. Evaluation of temperature changes in the pulp chamber during polymerization of light-cured pulp-capping materials by using a VALO LED light curing unit at different curing distances. Dent Mater J. 2014;33(6):764-9.

Advantages:

- Improved seal and bond to deep moist dentin
- Strong physical properties
- High radiopacity
- Significant calcium release

A Protective Liner: Liquid Apatite at Your Fingertips

- A protective liner for use under restorative materials, cements or other base materials.
- Pulp capping agent: TheraCal LC may be placed directly on pulpal exposures after hemostasis is obtained. It is indicated for any pulpal exposures, including carious exposures, mechanical exposures or exposures due to trauma.

* BISCO has, on file, the calcium release data for TheraCal LC.

¹ Gandolfi MG, Siboni F, Prati C. Chemical-physical properties of TheraCal, a novel light-curable MTA-like material for pulp capping. International Endodontic Journal. 2012 Jun;45(6):571-9.

² ADA definitions for direct and indirect pulp capping at <http://www.ada.org/en/publications/cdt/glossary-of-dental-clinical-and-administrative-ter>

³ Apatite-forming Ability of TheraCal Pulp-Capping Material, M.G. GANDOLFI, F. SIBONI, P. TADDEI, E. MODENA, and C. PRATI J Dent Res 90 (Spec Iss A):abstract number 2520, 2011 (www.dentalresearch.org)

⁴ Okabe T, Sakamoto M, Takeuchi H, Matsushima K (2006) Effects of pH on mineralization ability of human dental pulp cells. Journal of Endodontics 32, 198-201.

⁵ Sangwan P; Sangwan A; Duhan J; Rohilla A. Tertiary dentinogenesis with calcium hydroxide: a review of proposed mechanisms. Int Endod J. 2013; 46(1):3-19

⁶ Selcuk SAVAS, Murat S. BOTSALI, Ebru KUCUKYILMAZ, Tugrul SARI. Evaluation of temperature changes in the pulp chamber during polymerization of light-cured pulp-capping materials by using a VALO LED light curing unit at different curing distances. Dent Mater J. 2014;33(6):764-9.

⁷ Cantekin K. Bond strength of different restorative materials to light-curable mineral trioxide aggregate. J Clin Pediatr Dent. 2015 Winter;39(2):143-8.

⁸ Mechanical Properties of New Dental Pulp-Capping Materials Over Time. M. NIELSEN, R. VANDERWEELE, J. CASEY, and K. VANDEWALLE, USAF, JBSA-Lackland, TX, J Dent Res 93(Spec Iss A): 495, 2014 (www.dentalresearch.org)