

REUSE GUIDE

How to safely reuse Wesper patches.



STEP 1: CONFIRM ALL COMPONENTS

Every Wesper kit should come with two biometric patches, one charging pad, two charging cables, a quickstart guide, a pulse oximeter, and unused adhesives.



STEP 2: SANITIZE ALL COMPONENTS

- ✓ Use an alcohol wipe to clean the top and bottom of the patch, and the sensor window until all the surfaces are clean from any stains or debris.
- ✓ Wipe down the whole device with an alcohol wipe until it's clean from any stains or debris
- ✓ Use an alcohol wipe to clean the top and bottom of the charger pad until it's clean from any stains or debris.
- ✓ Wipe down the charger cable twice with an alcohol wipe until it's clean from any stains or debris
- ✓ Wipe down with an alcohol wipe until the linear is clean from any stains or debris.

STEP 3: PERFORM VISUAL INSPECTION


Once cleaned, perform a visual inspection. All components should be free of stain, deformation, scrapes and damages..

For the patches, pay extra attention to the fabric top, silicone bottom, sealing edges, and sensor window on the bottom.

If you have any stains, deformations, scrapes, or damages do not give the kit to a patient. Instead, order a new kit as soon as possible.

STEP 4: CHARGE THE KIT FOR NEXT PATIENT

It's important to set up the kit for your next patient so they can have the best experience possible.



Charge the patches for up at least three hours.



Make sure the patch is dry and clean.



Apply a new adhesive to both patches.



Remember to charge your CheckMe pulse oximeter until the display screen shows a fully charged icon. Remove the charging cable and connect the finger probe back on the device.

STEP 5: REPACK ALL COMPONENTS

Repack each component in the case for patients. Remember to select the number of adhesives (two per night of testing).



STEP 6: GIVE TO A PATIENT

Once all the above steps have been completed, it is now time to give it to another patient.

THANK YOU FOR CHOOSING WESPER FOR YOUR CLINIC.
WE LOOK FORWARD TO CONTINUING TO WORK
TOGETHER TO HELP THE WORLD SLEEP BETTER.

QUESTIONS?

Our team is standing by to answer your questions.
Reach out at providersupport@wesper.co