

## ECG:

## CABLE NOT FOUND

Office Medic Software is unable to detect the ECG whilst plugged in or there is a security policy in place blocking the QRS ECG from being detected.

## Solution -

- 2. If using a USB extender, disconnect and plug QRS ECG directly into PC.
- 3. Unplug the ECG from the computer.
  - a. Close Office Medic / Best Practice / Medical Director.
  - b. Plug the ECG into the computer.
  - c. Open Office Medic / Best Practice / Medical Director.
  - d. Test again.
- 4. Search 'Core Isolation' in Windows Search.
  - a. Turn off Memory Integrity.
- 5. Check Anti Virus or Windows Defender to see if anything is blocking the USB Device.
- 6. Uninstall Office Medic and all QRS Vectracor Drivers.
  - a. Control Panel > Programs > Uninstall Programs.
  - b. Uninstall Vectraplex\Office Medic.
  - c. Uninstall all 'Windows Driver Package Vectracor' and 'Windows Driver Package QRS Diagnostic' drivers.
  - d. Uninstall any ATEN and FTDI drivers as well. These are old pre 2012 QRS Drivers that may cause conflict.
- 7. Goto Windows Device Manager and remove QRS Devices from USB Controllers.
  - a. Goto Windows Start Menu, type 'Device Manager' and select it.
  - b. Goto View > Show Hidden Devices.
  - c. Goto Universal Serial Bus Controllers, select all 'Universal ECG Bus Controller' and 'QRS Spirometry' devices and uninstall.
  - d. If 'Unknown USB Device (Device Descriptor Request Failed) is there, then uninstall that as well.
- 8. Restart computer.
- 9. Reinstall Office Medic.
- 10. Plug the QRS ECG into the computer. Ensure it is plugged in before opening Office Medic.
- 11. Open Office Medic and test the ECG.
- 12. If still getting errors, email <a href="mailto:support@medilogic.com.au">support@medilogic.com.au</a> for further assistance.

## Note:

The blue QRS Unismart ECG is only compatible with Office Medic Version 6 and up.

The black QRS Universal ECG is only compatible with Office Medic Version 5 to 5.8.2.