



# Versana Premier<sup>TM</sup> Ultrasound

**Powerful imaging.  
Easy cleaning.**

**EASY CONSOLE AND PROBE CLEANING SUPPORTS  
COMPLIANCE WITH YOUR INFECTION-CONTROL POLICIES.**

Regular cleaning and disinfection of your Versana Premier ultrasound system can help prevent hospital-acquired infections (HAIs), a persistent problem that affects thousands of patients each year. Easy-to-clean equipment helps practitioners meet infection-control goals while maintaining workflow efficiency.

Versana Premier ultrasound systems and probes can withstand challenging healthcare environments and are compatible with multiple cleaning and disinfection agents.

## Save time and effort when cleaning

Versana Premier design features simplify cleaning:

- Removable trackball and easy-lock controls helps reduce cleaning time
- Washable keyboard film cleans easily and excludes dust from keyboard spaces
- Removable air filter simplifies cleaning



## Enjoy diverse options for probe reprocessing

You can choose from multiple branded agents to clean and disinfect probes for your Versana Premier system. They include:

- 27 agents for cleaning
- 88 agents for lower-level disinfection
- 27 agents for high-level disinfection

For a complete list of compatible reprocessing agents, consult your user manual or visit <https://www.gehealthcare.com/en/products/ultrasound/ultrasound-transducers>.

In addition, many Versana Premier probes\* are compatible with the GE trophon EPR® high-level disinfection system. The following probes can be used with trophon EPR and trophon2:

- |              |             |             |
|--------------|-------------|-------------|
| • 3Sc-RS     | • 6S-RS     | • 12S-RS    |
| • 12L-RS     | • L8-18i-RS | • 4C-RS     |
| • E8C-RS     | • E8Cs-RS   | • RAB2-6-RS |
| • RIC5-9A-RS | • 8C-RS     | • L6-12-RS  |

Sterile/sanitary sheaths are to be used on the probe during its actual use with patients. Wearing gloves protects the patient and operator.

## Select agents with ease

Testing confirms that the Versana Premier console is compatible with a variety of agents used for routine reprocessing. Choose from 10 agents that will either clean or disinfect your console:\*\*

### Disinfectants

- |                              |                         |
|------------------------------|-------------------------|
| • Tristel Wipes              | • CaviWipes™            |
| • Protex                     | • Sono Ultrasound Wipes |
| • PDI Super Sani-Cloth® Plus | • PDI Sani-Cloth® HB    |
| • PDI Sani-Cloth® Plus       | • Clorox® wipes         |

### Cleaners

- PDI Easy Screen Cleaning Wipe
- ALCOHOL (ISOPROPANOL) 70%

\* Product and probes may not be available in all countries and regions. Full product technical specifications is available upon request. Contact a GE Healthcare Representative for more information.

\*\* Cleaners and disinfectants listed for reprocessing of console include approval to use on the following areas: system cabinet, operator control panel, main display, touch screen and probe holder.

The pictured product may not reflect the actual physical appearance (e.g., color and design) of the commercially available device.

## Imagination at work

© 2019 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Versana Premier are trademarks of General Electric Company. CaviWipes is a trademark of Metrex. Sani-Cloth is a trademark of PDI. Clorox Healthcare is a registered trademark of The Clorox Company. trophon is a trademark of Nanosonics, Ltd. GE Healthcare is a distributor of trophon. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

August 2019  
JB69553XX

