


IMPORTANT TO READ BEFORE USING YOUR NEW MAGNIFYING GLASSES

Magnifying glasses are devices designed to enhance the vision of individuals with low vision. However, it is crucial to take certain precautions before using them.



It is important to note that they are not glasses for tired eyes or presbyopia; these glasses are visual aids equipped with high diopter lenses, allowing you to see letters or details in an enlarged size.

One disadvantage of these glasses is the reading distance; the higher the power of the diopters, the closer you will need to be to the text. For example:

+8.00 requires an optimal reading distance of 12 cm.

+10.00 requires an optimal reading distance of 10 cm.

In summary, when using them, you will need to get very close to the paper, similar to what watchmakers do when they work.

Our recommendation is to start using the glasses for short periods during the first few days and then gradually increase the wearing time as you feel more comfortable.

Thank you for your purchase! We appreciate your choosing us.