




# EasyAngle

## Product description

**easyangle** is a new type of goniometer that enables physiotherapists to assess range of motion of all joints with accuracy and ease. The tool addresses the limitations of existing solutions introducing a new standard for joint assessment.

### easyangle

- Designed for one hand operation - enabling support of patient
- Easy to read display
- Saves the last 5 measurements - no more scribbling angles on paper
- Can be used for all joint measurements, including rotations



The **easyangle** workflow illustrated in the image below.



*Align the device with the first limb and click the button.*



*Align the device with the second limb and click the button*



*The joint angle is calculated and displayed.*

**easyangle** is the only tool a physiotherapist need for high quality joint assessment.

## Technical specifications

Product name	EasyAngle
Operating environment	Not for home use
Operating temperature	0°C to 45°C
Enclosure	IP4X
Operating time*	2 weeks
Stand-by time	11 weeks
Charging time	2 hours
Sensor accuracy	± 1° within 180°
Certification	Certified as a Class 1 medical device according to MDD 93/42/EEC (LVFS 2003:11).

\*Assuming 12 minutes daily usage

## EasyAngle kit includes

- 1x EasyAngle device
- 2x Alignment guides, 280 and 195 mm
- 1x USB wall charger
- 1x USB cable for charging
- 1x Multilingual manual (SE/EN/FI/DK/DE/NL)


## Product video



<https://youtu.be/dsfLLOLvHsQ>




## Mechanical dimensions



All dimensions are in mm.  
General tolerance 1 mm


## Alignment guide

The alignment guide can be detached or changed as shown in the image below.



## Measurement guidelines

The EasyAngle can measure rotations in all planes. However, the device must be moved as shown in the image below to achieve optimal measurement results.



## Spherical alignment

If the alignment guide is detached the sensor house is designed for rigid 3 point contact on a spherical surface, such as a head. This can be used for continuous movement measurements.

An example can be seen in the link below where a measurement of cervical rotation is performed.

