# CERTIFICAZIONI CERTIFICATION

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*The developed products, studies, tests, manufacturing and assistance to customer are regulated and monitored by a certified Quality System according to ISO 13485-2016 for medical components by TÜV Rheinland and are entirely made in GIESSE TECHNOLOGY without the presence of external suppliers except for the sterilization process.*

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## Certification


*The continuous attention GiEsse Technology has towards the quality of its products is testified by constant laboratory tests and analyses, conducted to verify their safety towards patients and the ability to achieve the declared performances. To carry out these tests, GiEsse Technology avails itself of the collaboration of prestigious University Research Centers, such as the scientific-technological university Politecnico Milano and the University of Trieste, as well as laboratories adopting a certified Quality Management System for the execution of technical consultancy services in the biomedical field and laboratory services for microbiological and physical analyses. Tests and analyses are conducted according to the European standards (EN standards) and international standards (ISO standards), and according to the FDA (Food and Drug Administration) guidelines, in order to demonstrate the compliance of the products with the highest safety and performance standards.*

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
*... has also successfully started up and completed, the CE marking in accordance of its Medical Devices in accordance with Directives and EU Regulations.*

# TEST A FATICA FATIGUE TEST

Test con moncone dritto - Test with straight abutment

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GiEsse Technology's dental implants have been subjected to mechanical tests to evaluate their resistance to operating stresses.

The tests were performed on the complete implant system, including the endosseous implant and the abutment, according to the load scheme defined in UNI EN ISO 14801: 2016 "Dentistry - Implants - Dynamic loading test for endosseous dental implants". During the tests, which were conducted at the Materials' Testing Laboratory of Politecnico Milano and at the laboratory of CERMET, a certification body, the implant system was subjected to the application of 5,000,000 loading cycles. All tested GiEsse Technology dental implants have shown to have a practically infinite operating life (conventionally set at 5,000,000 cycles) under normal mechanical stresses produced by chewing.

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## MORPHOLOGICAL ANALYSIS OF SURFACE

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### Absence of Superficial Contamination

*The chemical analysis of the surface of the HSS treated implant highlights the absence of contaminating elements, and in particular, of the alumina used for sandblasting (present on the surfaces treated only with sandblasting), and the presence of titanium oxide phosphate hydrate, introduced through HSS treatment. This type of oxide promotes cell adhesion by osteoblasts, with a consequent increase in the quality of osseointegration.*

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
*HSS treatment also has a positive influence on the wettability of the implant surface.*

*In the following images, it is possible to notice the greater wettability of the surface of an HSS-treated implant compared to the only sandblasted surface.*

*The contact angle of the cell culture serum passes from 55 °C for the sandblasted surface to only 38 °C for the HSS treated surface. This parameter is an index of the speed of protein adsorption which will then influence the cell adhesion and consequently the osseointegration.*

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**TRATTAMENTO HSS**  
**HSS TREATMENT**  
High Speed Surface

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
*At the University of Trieste, the surface of GiEsse Technology dental implants treated with HSS was subjected to cell proliferation analysis to evaluate its biological response.*

*Cell proliferation on the surface of the HSS implant was extremely high in the first nine days of cell application compared to an identical implant treated only with sandblasting, resulting in a faster cell differentiation. After the ninth day, growth trends are similar. It is important to pay attention to the fact that the cell proliferation of the HSS system is similar to that of the cells present on the control. Below, we see images obtained through the epifluorescence microscope. Cells show a typical morphology of fibroblasts derived from human osteosarcoma (MG-63); we observe a massive concentration of cells highly adherent to the substrate (presence of many filopodia), organized on several layers and a massive amount of extracellular matrix.*

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*Epifluorescence microscope images. Cellular cytoskeleton, labelled with FITC - phalloidin (green) and the nucleus labelled with DAPI (blue).*

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## FITTING EVALUATION FIXTURE-ABUTMENT

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

### **Abutment-Implant Coupling**


*The quality of the implant-abutment coupling was assessed by micro computerized tomography, studying the case of the conical-hexagonal connection.*

*The results showed that the abutment-implant connection is highly coupled with few and very small spaces, among abutments, implant collar and connection screw. It is believed that these spaces, smaller than 25 microns and present only within the connection, are not sufficient to allow the penetration of bacterial agents.*

*The three-dimensional reconstructions shown in the images below underline the presence of several points of contact both between the healing screw and the implant and between the abutment and the internal part of the implant collar.*

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*The decontamination process, validated and implemented over the years, has been renewed by GiEsse Technology thanks to the collaboration with the University of Trieste.*


*In fact, a new CLEANING WITH PLASMA REACTOR process has been introduced. This, carried out at the end of the classic ultrasonic washes, allows an increase in the decontamination of the implant surface, reducing the risk of peri-implantitis and favoring the osseointegration process.*

*The implants to be treated are placed in a reactor where they are exposed to the flow of an high-power plasma inert gas (Argon). Thanks to the kinetic impact energy of the ions, atoms and molecules of contaminants and organic components present on the surface of the implant are expelled and then removed thanks to the vacuum pump.*

*Compared to traditional cleaning treatments, plasma decontamination does not require the use of cleaning liquids that could leave residues or impurities that could alter the osseointegration process; moreover, the flow of plasma is able to penetrate into slots or areas (e.g. etching slots) which are difficult to reach by traditional cleaning liquids.*

*The tests of cellular contamination and proliferation conducted in collaboration with the University of Trieste have highlighted the excellent level of cleanliness of the surface and the reduced presence of pro-inflammatory contaminants obtained thanks to the adoption of this new decontamination technology. Furthermore, the surfaces treated with plasma have a high wettability and are more prone to promote cell adhesion.*

*The reproducibility of the treatment and the considerable possibility to control the process variables make it possible to exploit this type of treatment even with products of complex geometry, offering a very high consistency in terms of quality.*

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*Cell activation markers production by U937 human macrophages.*

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**DOPPIA FIALA IN PETG BIOCOMPATIBILE**  
**DOUBLE BOX PETG BIOCOMPATIBILE**

*Dental implants are manufactured under conditions of absolute cleanliness and hygiene. Analyses conducted by external laboratories have confirmed that the level of pre-sterilization microbial contamination is very low, if not negligible, both before and after washing. The washing cycles adopted are effective both in the actual washing phase, completely eliminating the residues of processing oils and greases, and in the rinsing phase, removing all traces of the detergents used.*

*The packaging of GiEsse dental implants is carried out with an implant supported by a medical-grade titanium millable mounter and a retaining screw. The plastic materials used for packaging, which consists of an internal and external vial, with the presence of two caps and a security seal, are phthalate-free, comply with current regulations, and are suitable for gamma-ray sterilization.*

*The packaging of dental implants is performed within a controlled biological and particle contamination environment (clean room). The environmental contamination class has been determined so as to guarantee a constant and sufficiently low pre-sterilization bioburden on the product to ensure their sterilization.*

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*The gamma-ray sterilization treatment has been validated according to the UNI EN 11137 standard, and is periodically subjected to control tests to verify that it always remains effective and that the sterilization parameters do not present variations from what is established in the initial validation protocols.*

*Once the implants have been sterilized, their standard packaging is suitable to preserve the integrity of the sterile barrier system (i.e. to maintain the sterility of the contained implant) for a period of 5 years.*

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