

Upgrading to a New Version of QRlogix Software

For V1.5 and earlier.


1. Export any current templates to the Templates folder or another safe location.
 - a. C:\Users\Public\Documents\QRlogix\Templates
2. Take screenshots of your template forms in case they need to be rebuilt.

The screenshot displays the QRlogix Software interface. On the left is a sidebar with navigation options: Dashboard, Template 1, Cassette DS1, My template (selected), EB DS1, and File Drop 1. Below the sidebar are links for Import Template, Administrator, and Setup and Preferences. The main area is titled 'My template' and contains a form with the following fields: Prefix (SL23), Accession / Quantity (400), Specimen / Quantity (A), Block / Quantity (1), Name (PAN, PETER), and Stain (H&E). A blue 'PRINT' button is at the bottom of the form. To the right of the form is a preview of the label, which reads: 'SL23-400-A1', 'Lev 1 H&E', 'PAN, PETER', and '03/14/2024'. Below the preview is a 'Print History' section with a 'Show All' button and a list of ten print records, each stating 'Printed by: Administrator' on 'Mar/14/2024 at 15:28 PM'.


3. Copy/save the label designer code from any active templates in case they need to be rebuilt.

```
^XA
^CIO
^FT5,50
^A0N,33,31^FD@prefix@-@case@-@sample@@@block@-@slide2@^FS
^FT5,90
^A0N,33,31^FD@name@^FS
^FT5,130
^A0N,33,31^FD@panelname@^FS
^FT5,170
^A0N,33,31^FD@panelstain@^FS
^FT5,210
^A0N,33,31^FD@pathologist@^FS
^PQ1,1,1,Y
^XZ
```

4. Backup contents of your Sqlite folder located at: C:\ProgramData\QRlogixLMS\Sqlite (Or it could be located at C:\ProgramData\QRlogixSoftware\Sqlite)
 - a. Create new folder = Old
 - b. Place all contents of Sqlite folder into "Old" folder




5. Uninstall software in Control Panel



6. If PrintShare (File Drop) is installed, uninstall the AutomationService (File Drop Service) application in Control Panel.

Confirm the following folder was removed during uninstallation:
C:\AutomationService or C:\FileDropService
7. Download the new version of the software from the link provided.
8. Download the new version of File Drop if you are printing from CSV/TXT or connecting to LIS.
9. Follow the software installation instructions found in the User Manual.
10. If connecting to File Drop (formerly PrintShare), follow the installation instructions found in the User Manual.

11. Use the Import Template button to load previous templates.



12. If you have any problems loading old templates, email the exported template files, along with the form screenshot and label designer code to support@QRlogix.com and we'll rebuild the templates in the current version and send you the files.

Support:

Phone: 800-650-0632

Email: Support@QRlogix.com

<https://qrlogix.com/pages/qrlogix-software-support>