## **Technical Data Package**

C6155.006: cGMP SUITABLE HEAVY-DUTY
WASHDOWN LARGE CAPACITY ORBITAL
SHAKER WITH CUSTOM TAILGATE TOP
2" STROKE, 1-300 RPM, 500LB LOAD CAPACITY
230V-240V, 50/60Hz
SINGLE PHASE, PLUG TYPE F



5900 Schooner St, Belleville MI 48111 – USA TEL (800) 422-2558 – FAX (734) 665-9099 www.Eberbachlabtools.com

## **TDP Table of Content**

| I.   | Machi | ine Specifications                               | 3  |
|------|-------|--|----|
|      | a.    | Mechanical Specifications                        | 4  |
|      | b.    | Electrical specifications                        | 4  |
|      | c.    | Maintenance                                      | 4  |
|      | d.    | Load Capacity Chart                              | 5  |
|      | e.    | Controls   | 6  |
|      | f.    | Alarms & Faults                                  | 9  |
| II.  | Machi | ine Testing Qualification & validation           | 10 |
|      | a.    | Factory Acceptance Test                          | 12 |
|      | b.    | Installation Qualification                       | 13 |
|      | c.    | Operation Qualification                          | 14 |
|      | d.    | Site Acceptance Test                             | 15 |
|      | e.    | Performance Qualification                        | 16 |
| III. | Machi | ine Material of Construction (MOC) Documentation | 17 |
|      | a.    | Material Certification                           | 18 |
|      | b.    | Electronics Data                                 | 19 |
|      | C.    | Assembly Prints                                  | 23 |
|      | d.    | Wiring Diagrams                                  | 27 |
|      | e.    | Material Information                             | 31 |
|      | f.    | Spare Parts List                                 | 35 |
| V.   | Custo | omer Signoff                                     | 36 |

## I. Machine Specifications



5900 Schooner St, Belleville MI 48111 – USA TEL (800) 422-2558 – FAX (734) 665-9099 www. Eberbachlabtools.com

#### I.a <u>Mechanical Specifications</u>

- Stroke: 2.00" +/-0.04"
- Speed: Set to 1-300 RPM, set speed will be +/-2 RPM
   \*speeds over 150 RPM are not recommended without specifically engineered tops
- Length: 43.12" +/-0.50"Width: 43.12" +/-0.50"
- Height Unloaded: 33.50" +/-1.00"Weight Unloaded: 2350lb +/-10lb
- Max Capacity: 500 lbs (See load chart for speed specific limits)
- Tailgate Loading Capacity: 150lbs
- Base: Features included for use of narrow (21") forklift or pallet jack
- CFD engineering services available upon request

#### I.b <u>Electrical Specifications</u>


- 230-240V Power Cord (Plug Type F)
- 750W, IP65 Servo Motor with 10:1 Gear Box
- 13.2A Max amperage draw. This calculation is based on a motor overload (maximum amps, not normal running conditions) and max outlet amperage draw. Eberbach recommends a 15A outlet
- 1.6kWh max. This calculation is based on a motor overload (maximum wattage, not normal running conditions) and max outlet amperage draw
- Emergency Stop Button
- RJ45 Ethernet Port for accessing ProfiNet system
- USB Port for accessing Data Logs on HMI
- Washdown 4" HMI touch panel for controlling machine
- IP66 Enclosures for the drive & low voltage electronics
- Remote diagnostic functionality available upon request

#### I.c Maintenance

- Wipe exposed surfaces of machine with isopropyl alcohol twice per month or as needed.
- Keep touch screen free from dirt/debris.
- Verify E-Stop functionality once every three months unless local safety ordinances require otherwise.
- Ensure Cord Jacketing is not frayed or Broken.
- Replace the USB port caps when the ports are not in use.

## I.d Loading Capacity Chart

Always place containers in such a way as to balance the load in the center of the carrier.




#### I.e Controls



WELCOME SCREEN. Welcome screen is accessed by pressing F4. You will be able to end Simatic RUN TIME by pressing the Power button in the bottom left of the screen. 'Run' button goes to RUN SCREEN. 'Timed Run Button' goes to TIMED RUN SCREEN. 'Alarm' button goes to PLC ALARM BUFFER SCREEN. 'Info' button will take users to Machine specs.

#### **RUN SCREEN**



- 1) Set Speed Input Box
- 2) Set Speed Adjustment Arrows (1 RPM increments)
- 3) Set Speed Indicator
- 4) Elapsed Run Time (999:59:59) max value
- 5) Start Button (starts rotation)
- 6) Stop Button (stops rotation)
- 7) Elapsed Run Time Reset Button. Resets Time to (0:00:00)
- 8) Tachometer Box (displays actual speed of machine)
- 9) Load Meter (monitors torque load as percentage)

### **Programing Software**

## Installed software

Totally Integrated Automation Portal

Version V14 SP1 Update 9


▼ STEP 7 Professional

Version V14 SP1 Update 9

▼ WinCC Advanced

Version V14 SP1 Update 9

#### I.f Alarms & Faults




PLC ALARM BUFFER SCREEN. Shows PLC Alarms. Press the Home Icon to go to 'WELCOME SCREEN.' Press the 'Drive Faults' button to go to 'DRIVE ALARMS' screen.

## II. Machine Testing Qualification & Validation



5900 Schooner St, Belleville MI 48111 – USA TEL (800) 422-2558 – FAX (734) 665-9099 www.Eberbachlabtools.com



#### 1. Objective

The purpose of this protocol is to provide comprehensive document that validates the C6155.006 conforms with the quoted specifications; after manufacturing (FAT), prior to shipping (IQ), upon receipt at client's facility(OQ), and after installation at client's facility (SAT).

#### 2. Scope

This technical data package includes inspection procedures and acceptance criteria of tests to be performed. The scope of this document is limited to the C6155.006.

#### 3. Responsibilities

Responsibilities between Eberbach and Client are defined as follows:

| Task                                | Responsible Group |
|-------------------------------------|-------------------|
| Prepare FAT/IQ/OQ/SAT Documentation | Eberbach          |
| Perform FAT                         | Eberbach          |
| Provide FAT Data to Customer        | Eberbach          |
| Perform IQ/OQ/SAT                   | Client            |
| Review & Verify the IQ/OQ/SAT       | Client            |
| Approve the IQ/OQ/FAT/SAT           | Client            |

#### Eberbach is responsible for:

- Preparing the IQ/OQ/FAT/SAT documentation
- Performing FAT
- Provide FAT documentation to customer
- o Installation support
- Use and Care manuals

#### • Client is responsible for:

- o Reviewing and approving this IQ/OQ/FAT/SAT documentation
- Performing the IQ/OQ/SAT
- All rigging activities required
- All electrical installation requirements

#### 4. Reference

o C6155.006 Use and Care

## Factory Acceptance Test (Eberbach's Responsibility) II.a

C6155.006 230-240V 50/60Hz

| Serial#:  |                                       | Production Order#:   |                  |  |
|-----------|---------------------------------------|--|------------------|--|
| Step      | Description                           | Specification or Acceptance Criteria   | Result Pass Fail |  |
| 1         | Validate Run,<br>unloaded             | Set machine to run at 50 rpm. Speed must be within +/- 2 rpm for 5 min                       | - Fass           |  |
| 2         | Validate Run,<br>unloaded             | Set machine to run at 100 rpm. Speed must be within +/- 2 rpm for 5 min                      |                  |  |
| 3         | Validate Run,<br>unloaded             | Set machine to run at 300 rpm. Speed must be within +/- 2 rpm for 5 min                      |                  |  |
| 4         | Validate Timed<br>Run, unloaded       | Set Machine for a 5 minute timed run @ 50rpm.  Verify the machine stops after 5 minutes.     |                  |  |
| 5         | Verify E-stop<br>(if applicable)      | During machine operation, verify E-stop functionality (at all locations if multiples exist). |                  |  |
| 6         | Confirm machine finish                | Machine panels and base are free from scratches. (cGMP compliant SS if applicable).          |                  |  |
| 7         | Confirm machine finish                | Confirm that there are no sharp corners on machine / accessories.                            |                  |  |
| 8         | Validate stroke                       | Does the machine provide a 2" displacement?  Must be within +/04"                            |                  |  |
| 9         | Verify<br>firmware<br>version         | Verify the firmware version  |                  |  |
| 10        | Cumulative<br>Run Time<br>Validation  | Machine Run Off Time, min 8hrs:  |                  |  |
| 11        | Vibration<br>Analysis                 | Vibration Analysis Completed:  |                  |  |
| 12        | Machine<br>documentation<br>Completed | Verify the Use and Care for the machine is completed and packaged with the machine.          |                  |  |
| Checked b | oy:                                   |  |                  |  |

## **II.b** Installation Qualification

Eberbach will provide on-line execution of the IQ & OQ with customer team for 1 day.

| Chara            | Description  | 6 16 11 A 1 6 11 1   | Result |      |
|------------------|--|--|--------|------|
| Step Description |  | Specification or Acceptance Criteria   | Pass   | Fail |
| 1                | Space<br>Requirements  | Measure installation space for placement of equipment. 45" x 45" (Width x Depth) |        |      |
| 2                | 2 Electrical Check for electrical requirements. 2 20-240V, 50/60 Hz, 1-Phase   |  |        |      |
| 3                | Inspect shipped products for signs of damage that may have occurred during shipping.  Equipment Keep packaging material until Operation Qualifications has been complete.  Contact Eberbach if any shipping damage has occurred. |  |        |      |
| 4                | Equipment documentation provided by Eberbach   | Verify receipt of Use and Care from Eberbach                                     |        |      |
| 5                | Review FAT   | Verify that FAT has been completed by<br>Eberbach                                |        |      |

## **II.c** Operation Qualification

Eberbach will provide on-line execution of the IQ & OQ with customer team for 1 day.

| Ston                                       | Description  | Specification or Assentance Criteria              | Result |      |
|--|--|---|--------|------|
| Step Description                           |  | Specification or Acceptance Criteria              | Pass   | Fail |
| 1  | IQ completed Has the IQ been performed by the client?  |   |        |      |
| 2  | Inspect Validate that the equipment was received in good condition and that there is no damage to upon receipt the equipment or packaging. |   |        |      |
| 3 Equipment set up Unpa                    |  | Unpack and assemble the equipment per IQ.         |        |      |
| 4 Equipment set up Verify the e completed. |  | Verify the equipment/machine set up is completed. |        |      |
| 5 Begin SAT Client to conduct SAT          |  | Client to conduct SAT                             |        |      |

## Site Acceptance Test (Clients Responsibility) II.d

Serial#:\_\_\_\_\_

| Sten Description |                                       | Consideration on Assentance Criteria   | Result |      |
|------------------|---------------------------------------|--|--------|------|
| Step             | Description                           | Specification or Acceptance Criteria   | Pass   | Fail |
| 1                | Validate Run,<br>unloaded             | Set machine to run at 50 rpm. Speed must be within +/- 2 rpm for 5 min                       |        |      |
| 2                | Validate Run,<br>unloaded             | Set machine to run at 100 rpm. Speed must be within +/- 2 rpm for 5 min                      |        |      |
| 3                | Validate Run,<br>unloaded             | Set machine to run at 300 rpm. Speed must be within +/- 2 rpm for 5 min                      |        |      |
| 4                | Validate Timed<br>Run, unloaded       | Set Machine for a 5 minute timed run @ 50rpm.  Verify the machine stops after 5 minutes.     |        |      |
| 5                | Verify E-stop<br>(if applicable)      | During machine operation, verify E-stop functionality (at all locations if multiples exist). |        |      |
| 6                | Confirm machine finish                | Machine panels and base are free from scratches. (cGMP compliant SS if applicable).          |        |      |
| 7                | Confirm machine finish                | Confirm that there are no sharp corners on machine / accessories.                            |        |      |
| 8                | Validate stroke                       | Does the machine provide a 2" displacement?<br>Must be within +/04"                          |        |      |
| 9                | Verify<br>firmware<br>version         | Verify the firmware version  |        |      |
| 10               | Cumulative<br>Run Time<br>Validation  | Machine Run Off Time, min 8hrs:  |        |      |
| 11               | Vibration<br>Analysis                 | Vibration Analysis Completed:  |        |      |
| 12               | Machine<br>documentation<br>Completed | Verify the Use and Care for the machine is completed and packaged with the machine.          |        |      |

## **II.e** Performance Qualification

Eberbach provides custom engineering solutions. Please contact Eberbach for custom IQ/OQ/PQ for specific work.

| Cton        | Description  | Specification or Accontance Criteria   | Result |      |
|-------------|--|--|--------|------|
| Step        | Description  | Specification or Acceptance Criteria   | Pass   | Fail |
| 1           | Validate Run,<br>unloaded  | Set machine to run at 50 rpm. Speed must be within +/- 2 rpm for 5 min                       |        |      |
| 2           | Validate Run,<br>unloaded  | Set machine to run at 100 rpm. Speed must be within +/- 2 rpm for 5 min                      |        |      |
| 3           | Validate Run,<br>unloaded  | Set machine to run at 300 rpm. Speed must be within +/- 2 rpm for 5 min                      |        |      |
| 4           | Validate Timed<br>Run, unloaded  | Set Machine for a 5 minute timed run @ 50rpm.  Verify the machine stops after 5 minutes.     |        |      |
| 5           | Verify E-stop<br>(if applicable)   | During machine operation, verify E-stop functionality (at all locations if multiples exist). |        |      |
| 6           | Confirm machine finish   | Machine panels and base are free from scratches. (cGMP compliant SS if applicable).          |        |      |
| 7           | Confirm machine finish   | Confirm that there are no sharp corners on machine / accessories.                            |        |      |
| 8           | Validate stroke  | Does the machine provide a 2" displacement?  Must be within +/04"                            |        |      |
| 9           | Verify<br>firmware<br>version  | Verify the firmware version  |        |      |
| 10          | Cumulative<br>Run Time<br>Validation   | Machine Run Off Time, min 8hrs:  |        |      |
| 11          | Vibration<br>Analysis  | Vibration Analysis Completed:  |        |      |
| 12          | Machine documentation Completed  Verify the Use and Care for the machine completed and packaged with the machi |  |        |      |
| Checked by: |  |  |        |      |

# III. Machine Material of Construction (MOC) Documentation



5900 Schooner St, Belleville MI 48111 – USA TEL (800) 422-2558 – FAX (734) 665-9099 www.Eberbachlabtools.com

## III.a Material Certification

Material can be provided for the following items upon request:

| PartNo                             | DESCRIPTION                      | Material  |
|------------------------------------|----------------------------------|-----------|
| • CP6155.003.001                   | DRIVE ENCLOSURE STANDOFF         | 304SS     |
| • CP6155.003.002                   | BASE                             | 304SS     |
| • CP6155.003.313                   | MOTOR MOUNT                      | 316SS     |
| <ul> <li>CP6155.003.224</li> </ul> | FRONT PANEL                      | 316SS     |
| • CP6155.003.225                   | BACK PANEL                       | 316SS     |
| • EP6155.S.201                     | ROUND SUPPORT                    | 303SS     |
| • EP6155.SS.220                    | TOP COVER                        | 316SS     |
| • EP6155.SS.221                    | RIGHT SIDE PANEL                 | 316SS     |
| • EP6155.SS.222                    | LEFT SIDE PANEL                  | 316SS     |
| • EP6155.SS.226                    | FRONT & BACK MOUNTING<br>BRACKET | 316SS     |
| • EP6155.SS.227                    | SIDE MOUNTING BRACKET            | 316SS     |
| • EP6155.025                       | LOWER BEARING LADDER             | A36 STEEL |
| • EP6155.S.015                     | BEARING RETAINER                 | 303SS     |
| • EP6155.S.101                     | ECCENTRIC BLOCK                  | 303SS     |
| • EP6155.S.213                     | SPINDLE                          | 17-4SS    |
| • EP6155.S.253                     | RETAINING WASHER                 | 303SS     |
| • EP6155.S.252                     | CRANK PIN                        | 17-4SS    |
| • EP6155.S.202                     | UTILITY BOX SUPPORT              | 303SS     |
| • EP6155.205                       | UPPER BEARING LADDER             | A36 STEEL |

<sup>\*</sup>Non-Eberbach manufactured parts will not have material certifications available

## III.b Electronics data

Technical data sheets can be provided upon request:

| PartNo    | DESCRIPTION                     |
|-----------|---------------------------------|
| • 6544.1  | <b>EMERGANCY STOP BUTTON</b>    |
| • 6288.2  | IP68 USB PORT PANEL MOUNT       |
| • 6288.7  | IP68 ETHERNET PANEL MOUNT       |
| • 6506    | ROCKER SWITCH DPST ON-OFF 16A   |
| • 6760.66 | IP65 5A CIRCUIT BREAKER SWITCH  |
| • 6760.46 | IP65 15A CIRCUIT BREAKER SWITCH |
| • 6269    | 125V AC PANEL MOUNT RECEPTACLE  |
| • 6272.5  | AC INLET C20                    |
| • 7107.1  | DRIVE                           |

### **SIEMENS**

#### Data sheet

#### 6AV2124-2DC01-0AX0

SIMATIC HMI KTP400 Comfort, Comfort Panel, Key/touch operation, 4" widescreen TFT display, 16 million colors, PROFINET interface, MPI/PROFIBUS DP interface, 4 MB configuration memory, Windows CE 6.0, configurable from WinCC Comfort V11



| General information                              |                |  |  |
|--|----------------|--|--|
| Product type designation                         | KTP400 Comfort |  |  |
| Display  |                |  |  |
| Design of display                                | TFT            |  |  |
| Screen diagonal                                  | 4.3 in         |  |  |
| Display width                                    | 95 mm          |  |  |
| Display height                                   | 53.8 mm        |  |  |
| Number of colors                                 | 16 777 216     |  |  |
| Resolution (pixels)                              |                |  |  |
| <ul> <li>Horizontal image resolution</li> </ul>  | 480 Pixel      |  |  |
| <ul> <li>Vertical image resolution</li> </ul>    | 272 Pixel      |  |  |
| Backlighting                                     |                |  |  |
| <ul> <li>MTBF backlighting (at 25 °C)</li> </ul> | 80 000 h       |  |  |
| Backlight dimmable                               | Yes; 0-100 %   |  |  |
| Control elements                                 |                |  |  |
| Keyboard fonts                                   |                |  |  |
| <ul> <li>Function keys</li> </ul>                |                |  |  |
| <ul> <li>Number of function keys</li> </ul>      | 4              |  |  |
|  |                |  |  |

6AV2124-2DC01-0AX0 Page 1/10

01/29/2018

Subject to change without notice © Copyright Siemens

## **SIEMENS**

#### Data sheet

#### 6ES7211-1HE40-0XB0

SIMATIC S7-1200, CPU 1211C, COMPACT CPU, DC/DC/RELAY, ONBOARD I/O: 6 DI 24V DC; 4 DO RELAY 2A; 2 AI 0 - 10V DC, POWER SUPPLY: DC 20.4 - 28.8 V DC, PROGRAM/DATA MEMORY: 50 KB



| General information                                     |                          |
|---|--------------------------|
| Product type designation                                | CPU 1211C DC/DC/Relay    |
| Firmware version  | V4.1                     |
| Engineering with  |                          |
| Programming package                                     | STEP 7 V13 SP1 or higher |
| Display   |                          |
| with display  | No                       |
| Supply voltage  |                          |
| Rated value (DC)  |                          |
| • 24 V DC   | Yes                      |
| permissible range, lower limit (DC)                     | 20.4 V                   |
| permissible range, upper limit (DC)                     | 28.8 V                   |
| Reverse polarity protection                             | Yes                      |
| Load voltage L+   |                          |
| Rated value (DC)  | 24 V                     |
| <ul> <li>permissible range, lower limit (DC)</li> </ul> | 20.4 V                   |
| • permissible range, upper limit (DC)                   | 28.8 V                   |

 6ES7211-1HE40-0XB0
 Changes preserved


 Page 1/8
 19.01.2016
 © Copyright Siemens AG

#### SGM7J Gear Motors

The SGM7J gear motor product family pairs SGM7J servo motors with high precision, low backlash inline planetary gear heads resulting in a portfolio of rotary actuators fit for a wide range of applications. The family of gear motors has been thoroughly tested and adheres to the high levels of quality and performance expected from Yaskawa.

The high precision gear heads offer a variety of application advantages:

- Quiet operation helical cut gears contribute toward reduced vibration and noise
- @ High precision a standard backlash of 5 arc-min make this gear head ideal for the most accurate applications
- 9 High rigidity and torque capacity achieved with a design which incorporates uncaged needle roller bearings
- ® Optimized adapter bushing minimizes inertia allowing for more output torque to be realized
- O No leakage through the seal high viscosity, anti-separation grease does not liquefy and does not migrate away from the gears
- Maintenance-free no need to replace the grease for the life of the unit. The reducer can be positioned in any orientation



#### Model Designations





| Code | Specification |
|------|---------------|
| 01   | 100 W         |
| 02   | 200 W         |
| 04   | 400 W         |
| 08   | 750 W         |
| 15   | 1.5 kW        |

| 3rd digit | Power | Supply | Voltage |
|-----------|-------|--------|---------|
|-----------|-------|--------|---------|

| Code | Specification            |
|------|--------------------------|
| Α    | 200 VAC battery type     |
| В    | 200 VAC batteryless type |
| D    | 400 VAC battery type     |

4th digit Brake Option

| Code  | Specification |
|-------|---------------|
| Blank | No brake      |
| С     | 24 V Brake    |
|       |               |

| 5th digit | Gear box | backlash |
|-----------|----------|----------|
|-----------|----------|----------|


| Code | Specification      |
|------|--------------------|
| VL   | 5 arc-min backlash |


6th digit Gear head frame size


|   | Code | Specification |
|---|------|---------------|
| 1 | 050  | 50 mm         |
|   | 070  | 70 mm         |
|   | 090  | 90 mm         |
|   | 120  | 120 mm        |
|   | 155  | 155 mm        |


| 7th digit    | 0   | D-ti-  |
|--------------|-----|--------|
| TO THE COLOR | GBB | Hallic |

| 1 | Code | Specification |
|---|------|---------------|
|   | 03   | 3:1 Ratio     |
|   | 05   | 5:1 Ratio     |
|   | 10   | 10:1 Ratio    |
|   | 25   | 25:1 Ratio    |
|   | 50   | 50:1 Ratio    |













## **III.d Wiring Diagrams**













## **III.e Material Information**









## **II.f** Spare Parts List

| PartNo           | DESCRIPTION    | QTY |
|------------------|----------------|-----|
| • CP6155.006.002 | TAILGATE CABLE | 2   |
| • 6668.2         | POWER CORD     | 1   |

## **IV.** Customer Signoff

| Name:      |  |  |
|------------|--|--|
|            |  |  |
| Signaturo: |  |  |
| Signature: |  |  |
|            |  |  |
| Date:      |  |  |