Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung

REC Heating Plate

F20700081, F20710081

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch

Caution, hot surface! / Attenzione, superficie calda! / Attention, surface chaude! / Prudencia, superficie caliente! /
Vorsicht, heiße Oberfläche!

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.

This unit must be used for laboratory applications only. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The unit must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The minimum fire point of flammable solution is 750 °C. Only small amounts (< 50 ml) of flammable liquid can be used with the device.

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The unit is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xT8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it. It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The REC heating plate with ceramic plate, is used in the laboratory for general use and for all those applications that require a precise regulation of the heating plate temperature.

The ceramic heating plate (1) means that higher temperatures can be obtained compared to the traditional aluminium heating plate and it takes less time to heat the sample. Ceramic is an inert and very hard material and is resistant to almost any type of chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The instrument is turned on using the ON/OFF button (2). The temperature (4) of the heating plate is adjusted using the relative knob on the front panel, the value is indicated on the digital display (3).

NOTE: If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

Fig. 1

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- REC Heating Plate
- Power supply cord
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. The equipment shall be connected to a power socket with protective earth connection, using only the power supply cord provided with the instrument.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

NOTE: the mains cable must remain far away from the hot plate.

2.2 Start-up

Rotate the temperature knob completely to the left. Place the flask containing the sample. Then, set the temperature by turning the knob.

3. Operating controls

3.1 Turning on and regulating

Connect the unit to mains. In the next 5 seconds, displays shows the software version and the restarting mode (A or b, chapter "Setting mode") and turn it on using the on-off button. The green led indicates that the instrument is On.

To start the heating function turn the heating knob "Heating °C" on the front panel. Temperatures of from 5 to 550 °C can be selected at 5 °C intervals. The temperature setting is shown on the display. When heating is turned Off (Heating knob OFF) and the heating plate temperature exceeds 50 °C, the message "Hot" flashes on the display until the temperature of the heating plate falls to below 50 °C (or until the instrument is turned off using the ON/OFF).

NOTE: this warning is not active if the instrument is not connected to the power supply.

NOTE: in case of black-out, once the power is back the device will restart in the restarting mode set (OPt tyP), selected in chapter "Setting mode".

The knob on the front panel is easily accessible and is well away from the heating plate in order to ensure maximum operator safety as well as safeguarding the electronic components inside the unit.

3.2 Setting mode

Turn the knob twice from the minimum to the maximum during the visualization of the software version in order to enter into the setting mode.

Once there, by turning the knob, it is possible to pass from one of the following parameter to the next ones:

Parameter shown	Default value	Range	Description	
FS	tCJ	550	0 ÷ 550°C	Maximum temperature settable
OPt	tyP	b	A – b	(*) Restarting mode

Press on-off button in order to pass to the parameter you want to modify. The display will show the value saved for 3 seconds, and then the new value according to the position of the knob.

The new value is automatically saved after 4 seconds, if during this time the knob is not moved. "End" is shown and then the system goes back to the main menu.

Once the new values are set, switch off and restart the instrument.

(*) When **A** is selected: in case of black-out, once the power is back the device will remain off

When **b** is selected: in case of black-out, once the power is back the device will automatically restart with the latest settings

3.3 Error messages

If an error message appears on the display, please contact VELD Scientifica's technical service department.

When the display shows an error message, heating stops automatically.

NOTE: To remove the error message disconnect the instrument from the power supply.

AL1 Thermocouple is out of range

AL2 Excessive heating time

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized). It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Power supply	230 V / 50-60 Hz or 115 V / 50-60 Hz
Dimensions WxHxD	203x93.5x344 mm (8x3.68x13.5 in)
Weight	3.5 Kg (7.7 lbs)
Overall power	780 W
Dimension of the heating plate	180 x 180 mm
Programmable temperature range	5 – 550 °C
Type of temperature control	Digital
Overtemperature protection	Yes
Noise level	<< 80 dBa
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Overtoltage category	II
Pollution degree CEI EN61010-1	2
Max altitude	4000 m

6. Accessories / Spare parts

A00001071 Hemispheric bowl for 250ml flasks
 A00001072 Hemispheric bowl for 500ml flasks
 A00001073 Hemispheric bowl for 1000ml flasks
 A00001069 Support rod

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-2-010/2003

EN 61326-1/2013

and satisfies the essential requirements of the following directives:

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Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl