




Product Line-up

Amber [®] LiSi-POZ		Dimensions (mm)	pcs / Pack
	R10	Ø12.7 × T 10	5 Ingots
	R15	Ø12.7 × T 15	5 Ingots
	R20	Ø12.7 × T 20	3 Ingots

Available Shades



Indications with zirconia frameworks

- Crowns
- 3-unit anterior and premolar bridges
- Long-span and curved bridges
- Cantilever bridges
- Maryland bridges
- Implant supported crowns and bridges

Pressing Schedules

	Translucency	Size	Shade	Investment Ring	Start Temp.	Heating Rate	Max Temp.	Holding Time	Vacuum On	Vacuum Off
Amber LiSi-POZ	HT	R10 / R15	A1, A2, A3, A3.5, W4	Small (100g)	700°C	45°C/min	915°C	15 Min	700°C	915°C
	LT									
	HT	R20		Large (200g)				30 Min		
	LT									

※ Note : 1. There may be a difference between the displayed temperature and the real temperature of each furnace. When you use the Amber ingots, please verify the above standard schedule is suitable for your press furnace. If it is not, please try to find the optimum temperature through the following process.

- 1) If there are some traces of tiny bubble on the surface of the restoration ⇒ Please reduce the maximum temperature by 5~10°C and try pressing again.
- 2) If the marginal area of the restoration is not formed completely ⇒ Please increase the maximum temperature by 5~10°C and try pressing again.


Lithium Disilicate Press on Zirconia

Amber[®] LiSi-POZ

Most Innovative and Exciting

Easy Esthetics & Superior Strength

Opening Up New Era of Dental Restorations



More Lifelike

Amber® LiSi-POZ veneering has similar translucency to the enamel layer of natural teeth and the translucency of zirconia framework is similar to that of dentin of natural teeth. The high esthetics of Amber® LiSi-POZ enables it to replace a damaged natural teeth perfectly.




Restoration courtesy of CDT, Young Soo Kim

Procedures





Superior Strength

The tensile bond strength between zirconia frame work and Amber® LiSi-POZ is over 45 MPa after pressing.



Amber LiSi-POZ offers three times higher flexural strength than conventional veneering materials for Zirconia. After pressing the flexural strength is over 380 MPa.

The fatigue fracture strength of restorations made from Zirconia framework and Amber® LiSi-POZ is as high as monolithic zirconia crown.



1) 2) Cited from Pusan National University Dental Hospital, International Journal of Prosthodontics. Fracture strength after fatigue loading of lithium disilicate pressed zirconia crowns, 2015.
* Not a registered trademark of HASS Corp.