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DE

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ab Seite 3

EN

## **Operating Instructions**

starting on page 20

FR

## **Mode d'emploi**

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NL

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vanaf pagina 54

IT

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da pagina 71

## Inhalt

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## EMC Information

The appliance is complied with the following emission and immunity standard, and test level:

Emissions Test Standard	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Comply

Immunity Test standard	IEC 60601-1-2 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air
Electrical fast transient/burst IEC 61000-4-4	± 2kV, for power supply lines ± 1kV signal input/output 100kHz repetition frequency	± 2kV, for power supply lines ± 1kV signal input/output 100kHz repetition frequency
Surge IEC 61000-4-5	± 0.5kV, ± 1kV differential lines ± 0.5kV, ± 1kV, ± 2kV common mode	± 0.5kV, ± 1kV differential lines ± 0.5kV, ± 1kV, ± 2kV common mode
Voltage dips, interruptions and variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz

Note: UT is the AC mains voltage prior to application of the test level.

## EMC Information

### Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specification for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 – 390	TETRA 400	Pulse modulation 18Hz	1,8	0.3	27
	450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
	870						
	930						
	1720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2 400 - 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
	5500						
5785							

## Contents

Intended Use _____	22
Safety Instructions _____	22
The VenenWalker® » basic « _____	23
The History of IPC (Intermittent Pneumatic Compression) _____	24
The IPC principle _____	25
Benefits of the IPC _____	25
Applications _____	26
Items Supplied and Device Overview _____	27
Putting Into Operation _____	28
Use _____	29
Cleaning and Care _____	30
Additional Information on the VenenWalker® » basic « _____	31
Technical Data _____	32
Troubleshooting _____	33
Disposal _____	33
Customer Service _____	33
Other Information _____	34
EMC Information _____	34

## IMPORTANT NOTICES! PLEASE KEEP FOR REFERENCE!

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device.

The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

The illustrations in these instructions may differ slightly from the actual device.

We reserve the right to modify the device, packaging or enclosed documentation at any time in conjunction with our policy of continuous development.

## Dear Customer,

Congratulations on deciding to purchase the VenenWalker® » basic «!

The VenenWalker® » basic « is a certified medical device which can help you prevent venous diseases and reduce the risk of spider veins and varicose veins. If you have a medical indication, please consult your GP before use!

If you have any questions regarding the VenenWalker® » basic «, contact our website: **www.venenwalker.com**

or contact us from Mon – Fri by calling the following service number:

+49 (0) 38851 314337\*

\* Calls to German landlines are subject to your provider's charges.

We wish you every success when using the VenenWalker® » basic «.

## Meaning of the Symbols

Read operating instructions before use! This symbol appears in blue on the device.

Safety instructions are identified with this symbol. Please read them carefully and comply with the safety instructions in order to prevent personal injury and damage to property.

Tips, notes and recommendations are marked with this symbol.

Symbol for devices classified as "Type BF" (level of protection against electric shock)

**IP22**

Protected from touch by fingers and objects greater than 12mm and from water spray less than 15° from vertical.

Protection class II

Energy efficiency class V

Polarity of the coaxial connector (mains adapter)

Manufacturer

EU Representative

Serial number (Date of production / consecutive number)

With the CE mark, the manufacturer confirms that the product meets the requirements of the applicable EC directive and complies with the "general requirements" stipulated therein. The number 0123 alludes to the notified body (certifier). Together these certify that the German Medical Devices Act (MPG) has been observed by the distributor.

## Intended Use

The VenenWalker® »basic« may only be used for decongestion of the extremities in the foot and lower leg area. If you have a medical indication, please consult your doctor before use.

VenenWalker basic is a IPC (Intermittent Pneumatic Compression) device for private use for the treatment of the following symptoms:

1. Leg swelling caused by lack of movement
2. Leg swelling prophylaxis
3. Venous diseases such as varicose veins and spider veins
4. Lipoedemas
5. Swelling following a pulled muscle or sprain

The leg cuff is treated as the applied part.

## Safety Instructions

### The Proper Handling of Batteries

- Before inserting the batteries, check whether the contacts in the device and on the batteries are clean; clean them if necessary.
- Use only the battery type specified in the technical data.
- Always replace all the batteries at the same time. Do not use batteries of different types or brands, or batteries with different capacities. When replacing batteries, note the correct polarity (+/-).
- Take the batteries out of the device when they are discharged or if the device is not going to be used for a prolonged period. This will prevent any damage that might occur if a battery leaks.
- Remove dead batteries from the device immediately. Clean the contacts before you insert a new battery.
- The batteries must not be recharged or reactivated by other means, dismantled, thrown into a fire, immersed in fluids, or short-circuited.

### Danger of Injury

- Keep the device and packaging material away from children and animals. There is a risk of injury and suffocation!
- This device should not be used by people (including children) with impaired sensory or mental faculties or who do not have adequate experience and/or knowledge of such devices, unless they are supervised by a person responsible for their safety,

or have received instructions from such a person on how to use the device. Children must be supervised to ensure that they do not play with the device.

- The device should not be used if it is damaged.
- Use only accessories that are supplied with the device or recommended by the original manufacturer!
- Use the device only in dry indoor rooms.
- Use the device away from moist and wet areas.
- The main plug is used as disconnect to the mains supply, do not position the ME EQUIPMENT so that it is difficult to operate the disconnection device.
- The device contains no user serviceable parts.
- No modification of this equipment is allowed.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

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The comprehensive massage device for the prevention of venous complaints - at a very low price.

Too much sitting, too much standing, poor nutrition: According to consistent estimates, 20 % of all Europeans have venous leg complaints and more than half is affected by advanced vein problems. This global lifestyle disease starts off harmlessly with a tingling sensation in the legs. The legs then become heavy, swollen and tired and spider veins and varicose veins appear. The entire cardiovascular system is in danger and there is a risk of ulcers, thrombosis, a pulmonary embolism and an inability to work.

How does a venous disorder occur? Veins collect the used blood in the body's tissue and transport it back to the heart. Leg muscles and venous valves tense up and relax, thus acting like a pump to transport the blood forwards into the veins. However, if a person stands or sits for a long time, or if excess strain due to surplus weight leads to congestion, the blood cannot flow away properly and the venous walls become overstretched. Swollen legs are the result. The swollen leg forms an oedema due to the fact that liquid and protein pass through the slackened venous wall and into the tissue. A doctor should be consulted as a matter of urgency at the latest when such a venous congestion becomes a continual complaint. Medicines are not able to prevent the congestion, nor reverse the abnormal slackening of the veins. Surgical stockings are a tried and tested means of supporting the weakened tissue from the exterior and preventing swollen legs. However, they should be worn at all times. "Vein walking", a special type of sport which is very popular in the USA, specifically stimulates the circulation in the veins and trains all important muscles. Again, this sport is only effective if practised consistently.

Special massages using alternating-compression devices have proven to be particularly effective and time-saving in the treatment of venous complaints. The VenenWalker®

» basic « is the first ever comprehensively effective massage device which can efficiently prevent venous complaints - and this at a surprisingly low price.

The VenenWalker® » basic « comprises a control device and two double-walled leg cuffs which also cover the entire feet including the soles. The legs are lightly massaged by the VenenWalker® » basic « using an inflating and deflating air cushion which simulates the leg muscle pumping effect experienced during walking exercises. Venous reflux is improved, collected fluid is mobilised and can be flushed out.

The VenenWalker® » basic « was developed in close collaboration with medical specialists (phlebotomists), is easy to use and, if used properly, is absolutely safe and without any undesirable side effects. The entire kit contains the massage device and two leg cuffs including this detailed and easy-to-understand user operating instructions.

### A Specialist Comments

Dr. med. Frank Sömmmer, general practitioner and phlebotomist (doctor with specialist knowledge in venous disorders), with many years of experience in the treatment and prevention of venous disorders and Medical Director of the Institute for Compression Therapy in Lubeck, had the following to say about the VenenWalker®:

*“The VenenWalker® » basic « with short leg cuffs is very easy to use and can lead to a good increase in the venous return movement from the lower legs as it stimulates the important calf muscle pump.”*

## The History of IPC (Intermittent Pneumatic Compression)

The principle of intermittent pneumatic compression can be traced back to “quicksilver massage”. In 1899, a report was already circulating about a special treatment for lymphoedemas in the upper extremities (arms). Here, the extremities were immersed in an iron cylinder filled with quicksilver.

Back in 1917, pressure increases were used in a pneumatic chamber to “expel stagnant vein blood”. In the year 1929, “...intermittently inflating and deflating air cushions” were used for the first time in the treatment in bedridden patients with trophic ulcers. At the University Clinic of Leipzig, chronic venous circulation disorders with indurations (hardenings) and ulcers have been treated with the rhythmic alteration of the arterial blood flow since the beginning of the 1930s.

At around the same time, attempts were also being made in the USA to treat the consequences of artery occlusive diseases using intermittent compression.

In parallel with an increasing clinical use of the IPC, equivalent devices were further developed on a technical level right up to today's multi-chamber system for the treatment of severe lymphoedema.

## The IPC principle

Intermittent (alternating) compression, which is not generated through manual massage but rather by a device, is known in medicine as intermittent pneumatic compression, or IPC. It uses the principle of pneumatic alternating pressures. It is prescribed by doctors for the treatment of venous or lymphatic diseases. When used to prevent thrombosis and for sports injuries, it does not carry the risks found in medication-based treatment.

The IPC simulates the function of the natural muscle pumps found in the human body. Air is pumped with alternating pressure into double-walled leg cuffs for example. The pressure is increased and decreased in defined time segments. It should be possible to adjust the pressure level to suit the user and the pressure must never exceed 100 mmHg. The treatment duration and frequency depend on the person's condition when used for cosmetic purposes, or as recommended by the doctor when used for medical reasons. However, the treatment should last for at least 20 minutes. Any medical use should, of course, be supervised by a doctor.

The shape and type of cuff used is important. Due to the proven importance of the draining of the network of veins in the sole, the cuff should cover both the entire leg and the foot.

The rhythmic increase and decrease in pressure on the tissue leads to water and protein being pumped out of the interstitial tissue. If the patient has a medical indication, a compression dressing should be applied or a surgical stocking should be worn after using the IPC as the residual protein causes blood serum to flow from the capillaries into the tissue after an IPC treatment session. This measure prevents this flow occurring. These measures do not need to be taken if you simply remain lying down after treatment. When using the equipment at home it therefore makes sense to use it before bedtime.

## Benefits of the IPC

- Improvement of venous and lymphatic macrocirculation through the passive activation of the natural muscle pump
- Improvement of blood flow characteristics
- Positive flow-mechanical effects
- Considerable boost to the wound healing process
- Can promote the healing of lower-leg ulcers
- Favourable influence on lymphatic fistulas
- Carries a lower risk than other methods
- Pain-free, comfortable, uncomplicated and excellent tolerability

## Applications

### For Cosmetic Purposes

- Leg swelling caused by lack of movement, for example due standing for long periods of time in the case of sales personnel, due to sitting for longer periods for office work, or due to long-haul travel.
- The unnatural reproduction of connective tissue (fibrositis)
- Swollen legs through bedriddenness or pregnancy\*
- For varicose veins
- For muscle ache

### For medical purposes

- To treat lymphedema
- To treat vein thromboembolic
- To prevent embolisms due to bedriddenness
- To prevent venous complaints such as varicose veins
- To reduce swelling and heavy legs

## Contraindications

Do not use the IPC - or only use it only after consulting your doctor - if you have:

- Acute inflammations
- Open injuries
- Cardiac output problems with fluid retention in your legs
- Advanced arterial disease
- Suspected thrombosis
- Malignant tumors
- Pulmonary oedemas
- High blood pressure


## Items Supplied and Device Overview

- 1 Leg cuffs (2x)
- 2 Control unit
- 3 Storage bag
- 4 Air tubes (2x)
- 5 Mains adapter

Not shown:  
4x 1.5 volt demo batteries (AA type)

## Overview of Control Unit

- 6 Air tube connection for inflating the leg cuffs
- 7 Pressure adjuster
- 8 On/Off switch
- 9 On/Off control lamp
- 10 Air tube connection for deflating the leg cuffs

## Putting Into Operation

The device can be operated either with batteries or with the mains adapter supplied.

- Before use, insert the batteries or connect the mains adapter.

### Inserting/Changing Batteries

The battery compartment is on the rear of the control unit. You will need four 1.5 volt AA batteries. Use only standard batteries, never rechargeable batteries!

The demo batteries supplied only have enough power to operate the device for approx. 20 minutes.

1. Open the battery compartment.
2. Insert four 1.5 volt AA type batteries. Ensure the correct polarity (+/-) when inserting the batteries.
3. Close the battery compartment.

Your VenenWalker® » basic « is now ready for use!

### Connecting the Mains Adapter

1. Connect the mains adapter to the control unit. The connection is located below the area for opening the battery compartment.
2. Connect the mains adapter to a properly installed socket which corresponds to the technical data of the device (see the “Technical Data” chapter).

Your VenenWalker® » basic « is now ready for use!


## Use

### CAUTION!

- To prevent damaging the air chambers, do not bend the leg cuffs! Carefully fold them after use before storing.
- Do not open the zips on the leg cuffs during use as otherwise the cuffs could get damaged.
- Please make sure that the device is always free and is not covered by a blanket or other objects.
- If the appliance does not pump the air in cuff normally or it is not working in alternating inflating, do not use it or contact customer service.

The device can be operated both with batteries and with the supplied mains adapter. To connect the mains adapter or insert the batteries, read the chapter "Putting into Operation".

1. Put on the leg cuffs with the connections for the air tubes pointing outwards.
2. Connect the air tubes to the connections on the leg cuffs. Make sure that the tubes are connected properly and cannot slip out.
3. Connect the plug for the air tubes to the connection on the control unit for inflating the leg cuffs.
4. Set the pressure adjuster to its maximum setting.
5. Switch on the device by pressing the On/Off button. The On/Off control lamp lights up as soon as the device is turned on.
6. The leg cuffs are now alternately pumped up. If the pressure on your leg is too strong or if it feels unpleasant, set the pressure adjuster to a lower setting. The ideal pressure setting should feel strong but not unpleasant.

## How it works

Your legs are massaged by the rhythmic and alternating inflating and deflating air cushions. This massaging of your legs has the same effect on your veins and muscles as a slightly strenuous walk. For the best results, your legs should be positioned horizontally or raised during use so that the blood can circulate unimpeded.

We recommend using the device for approximately 20 – 30 minutes.

Please note that your legs are massaged one after the other by the inflating and deflating air cushions, not simultaneously. The pressure is built up firstly in one leg cuff and then in the other. This is not a sign that the device is faulty!

## After Use

Once you have finished use, the air must be removed from the leg cuffs.

1. Set the pressure control to the highest level.
  2. Switch off the device by pressing the On/Off button. The control lamp goes out.
  3. Unplug the air-tube connection plug from the top socket on the control unit.
- 
4. Now connect the plug for the air tubes to the connection on the side of the control unit.
  5. Switch on the device by pressing the On/Off button. The On/Off control lamp lights up as soon as the device is turned on. The air is now pumped out of the leg cuffs.
  6. Once the air has been pumped out, open the zip and remove the leg cuffs.
  7. Switch off the device by pressing the On/Off button.
  8. If you used the mains adapter, unplug it from both the mains socket and the device.

## Cleaning and Care

### CAUTION!

- Keep the device in a place which is protected from direct sunlight and extreme temperature fluctuations.
  - Do not use any caustic or abrasive cleaning agents to clean the device. These could damage the device.
- 
- Before cleaning, always remove the air tubes from the leg cuffs and from the control unit.
  - The leg cuffs can be cleaned with a slightly damp cloth and some household detergent. If the leg cuffs are used for more than one person and if they come into direct contact with the skin, you should disinfect the cuffs with a standard disinfectant to avoid cross contamination. **It is imperative that you make sure that no moisture enters the air inlet and outlet openings.**
  - When used for a long time, the ends of the tubes may expand and air could leak out. In such case, you can cut off a short piece of the ends.
  - Only lightly fold the cuffs together. Bending the cuffs with too much force can lead to cracks.
  - We recommend that you only clean the control unit when it is absolutely necessary! Clean the control unit with a slightly damp cloth and, if necessary, a little household detergent.

## Additional Information on the VenenWalker® »basic«

### Benefits

- Problem-free and low-risk option of supporting decongestion treatment with compression stocking dressings
- Easy to use
- Cuffs stretching from the foot to the knee are rhythmically and alternately inflated and deflated
- The air is automatically removed from the cuffs when the plugs are switched over
- The leg is gently massaged
- The function of the leg's muscle pumps is simulated
- Not only for cosmetic problems or sports injuries
- Also for older or overweight patients
- Can be used comfortably at home on a daily basis
- Can also be used on holiday, after sitting for a long time
- Continuous acceleration of the venous outward flow for several hours
- Swellings, pains and the feeling of heaviness in the legs are reduced

### Use at home for the prevention of

- Sports injuries
- Muscle ache
- Contusions
- Sprains
- Haematomas
- Spider veins
- Varicose veins
- Simple venous congestion
- Swelling
- Feeling of heaviness
- Slight disturbances in the arterial blood flow

#### CAUTION!

If you have any serious illnesses, please consult your GP before use!

### The intermittent compression is well proven for

Venous, arterial, lymphatic and connective-tissue diseases

### What are the results?

- Leg swelling can be reduced.
- Congestion pains can subside.
- The pressure in the varicose veins/spider veins can be reduced, the width can be reduced.
- The function of the venous valves can be improved.
- In the case of lymphedema, the VenenWalker® »basic« can be used as a therapy supplement.

## Technical Data

Product number: 06019

### Device

Product name: VenenWalker® » basic «

Model number: HNS018839

### Voltage supply

- In battery mode: 6 V DC (4x 1.5 volt batteries, AA type)
- In mains power mode: See below “Main adapter”

### Cuff pressure

- Operation: Maximum 120 mmHg
- Storage and transportation: no pressure

### Temperature range

- Operation: 15 – 40 °C
- Storage and transportation: 0 – 40 °C

### Air humidity

- Operation: 30 – 75 %
- Storage and transportation: 10 – 80 %

### Atmospheric pressure

- Operation: 700 – 1060 hPa
- Storage and transportation: 700 – 1060 hPa

Dimensions: 80 mm (l) x 45 mm (w) x 190 mm (h)

Weight: approx. 460 g

Degree of protection against electric shock: Type BF 

Degree of protection against the ingress of water, particulate matter: IP22

Expected service life of equipment: 5 years

### Mains adapter (Accessory)

Model number: UES06WV-060100SPA/UES06WB-060100SPA

Input: 100-240 V~, 50/60 Hz, 200 mA

Output: 6.0 V DC, 1.0 A

Protection class: II 

Efficiency level: VI

Polarity: 

Maximum length of cable: 2000 mm

ID of operating instructions: Z 06019 M GM V1.1 0321 (02 – month; 21 – year)

This appliance meets the requirement of Medical Device Directive 93/42/EEC for medical appliances.

## Troubleshooting

If the device stops working properly, check whether you are able to correct the problem yourself based on the tips below. If this is not possible, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

Check the following:

- Battery mode: Are the batteries inserted wrong or empty? Try using new batteries.
- Mains power mode:
  - Is the mains adapter not inserted in the plug socket correctly or is its cable not connected to the device correctly?
  - Is the plug socket defective? Try a different plug socket.
  - Check the fuse of your mains connection.

## Disposal

Dispose of the packaging material in an environmentally friendly manner and take it to a recycling centre.

Batteries should be removed before the device is disposed of and should be disposed of separately from the device at suitable collection points.

This device is governed by the European Directive 2012/19/ EU on waste electrical and electronic equipment (WEEE). Do not dispose of the device as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

## Customer Service

If you have any questions on the VenenWalker® »basic«, contact us via our website [www.venenwalker.com](http://www.venenwalker.com) or by telephone on

+49 (0) 38851 314337\* (Mon – Fri, 8 am – 12 midday)

\* Calls to German landlines are subject to your provider's charges.

## Other Information

### The VenenWalker® » basic « is manufactured by:

Zero-Plus International Limited  
 Room 1004, 10/F.  
 71-75 container port road  
 Kwai Chung  
 N.T., Hong Kong

### European Authorised Representative and Importer:

Globalmind Consumer Electronics GmbH  
 Ernst-Mantius-Str. 11  
 21029 Hamburg  
 Germany

## EMC Information

- **Warning:** Use of this appliance adjacent to or stacked with other device should be avoided because it could result in improper operation. If such use is necessary, this appliance and the other device should be observed to verify that they are operating normally. If the appliance does not light on or flickering, it means that the appliance is not operating normally.
- **Warning:** Do not change any cable without authorization of the manufacturer. Use of accessories and cables other than those specified or provided by the manufacturer of this appliance could result in increased electromagnetic emissions or decreased electromagnetic immunity of this appliance and result in improper operation.
- **Warning:** Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.
- **Warning:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of this appliance (HNS018839), including power cable. Otherwise, degradation of the performance of this equipment could result.
- If the appliance does not pump the air in cuff normally or it is not working in alternating inflating and deflating due to the electromagnetic disturbances, do not use it and try to place the appliance in other location.
- For maintaining the basic safety and essential performance with regard to electromagnetic disturbances for the expected service life, please precaution to take if the use location is near (e.g. less than 1,5 km from) from AM, FM or TV broadcast antennas.

## EMC Information

The appliance is complied with the following emission and immunity standard, and test level:

Emissions Test Standard	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Comply

Immunity Test standard	IEC 60601-1-2 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air
Electrical fast transient/ burst IEC 61000-4-4	± 2kV, for power supply lines ± 1kV signal input/output 100kHz repetition frequency	± 2kV, for power supply lines ± 1kV signal input/output 100kHz repetition frequency
Surge IEC 61000-4-5	± 0.5kV, ± 1kV differential lines ± 0.5kV, ± 1kV, ± 2kV common mode	± 0.5kV, ± 1kV differential lines ± 0.5kV, ± 1kV, ± 2kV common mode
Voltage dips, interruptions and variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	3 V 0.15MHz – 80 MHz 6V in ISM and amateur radio bands between 0.15MHz and 80 MHz 80 % AM at 1 kHz	3 V 0.15MHz – 80 MHz 6V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	10V/m 80MHz – 2.7 GHz 80 % AM at 1 kHz	10V/m 80MHz – 2.7 GHz 80 % AM at 1 kHz

Note: UT is the AC mains voltage prior to application of the test level.

## EMC Information

### Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specification for EN-CLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 - 390	TETRA 400	Pulse modulation 18Hz	1,8	0.3	27
	450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1kHz sine	2	0.3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation 217Hz	0,2	0.3	9
	745						
	780						
	810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
	870						
	930						
	1720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
	1845						
	1970						
	2450	2 400 - 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
	5240	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation 217Hz	0,2	0.3	9
	5500						
5785							

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