




As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation

1 Remove the Smart-Tip[™] probe from test unit.



2 Sample suspected substance.

Powder

Touch the Smart-Tip[™] probe to the sample.




Tablet

Firmly press Smart-Tip[™] probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.



Test Preparation


3 Lock Smart-Tip[™] probe into test unit. Press firmly.



STOP
Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.

Identification

4 With thumb and forefinger, squeeze chamber 1 until capsule breaks.




5 Gently shake test unit for 2-3 seconds.




Identification

6 With thumb and forefinger, squeeze chamber 2 until capsule breaks.



7 Gently shake test unit for 2-3 seconds.




A blue, blue-violet or purple color will develop. Proceed to Step 8.

Identification

8 With thumb and forefinger, squeeze chamber 3 until capsule breaks.



9 Vigorously shake test unit for 20 seconds.



Hold test unit steady and keep upright, until 2 layers are visibly separated. This may take up to 2 minutes. The process may occur sooner.

Results



The appearance of a violet color in the upper layer of liquid within two minutes indicates the presence of Ephedrine or Pseudo-Ephedrine.



If the upper layer of liquid does not turn a violet color within two minutes, Ephedrine is not present.



A blue color does not indicate a positive result.



Disposal

The test kit is safe to dispose of in the trash.

NET 01-21

IMPORTANT INFORMATION

Important Information

- IDenta Touch&Know[™] Discreet Drug Test Kits** are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
- All results are presumptive.** If necessary, further testing by a competent drug testing lab can provide additional verification.
- A color change indicates a positive result for the listed drugs only if it occurs within two minutes.**

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After five minutes, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.