



**Read Before Use**

As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation

**1** Remove the Smart-Tip<sup>™</sup> probe from test unit.




**2** Sample suspected substance.

Liquid

Place the absorbent paper of Smart-Tip<sup>™</sup> probe into the liquid. Make sure the paper is wet.



**3** Lock Smart-Tip<sup>™</sup> probe into test unit. Press firmly.




**STOP**

**Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.**

Identification

**4** With **thumb and forefinger**, squeeze chamber 1 until capsule breaks.




**5** Gently shake test unit for 2-3 seconds.




Identification

**6** With **thumb and forefinger**, squeeze chamber 2 until capsule breaks.



**7** Gently shake test unit for 2-3 seconds.



Results



The appearance of a light orange color within one minute indicates the presence of **GBL** ( $\gamma$ -Butyrolactone).



If a light orange color does **not** develop within one minute or the liquid remains a yellow color, GBL is **not** present.



**Disposal**

The test kit is safe to dispose of in the trash.

NET 07-21

### IMPORTANT INFORMATION

- IDenta Touch&Know<sup>™</sup> Discreet Drug Test Kits** are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
- All results are presumptive.** If necessary, further testing by a competent drug testing lab can provide additional verification.
- A color change indicates a positive result for the listed drugs only if it occurs within the first minute.**

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.