



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip [™] probe from test unit.	Sample suspected substance.	
		2	Powder	Tablet
		Touch the Smart-Tip [™] probe to the sample.		Firmly press Smart-Tip [™] probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.

Test Preparation	3	Lock Smart-Tip [™] probe into test unit. Press firmly.	STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.	4	With thumb and forefinger, squeeze chamber 1 until capsule breaks.	5	Gently shake test unit for 20 seconds.
		Identification					

Results		The appearance of a red-purple color within one minute indicates the presence of Flunitrazepam .
		The appearance of a yellow color within one minute indicates the presence of Nitrazepam, Clonazepam or Nimetazepam .
		The appearance of a green color within one minute indicates the presence of Lorazepam or Temazepam .
		If none of the above colors develop within one minute, the above listed substances are not present.

Disposal The test kit is safe to dispose of in the trash.

IMPORTANT INFORMATION

- IDenta Touch&Know[™] Discreet Drug Test Kits** are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
- All results are presumptive.** If necessary, further testing by a competent drug testing lab can provide additional verification.
- A color change indicates a positive result for the listed drugs **only if it occurs within the first minute.**

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.