




As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation

1

Remove the Smart-Tip[™] probe from test unit.



2

Sample suspected substance.

Powder


Touch the Smart-Tip[™] probe to sample in different places on the sample. Make sure the tip is completely covered.



Test Preparation

3

Lock Smart-Tip[™] probe into test unit. Press firmly.




STOP

Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.

Identification

4

With thumb and forefinger, squeeze chamber 1 until capsule breaks.



5

Gently shake test unit for 2-3 seconds.



Results



The appearance of a red-purple color within one minute indicates the presence of **Fentanyl**.



The appearance of a purple color within one minute indicates the presence of an opiate such as **Heroin, Morphine** or **Codeine**.
* This does not exclude the possibility that Fentanyl may be present.



The appearance of a reddish-orange color within one minute indicates the presence of **Amphetamine** or **Methamphetamine**.
* This does not exclude the possibility that Fentanyl may be present.



If **none** of the above colors develop within one minute, the above listed drugs are **not** present in the tested sample.



Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.



Dilution

6

With thumb and forefinger, squeeze chamber 2 until capsule breaks.



7

Gently shake test unit for 2-3 seconds.




Wait 5 minutes.
The test window will warm up. This is normal.

Dilution

8

With thumb and forefinger, squeeze chamber 3 until capsule breaks.



9

Gently shake test unit for 2-3 seconds.



The test window will warm up. This is normal.



Disposal

After Dilution, the test kit is safe to dispose of in the trash.

NET 01-21

IMPORTANT INFORMATION

1. IDenta Touch&Know[™] Discreet Drug Test Kits are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
2. All results are **presumptive**. If necessary, further testing by a competent drug testing lab can provide additional verification.
3. A color change indicates a positive result for the listed drugs **only if it occurs within the first minute**.

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.