



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip™ probe from test unit.	2 Sample suspected substance.	
			Plant Matter Touch the Smart-Tip™ probe to sample. Make sure the tip is completely covered.	White Powder Touch the Smart-Tip™ probe to the sample.

Test Preparation	3	Lock Smart-Tip™ probe into test unit. Press firmly.	STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.	4	With thumb and forefinger, squeeze chamber 1 until capsule breaks.	5	Gently shake test unit for 2-3 seconds.
				Identification			

Results



The appearance of a yellow color within one minute indicates the presence of a **Synthetic Cannabinoid** (street names – **K2, Spice**), such as Naphthoylindoles, Benzoylindoles, Naphthoylpyrroles.



If **no** yellow color develops within one minute, the above listed drugs are **not** present.

*The test will **not** react to the following substances: Classical cannabinoids, Nonclassical cannabinoids, Phenylacetylindoles.



Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.



Dilution	6	With thumb and forefinger, squeeze chamber 2 until capsule breaks.	7	Gently shake test unit for 2-3 seconds.
				Wait 5 minutes. The test window will warm up. This is normal.

Dilution	8	With thumb and forefinger, squeeze chamber 3 until capsule breaks.	9	Gently shake test unit for 2-3 seconds.
				The test window will warm up. This is normal.

Disposal After Dilution, the test kit is safe to dispose of in the trash. NET 07-21

IMPORTANT INFORMATION

1. IDenta Touch&Know™ Discreet Drug Test Kits are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
2. All results are **presumptive**. If necessary, further testing by a competent drug testing lab can provide additional verification.
3. A color change indicates a positive result for the listed drugs **only if it occurs within the first minute**.

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.