



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip [™] probe from test unit.	2		Sample suspected substance.			
			Powder	Tablet	<p>Touch the Smart-Tip[™] probe to the sample.</p>	<p>Firmly press Smart-Tip[™] probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.</p>		
Test Preparation	3	Lock Smart-Tip [™] probe into test unit. Press firmly.	<p>STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.</p>	Identification	4	With thumb and forefinger, squeeze chamber 1 until capsule breaks.	5	Gently shake test unit for 2-3 seconds.

Results		The appearance of a yellow color within one minute indicates the presence of a Cathinone , or an analogue such as Buphedrone, Mephedrone, Methylcathinone, or N,N-Dimethylcathinone.
		If no yellow color develops with one minute, the above listed drugs are not present.



Proceed to Dilution (steps 6-9), regardless of whether or not a color appeared.



Dilution	6	With thumb and forefinger, squeeze chamber 2 until capsule breaks.	7	Gently shake test unit for 2-3 seconds.
				<p>Wait 5 minutes. The test window will warm up. This is normal.</p>

Dilution	8	With thumb and forefinger, squeeze chamber 3 until capsule breaks.	9	Gently shake test unit for 2-3 seconds.
				<p>The test window will warm up. This is normal.</p>

Disposal After Dilution, the test kit is safe to dispose of in the trash. NET 07-21

IMPORTANT INFORMATION

- IDenta Touch&Know[™] Discreet Drug Test Kits** are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
- All results are presumptive.** If necessary, further testing by a competent drug testing lab can provide additional verification.
- A color change indicates a positive result for the listed drugs **only if it occurs within the first minute.**

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.