



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

Test Preparation	1	Remove the Smart-Tip [™] probe from test unit.	2	Sample suspected substance.	
				Blotter (Stamp)	Tablet
				Roll and insert the blotter into the hole.	Firmly press Smart-Tip [™] probe on the tablet while scrubbing in a circular motion until the tip is completely covered with substance.

Test Preparation	3	Lock Smart-Tip [™] probe into test unit. Press firmly.	STOP Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit.	Identification	4	With thumb and forefinger , squeeze chamber 1 until capsule breaks.	5	Gently shake test unit for 20 seconds.

Identification	6	With thumb and forefinger , squeeze chamber 2 until capsule breaks.	7	Gently shake test unit for 2-3 seconds.


Results		The appearance of a purple-violet color within one minute indicates the presence of LSD (lysergic acid diethylamide) or DMT (dimethyltryptamine).
		If no violet-purple color develops within one minute, LSD and DMT are not present.

	Disposal	The test kit is safe to dispose of in the trash.
--	-----------------	--

IMPORTANT INFORMATION

- IDenta Touch&Know[™] Discreet Drug Test Kits** are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
- All results are presumptive.** If necessary, further testing by a competent drug testing lab can provide additional verification.
- A color change indicates a positive result for the listed drugs only if it occurs within the first minute.**

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.



Manufactured by IDenta Ltd
 www.touch-know.com
 info@touch-know.com