



As a precaution, gloves and eye protection should be worn when handling drug substances and using test units. Do not puncture any part of unit. Do not remove probe once a capsule is broken. In the unlikely event of a spill, rinse affected area thoroughly with water. Keep out of reach of children. All results are presumptive. If necessary, further testing by a competent drug testing lab can provide additional identification.

|                  |          |                                                                       |                                                                                                |        |
|------------------|----------|-----------------------------------------------------------------------|------------------------------------------------------------------------------------------------|--------|
| Test Preparation | <b>1</b> | Remove the Smart-Tip™ probe from test unit.                           | <b>2</b> Sample suspected substance.                                                           |        |
|                  |          |                                                                       | Powder                                                                                         | Liquid |
|                  |          | Swipe the paper of the Smart-Tip™ probe across the suspect substance. |                                                                                                |        |
|                  |          |                                                                       | Place the absorbent paper of the Smart-Tip™ probe into the liquid. Make sure the paper is wet. |        |

|                  |          |                                                     |                                                                                                                     |                |                                                                    |          |                                         |
|------------------|----------|-----------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|----------------|--------------------------------------------------------------------|----------|-----------------------------------------|
| Test Preparation | <b>3</b> | Lock Smart-Tip™ probe into test unit. Press firmly. | <b>STOP</b><br>Do not continue unless probe is locked below the Flex Pins and cannot be removed from the test unit. | <b>4</b>       | With thumb and forefinger, squeeze chamber 1 until capsule breaks. | <b>5</b> | Gently shake test unit for 2-3 seconds. |
|                  |          |                                                     |                                                                                                                     | Identification |                                                                    |          |                                         |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  |        |                                                                                                                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--------|---------------------------------------------------------------------------------------------------------------------|
| Results                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |  | Powder | The appearance of a dark blue color within one minute indicates the presence of <b>GHB</b> (γ-Hydroxybutyric Acid). |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  | Liquid | The appearance of an intense blue color on the paper within one minute indicates the presence of <b>GHB</b> .       |
| <p><b>Notes for liquid testing:</b></p> <ul style="list-style-type: none"> <li>* The test is more reliable for water, tea, black coffee, beer, vodka.</li> <li>* The detection limit is approximately 2.5 mg/ml of GHB.</li> <li>* The minimum detection limit for white and red wine is approximately 15mg/ml.</li> <li>* Water by itself with no GHB will produce a green color.</li> <li>* The kit is not reliable for use with milk or any drink with milk, because milk causes a blue color to appear.</li> <li>* The kit is not reliable for use with citrus drinks, such as orange or grapefruit juice.</li> </ul> |  |        |                                                                                                                     |
| <p> If <b>no</b> blue color develops within one minute or a color other than blue develops, such as green, orange or red, GHB is <b>not</b> present.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |  |        |                                                                                                                     |

**Disposal** The test kit is safe to dispose of in the trash. NET 07-21

### IMPORTANT INFORMATION

- IDenta Touch&Know™ Discreet Drug Test Kits** are used to determine the presence of certain drugs in substances such as powders, crystals, granules, flakes, pills, tablets, or liquids.
- All results are presumptive.** If necessary, further testing by a competent drug testing lab can provide additional verification.
- A color change indicates a positive result for the listed drugs only if it occurs within the first minute.**

If there is **NO** color change within the time indicated, the listed drugs are **NOT** present in the tested sample. After a minute, the sampler and the liquid may begin changing color as a normal part of the testing process and this color change does **NOT** indicate a positive result.