Dia Dent BIA-KONDENSOR

Featuring both a tapered Nickel-Titanium end and a stainless steel end, Dia-Kondensor can be used in any common obturation technique, including warm and cold condensation.

Purpose of use

Dia-Kondensor is that of vertically plugging and condensing the root canal filling material into the root canal.

Applications

- 1) Preparations prior to the use
- ① Before using Dia-Kondensor, check if it has been stored in a clean condition.
- ② Check for any noticeable external damage and make sure there is no dirt, dust or any foreign substance on the Product.
- 2) Application Method
- ① According to its intended use, choose the most appropriate endodontic instruments(type) and use it for reciprocating and vertical filling into the dental root canal.
- 2 Hold the product in your hand and apply the appropriate pressure to fill the canal with the filler.
- 3) Cleaning & Sterilization
- ① Cleaning

Soak the instrument in water to remove any chemical or other substances before wipping with alcohol. Repeat until the substance is completely removed.

② Sterilization

Before and after each use, steam autoclave the Dia-Kondensor. Recommended drying time after wet autoclaving is 20minutes at room temperature.

Storage & Maintenance

- 1) Store at room temperature
- 2) After use, Dia-Kondensor should be cleaned with alcohol to prevent foreign matters from drying out or to remove any dry matter.
- 3) No specific requirements as stainless steel or NiTi are used. However, it is preferable to store devices in a clean, dry place protected from damp.

Caution

- 1) Dia-Kondensor must be sterilized before use.
- 2) Do not apply too much force when inside a root canal
- 3) Do not use the Product for any other purpose.
- 4) It is the user's responsibility to sterilize products before first use and re-use.

Symbols



EC REP















Warranty and limitation of liability

DiaDent Group International will replace product that is proven to be defective. DiaDent Group International does not accept liability for any damage or loss direct or consequential, stemming from the use of or inability to use the product described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection

Contact Information DiaDent Group International 16,0songsaengmyeong 4-ro,0song-eup.

16, Osongsaengmyeong 4-ro, Osong-Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Korea

Tel: 82-43-266-2315 Fax: 82-43-262-8658 www.diadent.co.kr diadent@diadent.co.kr

DiaDent Group International Inc.

11-3871 North Fraser Way, Burnaby, B.C., Canada, V5J 5G6

Tel: 1-604-451-8851 Fax: 1-604-451-8865 www.diadent.com diadent@diadent.com

DiaDent Europe B.V.

www.diadenteurope.com diadent@diadenteurope.com



Copyright May 29,2019 [IFU : Dia-Kondensor (Rev.0)]