DIA ULTRASONIC
Ultrasonic Endodontic Activator

DiaDent DiaDent Group International

DiaDent Group International

16, Osongsaengmyeong 4-ro, Osong-eup,
Heungdeok-gu, Cheongju-si,
Chungcheongbuk-do, 28161, Korea
Tel : 82-43-266-2315
Fax : 82-43-262-8658
www.diadent.co.kr
diadent@diadent.co.kr

DiaDent Group International Inc.

11-3871 North Fraser Way, Burnaby
B.C., Canada, V5J 5G6
Tel : 1-604-451-8851
Fax : 1-604-451-8865
www.diadent.com
diadent@diadent.com

DiaDent Europe B.V.

Antennestraat 70, 1322AS
Almere, the Netherlands
Tel : 31-36-549-8607
Fax : 31-36-536-7317
www.diadent europe.com
diadent@diadent europe.com

REF Oct. 11, 2022(Rev.0)

DIA ULTRASONIC
Ultrasonic Endodontic Activator

Instruction Manual

DiaDent

©Copyright 2022 DiaDent Group International

DIA ULTRASONIC

Ultrasonic Endodontic Activator

Instruction Manual

English	4
français	20

DiaDent

Index

1. Overview	5
2. Precautions before use	6
2.1 Precautions for use	6
2.2 Storage Conditions	7
3. Components Checklist	8
4. Product Description	8
5. Components	9
6. Operation Method	10
6.1 Preparation before use	10
6.2 Instruction for use and operation procedures	10
6.3 Ultrasonic Tip Installation / Removal	11
6.4 Technique Instructions	11
7. Precautions	12
7.1 Handpiece	12
7.2 Charger and Battery	12
7.3 Ultrasonic Tip and Wrench	12
8. Cleaning, Sterilization and Maintenance	13
8.1 Maintenance	13
8.2 Cleaning and Sterilization Method	13
9. Disposal	14
10. Product Codes	14
11. Troubleshooting	15
12. Warranty	16
13. Specifications	16
14. Symbols	17

We sincerely thank our customers for choosing Dia-Ultrasonic. This product, with a simple design and a dedicated charger, is a dental Passive Ultrasonic Irrigation (PUI) device that cleans the root canal using ultrasonic. It is an innovative product that can effectively clean root canals through ultrasonic acoustic stream and cavitation. By applying dual frequency, powerful cleaning and detailed cleaning are all possible at once.

In addition, the dedicated ultrasonic tip improves durability and reduces the chance of fracture, allowing for longer use.

Special Features

- An ergonomic design for added convenience for use.
- Dual frequency : Powerful and thorough root canal irrigation
- A dedicated charger base for convenience and charging stability.
- A smart function that automatically selects the optimal frequency and soft start mode with a single push of a button.
- High-capacity battery (2,750mAh) : long operation time

This manual describes the detailed operation method and setting method. Please read this manual carefully before using Dia-Ultrasonic. Please keep this manual in an easily accessible place for quick and easy reference.

Indications for Use

Root canal irrigation plays a very important role in the success of root canal treatment. Dia-Ultrasonic is a dental PUI (Passive Ultrasonic Irrigation) device that cleans the root canal using ultrasonic energy (30KHz ~ 45KHz). It can prevent damage and deformation of the root canal, but also powerful enough to be able to effectively remove smear layer, dentin debris and bacteria using Acoustic Stream and Cavitation. When used together with a cleaning solution such as NaOCl, it disrupts biofilm and is antimicrobial and has antibacterial action.

2 Precautions before use

To reduce the risk of injury to the operator and patient, and to ensure an appropriate environment for product operation, please follow the manufacturer's instructions below.

CAUTION

Precautions informs you of information that can cause serious risk to the operator and patient and damage to the device if the device is not used correctly. The manufacturer is not responsible of any injury or damage if the user did not follow the user manual.

2.1 Precautions for use

1. Warning

- 1) Electromagnetic waves or other medical equipment used in hospitals may cause malfunction of the device or affect other medical equipment.
- 2) To reduce the risk of electric shock, do not open, disassemble or repair the handpiece or charger. Use only the dedicated charger and adapter provided by DiaDent.
- 3) Do not use in the presence of explosives and combustible materials.
- 4) Use only genuine components and parts provided by DiaDent. The use of non-genuine components and parts may cause product damage, accidents, or malfunction, and the manufacturer is not responsible.
- 5) If any malfunction or abnormality is found, stop using it immediately and contact the manufacturer.
- 6) Do not use this device on patients with pacemakers.

2. General Cautions

- 1) This product is intended for use only by trained and certified dentists or endodontic specialists for the intended purpose. If it is used by a user not mentioned in this manual or used for an unspecified purpose, the manufacturer is not responsible for any injury or damage.
- 2) Before use on a patient, the user must have a good understanding of this device and operating techniques.
- 3) To avoid the risk of patient injury, do not allow patients to touch or operate the product.
- 4) This product is not waterproof. Do not spray with or immerse the handpiece, charger, or adapter in any liquid.
- 5) For safe use, do not use continuously and follow the specified duty cycle (refer to 13. Specifications).
- 6) Keep solvents, flammable liquids, and sources of intense heat away from the device as they may damage the plastic housing of the device, the seals or the operating buttons.
- 7) Do not allow any cleaning agents to enter the device during cleaning as they could cause an electrical short or a dangerous malfunction.

3. Precautions for application

- 1) If the ultrasonic tip is damaged or fractured, discontinue use immediately.
- 2) Always use caution as there is a risk that fragments may enter the root canal if the ultrasonic tip is damaged.
- 3) Do not apply strong shocks or vibrations to the handpiece. Tighten the ultrasonic tip securely using the provided wrench. (Incorrect tightening may cause malfunctions and poor performance.)
- 4) Do not operate the device without any liquid. This will cause ultrasonic tip fracture.
- 5) Replace the ultrasonic tip with a new ultrasonic tip periodically. The recommended ultrasonic tip

replacement cycle is within 500 cycles or within 3 months.

- 6) Be sure to sterilize the ultrasonic tip before use to avoid contamination and cross infection. For sterilization method, refer to 8. Cleaning, Sterilization and Maintenance.
- 7) Do not touch the ultrasonic tip with your hands during operation. Heat can cause burns.
- 8) The tip of the ultrasonic tip is sharp. When moving the handpiece, be careful not to injure the body of the doctor, patient, or others, and be careful not to injure the patient's root canal during the procedure.

 Do not operate without liquid and do not touch the ultrasonic tip during operation. It may cause burns.

2.2 Storage conditions

- 1) For the best use of the equipment, it must be stored under the following environmental conditions. It may cause product malfunction or fire.
 - Store away from any liquid.
 - Do not store in a place where there is a risk of being adversely affected by sunlight, dust, salt, or air containing sulfur.

[Proper storage conditions]

▶ Atmospheric Pressure: 70 ~ 106KPa ▶ Temperature: -10 ~ 40℃ ▶ Humidity: 30 ~ 85%

- 2) Be sure to clean it thoroughly so that it does not interfere with the next use.
- 3) If any malfunction or abnormality is found, stop using it immediately and contact the manufacturer.
- 4) When not in use for a long time, do not store the battery in a discharged state. Store the charger and handpiece separated and unplug the power adapter.

This product has been designed to provide adequate protection from harmful interference to the operation and installation of general medical devices as a device that has been tested for electromagnetic compatibility in accordance with IEC60601-1-2. This product generates, uses, and can radiate radio frequency energy and may cause harmful interference to other devices in the vicinity.

If the use of this product does cause harmful interference to other devices, the user is encouraged to try to correct the interference by one or more of the following measures:


- Change the location of the receiver.
- Maintain a distance between this product and the device causing the interference.
- Contact the place of purchase or manufacturer.
- The battery of this product has completed the battery transportation test (UN DOT.38.3) and battery safety test (IEC 62133-2).

CAUTION

Potential electromagnetic waves or effects of other medical devices used in hospitals/clinics may cause malfunctions or adversely affect other medical devices.

3 Components Checklist

Make sure you have received all the products listed below.
If the product is missing or damaged, contact the place of purchase immediately.

- | | |
|---|---|
| 1 Handpiece 1EA | 5 Power Supply (adapter) 1EA |
| 2 Charger 1EA | 6 Power Cord 1EA |
| 3 Ultrasonic Tips 2EA(in the ultrasonic tip case) | 7 Battery pack 1EA (inside the handpiece) |
| 4 Wrench 1EA | 8 User Manual 1EA |


4 Product Description

Dia-Ultrasonic is luxurious, yet simple to use. It has an ergonomic design and a built-in rechargeable battery that can be recharged using the dedicated charger base and cradle. Various protection circuits (overcharge and over discharge protection circuits) are built-in for safety. LED light on the handpiece changes color to indicate the battery level.


Dia-Ultrasonic is used for root canal treatment using ultrasonic activation. Dia-Ultrasonic generates vibration with ultrasonic frequency required to create sufficient acoustic streaming and cavitation necessary for effective cleaning, penetration, and removal of vapor locks. A thoroughly cleaned root canal system provides improved treatment outcomes. This device should only be used by licensed dentists.

5 Components


■ Handpiece

■ Charger

■ Wrench

6 Operation Method

6.1 Preparation before use


- 1) Please read the instruction manual thoroughly.
- 2) Check the appearance of the product for abnormalities such as deformation or damage.
- 3) Check the battery charge.
- 4) Inspect the handpiece and components for defects.

6.2 Instruction for use and operation procedures

- 1) Connect the adapter and power cord, insert the plug into the outlet.
- 2) Connect the DC plug to the charger.
- 3) Charging starts when the handpiece is placed in the charger. Indicator lamp will light up in orange when charging, and green when charging is complete.
- 4) After removing the handpiece from the charger, secure the ultrasonic tip to the handpiece using the provided wrench. Tighten until the ultrasonic tip no longer spins. Do not over-tighten.
- 5) Press and hold the Operation Button for about 1 second to turn on the power. When the button is pressed in the standby mode, it turns on immediately.
- 6) Check that the Battery Level Indicator LED light is green.
- 7) Fill pulp chamber with NaOCl, EDTA, or other final rinse solution and insert the ultrasonic tip.
- 8) When the Operation Button is pressed, the operation status LED light turns on and operate for 15 seconds.
- 9) After operation, when the Operation Button is pressed or after 15 seconds, the operation status LED turns off and operation stops.
- 10) Charge when the orange LED on the battery level indicator lights up. Charge immediately when the orange LED flashes.
- 11) After use, remove the ultrasonic tip and clean it. For cleaning instructions, refer to 8. Cleaning, Sterilization and Maintenance.


⚠ For safety, do not use continuously, follow the specified duty cycle (refer to 13. Specifications).

[Adapter, Charger Connection]


6.3 Ultrasonic Tip Installation / Removal

⚠ Never screw or unscrew ultrasonic tip while the unit is activated.


- 1 Hand thread tip onto the handpiece. Insert the wrench into the ultrasonic tip joint.


- 2 Turn (clockwise) within the range of the wrench's maximum operating angle (about 25°). The wrench is designed to move only within a maximum of about 25°. Force should only be applied within the designed angle.

⚠ Do not overtighten.


- 3 When the force is released, wrench will automatically revert to its original position.
* Turn counterclockwise for ultrasonic tip removal.


6.4 Technique Instructions

- 1 Shape the root canal to be cleaned.

- 2 Fill the pulp cavity with cleaning agent / irrigation solutions (NaOCl, EDTA, etc.).

- 3 Insert the ultrasonic tip into the handpiece.

- 4 Tighten the ultrasonic tip with the provided wrench.

- 5 Press the start button to start root canal irrigation.

- 6 Move ultrasonic tip up and down using a small (2-3mm) vertical motion, maintaining a distance of 2mm from working length.

- 7 Activate the irrigation solution with the device 2-3 times for optimal result.

- 8 Suction out any remaining debris from the canals.

⚠ Avoid using the device in severely curved canals due to the risk of ultrasonic tip fracture.

⚠ Avoid hitting the ultrasonic tip against the root canal wall.

7 Precautions

7.1 Handpiece


- 1) Do not press the button while in use.
- 2) Do not operate in the atmosphere(air).
- 3) Be careful as the cleaning liquid may splash. If the ultrasonic tip is wet with cleaning liquid, wipe it off before use and storage.
- 4) Charge the device when the orange LED on the battery level indicator lights up. When the orange LED is blinking, the battery level is very low and the device needs to be charged immediately to restore normal functionality, as low battery level can adversely affect the performance.
- 5) The handpiece goes into sleep mode after 10 minutes when the device is not in use and turns off automatically after 60 minutes.
- 6) Do not short-circuit the handpiece's charging contacts. There is a risk of product damage and fire.
- 7) Do not overtighten the ultrasonic tip. The ultrasonic tip and body joint may be damaged.

7.2 Charger and Battery


- 1) Use only the power adapter supplied by DiaDent. If a third-party adapter is used, the internal parts may be damaged and may not operate normally and may cause a fire.
- 2) When the handpiece is removed from the charger, the LED turns off after a few seconds. If the handpiece is removed from the charger while charging, the green LED turns on and then turns off.
- 3) If the LED light alternates between green and orange during charging, the charging contacts are unstable. Remove the handpiece from the charger and put it in again to charge.
- 4) If the battery level is very low while charging, the yellow green LED light may light up. As the battery level increases, orange light will come on.
- 5) Do not short-circuit the charging contacts in the charger. There is a risk of product damage and fire.
- 6) When the device is not in use for a long period of time, separate the handpiece and charger and store them after charging.
- 7) Do not drop the product or subject it to impact. Doing so may cause material leakage, damage or fire in the battery installed inside the product.
- 8) Do not disassemble the handpiece case arbitrarily to expose the battery pack to the outside.

7.3 Ultrasonic Tip and wrench

- 1) After fastening the ultrasonic tip, the angle should be within 10° left and right. If it is not tightened enough, the ultrasonic output may be weakened. If it is tightened too much, it may damage the fastener and cause patient discomfort during the procedure.
- 2) If the angle exceeds 10°, remove the end Cover of the wrench and adjust the tension lightly(Use an M5 Hexagon/Allen Key, clockwise - strong tension, counterclockwise - weak tension).

- 3) When tightening the ultrasonic tip with the wrench, force should be applied only until the moment the head moves. If excessive force is applied, the clamping angle may be out of the standard range and the joint may be damaged.

8 Cleaning, Sterilization and Maintenance

Users should regularly check the device.

8.1 Maintenance

- 1) Before use
 - Check the handpiece and other parts for damage.
 - Inspect the ultrasonic tip for damage and, if present, replace it with a new one.
 - Check the ultrasonic tip for contamination, if any, clean and sterilize.
- 2) After use
 - Check that there are no loose parts in the connection part before tightening.
 - After removing the ultrasonic tip from the handpiece, clean the outer surface of the handpiece using a soft cloth / cotton ball / cotton pad soaked with alcohol.
- 3) Yearly Maintenance
 - You must carefully check that your device is safe and functioning properly. It is necessary to check that there are no abnormalities.

8.2 Cleaning and Sterilization Method

- 1) Handpiece
 - ① Remove the handpiece from the charger for cleaning.
 - ② Clean the outer surface of the handpiece using a soft cloth / cotton ball / cotton pad soaked with alcohol.
 - ③ Be careful not to immerse the handpiece in the cleaning solution.
- 2) Charger
 - ① Remove the handpiece from the charger.
 - ② Disconnect the charger from the adapter and outlet.
 - ③ Clean the outer surface of the charger with a soft cloth / cotton ball / cotton pad soaked with alcohol.
 - ④ Be careful not to immerse the charger in any cleaning solution.
- 3) Ultrasonic Tips
 - ① Remove the ultrasonic tip from the handpiece.
 - ② Use an alcohol swab to remove all foreign substances from the ultrasonic tip and autoclave. Autoclave at 121°C for 30 minutes. After autoclaving, dry at room temperature for 30 minutes. The ultrasonic tip can be reused up to 500 times or up to 3 months.
- 4) Wrench
 - ① Check for any damage.
 - ② If there is dirt on the surface of the wrench, wipe it with an alcohol swab.

9 Disposal

1) Handpiece and charger disposal

- Remove the battery pack before disposal. Dispose of handpieces and chargers as electronic waste. Please dispose of it in accordance with your local electronic product disposal regulations.

2) Battery Pack Disposal and Precautions

- Please dispose of the battery pack as chemical waste. Dispose of in accordance with appropriate local regulations.
- When disposing of the battery pack, do not incinerate or dispose of it in a fire. Heat can cause explosion and fire.
- Dispose of the battery in a discharged state to prevent heat generation of the battery due to short circuit.
- Do not disassemble, open or cut the battery.
- Do not expose the battery to heat or fire.
- Beware of short circuits in the battery.

3) Ultrasonic Tips

- Remove the ultrasonic tip from the handpiece before disposal.
- Dispose of ultrasonic tips as medical waste.
- Dispose of the ultrasonic tips in accordance with appropriate local regulations.

10 Product Codes

- Dia-Ultrasonic Regular Kit : 4009-1001
- Dia-Ultrasonic Hand-piece : 4009-1110
- Dia-Ultrasonic Charger : 4009-1120
- Dia-Ultrasonic Tip : 4009-1130
- Dia-Ultrasonic Wrench : 4009-1140
- Dia-Ultrasonic U-battery(Battery Pack) : 4009-1150
- Dia-Ultrasonic Power Supply (Adapter) : 4009-1160
- Dia-Ultrasonic Power Cord : 4009-1170

11 Troubleshooting

Here are some common problems you may encounter and how to solve them. For problems not presented here, please contact the place of purchase or the A/S Center of the head office.

Problem	Solution
1. Power does not turn on	<ul style="list-style-type: none"> - Make sure the charger and power adapter are working properly. - Charge the battery pack.
2. Handpiece does not work	<ul style="list-style-type: none"> - Make sure the battery pack is fully charged before use. - Check the charger and power adapter. - Request service from the place of purchase or manufacturer.
3. Handpiece is operating abnormally	<ul style="list-style-type: none"> - Check if the operation button is pressed. - Press and hold the operation button for more than 12 seconds to turn off the power, then turn it on again. - Request service from the place of purchase or manufacturer.
4. When the output is abnormally weak	<ul style="list-style-type: none"> - Make sure the battery is fully charged. - Make sure the ultrasonic tip is properly tightened. - Retighten using the provided wrench.
5. Charger does not work	<ul style="list-style-type: none"> - Make sure the power adapter is the one provided by the manufacturer. - Make sure the power adapter is plugged into an outlet with the correct voltage rating. - Check the charger and power adapter for damage. - Check that the handpiece is properly placed in the charger. - Do not use the device if the battery or charger of the handpiece is wet. Turn off the power immediately and dry it completely. - If the charger still does not work even after checking the above solutions, contact the place of purchase or A/S center.
6. When the battery drains quickly (Short battery life)	<ul style="list-style-type: none"> - Lithium-ion batteries have a long lifespan, however, performance can still deteriorate when used for a long period of time. - If a problem occurs during use, the battery should be replaced with a new one. - We recommend replacing the battery with a new one every two years. - The battery is not user replaceable. Request service from the place of purchase or manufacturer.
7. When the ultrasonic tip(needle part) is broken	<ul style="list-style-type: none"> - Inject root canal irrigation solution into the root canal and suck the ultrasonic tip with suction. - After adding the cleaning solution and installing a new ultrasonic tip, vibration to move the broken ultrasonic tip out of the root canal and remove or suck up with suction.









12 Warranty

- 1) DiaDent provides a 1-year free warranty and repair of this device (1 year from the date of purchase) for defects in performance and accessories when used in compliance with the instructions and precautions in the product manual. (However, ultrasonic tips and battery packs are consumables and are not covered under warranty.)
- 2) To obtain warranty coverage, please contact the place of purchase to ship the product to the manufacturer's service center.
- 3) Product shipments for product warranty benefits must be securely packaged and the shipping cost must be prepaid. Manufacturers do not accept products shipped on a postpaid basis. The manufacturer is responsible for repair and replacement of the product.
- 4) The warranty is void in the following cases :
 - Failure caused by not using the product in the manner described in the user manual
 - Failure caused by negligence in handling
 - Failure due to natural disaster
 - Failure caused by not using the specified power source (Do not use the charging adapter with other products)
 - If the product is disassembled and repaired

13 Specifications

Product Name		Ultrasonic Endodontic Activator
Model Name		Dia-Ultrasonic
Ultrasonic Frequency		30KHz ~ 45KHz
Charging Time		Approximately 3 hours 30 minutes
Duty Cycle		Operation Time : up to 15 seconds Break Time : at least 5 seconds
Operation Time per Full Charge		For 15 seconds up to 500 times
Dimensions	Hand-piece	178 x 25 x 25mm(L x W x H)
	Charger	91 x 91 x 44mm(L x W x H)
	Ultrasonic Tip	27 x 32 x 4mm(L x W x H)
	Wrench	74 x 12 x 12(L x W x H)
Weight	Hand-piece	101g
	Charger	135g
	Ultrasonic Tip	0.8g
	Wrench	15g
Power Supply (Adapter)	Input	100-240V~, 50/60Hz, 0.4A
	Output	5V  2A
Battery	Voltage, Capacity	Li-ion 3.65V  , 2750mAh / Rechargeable
Operating Conditions	Temperature(°C)	10 ~ 30 °C
	Ambient Pressure	70 ~ 106 (KPa)
	Relative Humidity	30 ~ 85 %
Transportation and Storage Conditions	Temperature(°C)	-10 ~ 40 °C
	Ambient Pressure	70 ~ 106 (KPa)
	Relative Humidity	30 ~ 85 %

14 Symbols

Symbol	Description	Symbol	Description
	Cautions indicate which, if not observed, could result in damage to the user and/or device		Direct current
	CE certificate: Product conformance with the applicable European Union Directives		Authorized EU Representative
	Class II Equipment		Type B Applied Part
	Manufacturer		Manufacture Date
	Serial No		Subject to electrical/electronic equipment waste regulations
	See instruction manual		Sitting Prohibited
	Humidity Limitation		Temperature Limitation
	Catalog Number		Atmospheric Pressure Limitation
	Caution, hot Surface Be careful to prevent burns		Keep Dry
	The Restriction of the use of certain Hazardous Substances in electrical and electronic equipment		Medical Device
	Refer to accompanying Documents		Refer Servicing to Qualified personal only
	Prescription Only		

EMC (Electromagnetic Compatibility) Information

Phenomenon	Basic EMC standard or test method	Port tested	Test level/requirement
Conducted Emissions (Mains terminal disturbance voltage)	EN 60601-1-2:2015 IEC 60601-1-2:2014 EN 55011:2016+A11:2020 CISPR11:2015+A1:2016+A2:2019	AC Mains	Group1, Class A
Radiated disturbance	EN 60601-1-2:2015 IEC 60601-1-2:2014 EN 55011:2016+A11:2020 CISPR11:2015+A1:2016+A2:2019	Enclosure	Group1, Class A
Harmonic Current Emission	EN IEC 61000-3-2:2019 IEC 61000-3-2:2018	AC Mains	Class A
Voltage change, Voltage fluctuations and Flicker Emission	EN 61000-3-3:2013+A1:2019 IEC 61000-3-3:2013+A1:2017	AC Mains	Max Pst : 1 Plt : 0.65 Max Tmax : 0.5 Max dmax : 4% Max dc : 3.3% General test conditions
Electrostatic Discharge	EN 61000-4-2:2009 IEC 61000-4-2:2008	Enclosure	Standard B Air Discharges (2, 4, 8, 15) kV Contact Discharges (2, 4, 6, 8) kV Positive(+) and Negative(-)
Radiated RF Electromagnetic Field	EN 61000-4-3:2006+A1:2008+A2:2010 IEC 61000-4-3:2006+A1:2007+A2:2010	Enclosure	Frequency Range : 80 MHz~2.7 GHz and 3 V/m Modulation : AM 80% with 1 KHz sine wave
Proximity Fields from RF wireless Communications Equipment	EN 61000-4-3:2006+A1:2008+A2:2010 IEC 61000-4-3:2006+A1:2007+A2:2010	Enclosure	Range of Frequency(MHz) : 385, 450, 710, 745, 780, 810, 870, 930, 1720, 1845, 1970, 2450, 5240, 5500, 5785 Polarity : Horizontality and Verticality
Electrical Fast Transient/Burst	EN 61000-4-4:2012 IEC 61000-4-4:2012	AC & DC Power, Signal In/ Out port	AC&DC Power : ± 2 kV, Signal In/Out : ± 1 kV, Repetition Frequency : 100 kHz
Surge Transient	EN 61000-4-5:2014+A1:2017 IEC 61000-4-5:2014+A1:2017	AC & DC	Line to Line : ± 0.5 kV, ± 1 kV Line to Earth : ± 0.5 kV, ± 1 kV, ± 2 kV
Conducted Disturbances Induced by RF fields	EN 61000-4-6:2014 IEC 61000-4-6:2013	AC Mains	Frequency Range : (0.15 ~ 80) MHz Test Voltage : 3.0 V(General), ISM and amateur radio bands : 6.0 V Modulation : AM 80% with 1 KHz sine wave

Power Frequency Magnetic Field	EN 61000-4-8:2010 IEC 61000-4-8:2009	Enclosure	30 A/m 50 Hz or 60 Hz
Voltage Dips and Interruptions	EN IEC 61000-4-11:2020 IEC 61000-4-11:2020	AC Mains	Voltage Reduction : 0.5 Cycle to 100% : C, 1 Cycle to 100% : C, 25/30 Cycles to 30% : C, 250/300 Cycles to 100% : C Number of pulses : 3 at each level Recovery time between pulses : 10 s

NOTE

- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals(CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- If there is any serious incident that occurs while using the device, the incident should be reported to DiaDent Group International, the related authorities and/or the dealer where the purchase was made.

Manufacturer :

DiaDent Group International

16, Osongsaengmyeong 4-ro, Osong-eup,
Heungdeok-gu, Cheongju-si,
Chungcheongbuk-do, 28161, Korea
Tel : 82-43-266-2315
Fax : 82-43-262-8658
www.diadent.co.kr
diadent@diadent.co.kr

Distributor :

DiaDent Group International Inc.

11-3871 North Fraser Way, Burnaby
B.C., Canada, V5J 5G6
Tel : 1-604-451-8851
Fax : 1-604-451-8865
www.diadent.com
diadent@diadent.com

DiaDent Europe B.V.

Antennestraat 70, 1322AS
Almere, the Netherlands
Tel : 31-36-549-8607
Fax : 31-36-536-7317
www.diadent europe.com
diadent@diadent europe.com

EC REP Authorized representative in the European Community

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